SOAP3-dp: Fast, Accurate and 1 Sensitive GPU-based Short Read 2

Aligner 3

- Ruibang Luo^{1,*}, Thomas Wong^{1,*}, Jianqiao Zhu^{1,4,*}, Chi-Man Liu^{1,*}, Edward Wu¹, 4
- Lap-Kei Lee¹, Haoxiang Lin², Wenjuan Zhu², David W. Cheung¹, Hing-Fung Ting¹, 5
- 6 Siu-Ming Yiu¹, Chang Yu², Yingrui Li², Ruigiang Li³, Tak-Wah Lam¹
- 7
- 8 ¹HKU-BGI Bioinformatics Algorithms and Core Technology Research Laboratory & Department
- ğ of Computer Science, University of Hong Kong, Hong Kong.
- 10 ² BGI-Shenzhen, Shenzhen, Guangdong 518083, China.
- 11 ³ Peking-Tsinghua Center for Life Sciences, Biodynamic Optical Imaging Center and School of Life
- 12 Sciences, Peking University, China.
- 13 ⁴ Department of Computer Sciences, University of Wisconsin-Madison, USA.

liyr@genomics.cn

lirq@pku.edu.cn

twlam@cs.hku.hk

- 14 * These authors contributed equally to this work. 15
- 16 Correspondence should be addressed to T. L. (twlam@cs.hku.hk).
- 17 18 **Ruibang Luo** rbluo@cs.hku.hk
- 19 Thomas Wong thomaskf@gmail.com
- 20 Jiangiao Zhu jiangiao@cs.wisc.edu
- 21 Chi-Man Liu cmliu@cs.hku.hk
- 22 Edward Wu edwardmkwu@gmail.com
- 23 Lap-Kei Lee lklee@cs.hku.hk
- 24 Haoxiang Lin linhaoxiang@genomics.cn
- 25 26 Wenjuan Zhu zhuwenjuan@genomics.cn dcheung@cs.hku.hk
- David W. Cheung
- 27 Hing-Fung Ting hfting@cs.hku.hk
- 28 Siu-Ming Yiu smyiu@cs.hku.hk yuchang@genomics.cn
- 29 Chang Yu
- 30 Yingrui Li
- 31 Ruiqiang Li
- 32 Tak-Wah Lam 33

Abstract: 1

2

3 To tackle the exponentially increasing throughput of Next-Generation

4 Sequencing (NGS), most of the existing short-read aligners can be configured to

- 5 favor speed in trade of accuracy and sensitivity. SOAP3-dp, through leveraging
- the computational power of both CPU and GPU with optimized algorithms, 6 7 delivers high speed and sensitivity simultaneously. Compared with widely
- 8 adopted aligners including BWA, Bowtie2, SeqAlto, GEM and GPU-based aligners
- 9 including BarraCUDA and CUSHAW, SOAP3-dp is two to tens of times faster,
- while maintaining the highest sensitivity and lowest false discovery rate (FDR) 10
- on Illumina reads with different lengths. Transcending its predecessor SOAP3, 11
- which does not allow gapped alignment, SOAP3-dp by default tolerates 12
- 13 alignment similarity as low as 60%. Real data evaluation using human genome
- demonstrates SOAP3-dp's power to enable more authentic variants and longer 14
- 15 Indels to be discovered. Fosmid sequencing shows a 9.1% FDR on newly
- 16 discovered deletions. SOAP3-dp natively supports BAM file format and provides
- 17 a scoring scheme same as BWA, which enables it to be integrated into existing
- 18 analysis pipelines. SOAP3-dp has been deployed on Amazon-EC2, NIH-Biowulf
- 19 and Tianhe-1A.
- 20

21 Keywords: Genome Alignment; GPU Acceleration; Burrows-Wheeler Transform;

22 **Dynamic Programming**

Maintext: 23

24

- 25 Introduction: 26 With the rapid advancement of Next-Generation Sequencing technologies, 27 modern sequencers like Illumina HiSeq 2500 can sequence a human genome into 28 600 million pairs of reads of 100bp in length (total 120 Gigabases) in merely 27 29 hours. The cost is also decreasing fast. By 2013 year's end, sequencing a human 30 genome is projected to cost less than \$1,000. Bioinformatics research using 31 sequencing data often starts with aligning the data onto a reference genome, 32 followed by various downstream analyses. Alignment is computationally 33 intensive; the 1000 genomes pilot paper[1] published in 2010 reported that a 34 1192-processor cluster was used to align the reads using MAO[2]. This kind of 35 computing resources is not available to most laboratories, let alone clinical 36 settings. Although considerable advances have been made on new aligners. 37 alignment still remains a bottleneck in bioinformatics analyses. Thus, ultra-fast 38 alignment tools without relying on extensive computing resources are needed. 39 40
- There are quite a few software tools for aligning short reads onto a reference
- 41 genome. The more popular ones include MAQ, Bowtie[3], BWA[4] and SOAP2[5].
- 42 The faster ones[3-5] index the reference genome using the Burrows-Wheeler
- 43 Transform (BWT), which is efficient for aligning short reads with limited
- mismatches, but inefficient for aligning with gaps. These tools (running on a 44
- 45 quad-core processor) take tens of hours to align 120 Gigabases with limited (or

1 even none for Bowtie and SOAP2) gapped alignment found. Alignment gaps can

2 result from insertions and deletions (Indels), which are thought to comprise over

3 20% of genetic variations[6] and contribute to human traits[7]. Hence, a

4 successor is expected to be faster and more sensitive to gaps.

5

SegAlto[8] and GEM[9] were published very recently. SegAlto is a hash-based 6 7 aligner that improves an earlier hash-based aligner SNAP[10] (reported to have 8 relatively poor sensitivity for real reads and provide no mapping quality[8]) 9 using additional global and local alignments. SegAlto is slower than SNAP, vet SeqAlto is still faster than the BWT-based aligners except Bowtie2. GEM mapper 10 leverages string matching with filtration to search the alignment space more 11 efficiently[9]. GEM is faster than comparable state-of-the-art aligners. Yet it 12 13 does not provide PHRED[11] compliant mapping quality score; this saves some 14 tedious computation, but prohibiting it from integrating into existing analysis 15 pipelines.

16

17 Nowadays, general-purpose computing on graphics processing units (GPUs) is 18 becoming popular. A GPU is a piece of low-cost hardware providing massive 19 parallelism but with limited memory and restricted usage. A number of GPU 20 bioinformatics tools emerged last year[12]. BarraCUDA[13] implements BWA to 21 align reads in parallel on a GPU; limited by the branch and bound trie algorithm 22 that requires extensive decisions making, BarraCUDA works sub-optimally on 23 GPU and gains a 4-time boost than a single-thread BWA. SOAP3[14] successfully 24 exploits the massive parallelism of a GPU with tailor-made GPU-BWT and read-25 characteristics sensitive load balancing to effectively align short reads. Albeit not 26 supporting gapped alignment, which makes it unsuitable for production, it is to 27 date the most competitive aligner for ungapped alignment.

28

29 Here we present a GPU-based alignment software SOAP3-dp that allows multiple 30 mismatches and gaps, which is well suited for production environments (real 31 data alignments) than the predecessor SOAP3. A simple approach to extend 32 mismatch alignment to gapped alignment is to first identify candidate regions by 33 exact or mismatch alignment of short substrings (seeds) in the reads, then use 34 dynamic programming to perform a detailed alignment of the read to the regions. Such an approach has been widely used (e.g., Bowtie2). The bottleneck occurs as 35 substring alignment often results in a large number of candidates, especially 36 37 when mismatches are allowed. As a result, reads with too many candidates are 38 often ignored due to time constraint. On the other hand, the parallelism of GPU 39 apparently would allow many candidates to be verified in parallel; yet dynamic 40 programming is memory consuming, and the limited-memory of GPU becomes a 41 prohibiting factor to fully utilize the parallelism. SOAP3-dp gives a pragmatic 42 realization of this approach (Figure 1, Supplementary Note). By exploiting 43 compressed indexing and memory-optimizing dynamic programming on a GPU, 44 SOAP3-dp can efficiently tackle a large number of candidates in parallel, and thus 45 can examine gapped alignments extensively and achieve a drastic improvement 46 in both speed and sensitivity over other tools.

47

48 See the Supplementary Note about the design of the dynamic programming 40 which attempts to minimize the memory usage for each candidate so as to let a

49 which attempts to minimize the memory usage for each candidate so as to let a

1 GPU to handle hundreds of candidates in parallel while using limited shared

2 memory. The note also shows the details of SOAP3-dp's intricate engineering

3 solution to finding the optimal way to align different reads using either the CPU

4 or the GPU.

56 Experiments & Performance:

We compared SOAP3-dp to other short-read alignment software on the speed,
sensitivity and accuracy. We used both real and simulated Illumina data.

Sensitivity and accuracy, we used both real and simulated illumina data.
 Furthermore, we tested out SOAP3-dp's alignment quality for variant calling

rui thermore, we tested out SOAP3-up's alignment quality for variant calling
 using real data. In particular, 41-fold of 100bp (PE100) and 77-fold of 150bp

11 (PE150) Illumina paired-end reads of YH[15] sample have been generated

12 (Supplementary Table 1) for the testing for variant calling.

13

Alignment Performance. We first used real data to test SOAP3-dp with BWA
(BWA-SW[16] for SE reads), Bowtie2, SeqAlto, GEM, BarraCUDA and CUSHAW.
The aim was to compare the time and alignment rate when each runs in the
default setting. Next, to assess the accuracy and sensitivity, we used simulated
reads whose correct alignments were known. We then considered more detailed
comparison with the software running in different settings. We also attempted to
compare SOAP3-dp against its predecessor SOAP3.

21

22 In our experiments, we assume input reads are plain text instead of in gzip-23 compressed format. This is because GEM (to the date of paper submission) does 24 not accept gzip-compressed FASTQ file. All other software can handle 25 compressed input, which is getting common nowadays. Regarding output format, 26 we require all software to use SAM format, which is mandatory for downstream 27 analysis software including GATK[17] and SAMTOOLS[18]. All software except 28 GEM can output directly in SAM format; GEM first outputs in a simple format and 29 then takes an extra step to convert to SAM format. To test GEM's efficiency when 30 using a simple format, we include a comparison to SOAP3-dp also using a simple 31 format (see the Remark section).

32

33 Real data. We used three real datasets for benchmarking of alignment 34 performance: (1) a lane (122.43M reads) from PE100, and (2) a lane (374.87M 35 reads) from PE150, and (3) SRR211279 (25.23M 100bp paired-end reads Generated by Illumina GAIIx) from the Washington University Genome 36 37 Sequencing Center. We tested SOAP3-dp and six other aligners (namely, BWA, 38 Bowite2, SeqAlto, GEM, Barracuda, and CUSHAW; see Supplementary Note for 39 receipts), all using 4 CPU threads and 1 GPU device (For GPU aligner). As shown 40 in Table 1, SOAP3-dp is much faster than all others (Supplementary Tables 2-9 41 for more details). It is at least 3.5 times faster than GEM, and 7 to 15 times faster 42 than the other five. SOAP3-dp also gave better alignment rate consistently. 43 SeqAlto comes closest, aligning 0.48% to 3.6% less reads than SOAP3-dp, and the 44 others are in the range of 2% to 8% less than SOAP3-dp. Notice that except 45 SOAP3-dp, aligners usually have an obvious drop in alignment rate for longer 46 reads (dataset 2). For SOAP3, its alignment rate, as expected, is much poorer 47 than SOAP3-dp and the others (this is due to lack of gapped alignment); furthermore, SOAP3 is slower than SOAP3-dp for 100bp reads. The two GPU-48 49 based aligner Barracuda and CUSHAW are rather primitive in optimizing GPU's

1 utilization and overheads, and their performance was dominated by new CPU-

2 based aligners like GEM and Bowtie2. We did not include Barracuda and

3 CUSHAW for further experiments.

4

Simulated data. To assess the accuracy and sensitivity of SOAP3-dp, we used the
short read simulator Mason[19] to obtain 5 sets of 6M Illumina-style paired-end
(PE) reads with 500bp insert size from GRCh37 major build, with length ranging
from 50-250bp.

9

Notably, Bowtie2, SegAlto and GEM were designed with switches to favor speed 10 at the expense of accuracy and sensitivity. We applied "very-fast", "sensitive", 11 and "very-sensitive" switches to Bowtie2, "fast (-f)" to SeqAlto and "fast adaptive 12 13 (--fast-mapping)", "fastest (--fast-mapping=0)" to GEM. For SOAP3-dp, we tested three versions whose indices are based on 1/4 sampled, 1/2 sampled and full 14 15 suffix array (SA), respectively. Different sized SAs still deliver identical alignment 16 results, but a smaller one consumes less memory at the expense of slightly 17 longer alignment time. All parameters of SOAP3-dp and SOAP3 remained as 18 default while parameters for other aligners were set to favor different read types

and lengths as suggested by previous studies. Detailed command lines were

- summarized as receipts in the Supplementary Note. In total, 15 sets of programsand parameters were compared.
- 22

23 In all datasets of simulated reads, SOAP3-dp gives consistent advantage. It is 24 faster and simultaneously has higher sensitivity and lower FDR over all other 25 tools (Table 2 for 100bp PE, Supplementary Tables 10-13 for all other simulated 26 datasets, Figure 2). For 100bp reads, SOAP3-dp with full SA takes 132 seconds 27 to align 6M read pairs, and it is 2.26 to 8.74 times faster than the others (others 28 using the fastest switches). SOAP3-dp's sensitivity is 99.96%, which is 0.13 -29 0.85% higher than the others (others using the sensitive switches), and SOAP3-30 dp's FDR is 0.34%, which is lower than the others by 0.13 - 0.85%. Apparently 31 the simulated data is easier to align than the real data due to recombination 32 hotspots with intensified variants in real genome [6]. SOAP3-dp consumes more 33 memory (9.3, 11.9, 17.2GB for 1/4, 1/2 and full SA in average, respectively) than 34 other software; Bowtie2 has the least (3.5GB). Nevertheless, workstations and 35 servers nowadays are equipped with at least 16GB or even 32GB of memory; 36 SOAP3-dp is designed to take advantage of the available memory to achieve

- 37 speed.
- 38

Mapping quality score is mandatory for most of the popular downstream
analysis tools such as GATK and SAMTOOLS. SOAP3-dp uses a scoring scheme
same as BWA so as to make its alignment results compatible to the expectations
of existing analysis tools. As shown in Figure 3, BWA, Bowtie2, SeqAlto and
SOAP3-dp provide mapping quality scores that can differentiate different
alignments properly, while GEM's scores are too rough, and too many incorrectly
aligned reads are given high quality scores, which makes it unsuitable for

- 46 downstream analysis.
- 47

48 **Remarks.** Note that GEM can output in a simple format to save time. When
49 compared to SOAP3-dp in its own simple format, GEM and SOAP3-dp can both

- 1 save a lot of time; for the 6M simulated paired-end data of length 100bp, the
- 2 alignment time of GEM and SOA3-dp is reduced to 90 seconds and 38 seconds,
- 3 respectively. Downstream programs, if redesigned to utilize these specific
- 4 formats, could save time.
- 5
- 6 The simulated dataset is relatively small (6M read pairs), thus when using a large
- 7 SA, the index loading time of SOAP3-dp dominated the total elapsed time.
- 8 Considering only the alignment time (time consumption after index loading,
- 9 including input of raw reads and output of alignment results), SOAP3-dp using
- 10 the full and 1/2 SA is 12 and 9 seconds faster than 1/4 SA, thus for large real
- 11 datasets, 1/2 and full SAs are suggested if memory permits.
- 12

13 Different generations of GPU device differ in speed. We compared the

- 14 performance of SOAP3-dp between the latest GPU "GTX680" and a previous
- 15 generation "Tesla C2070" using simulated datasets. The alignment time extended
- 16 about 10% for each dataset (Figure 4) using "Tesla C2070". Furthermore, a large
- 17 real set of 150bp paired-end reads was used. The alignment using the "GTX680"
- consumed 6,835 seconds, which is 4,658 seconds (1.68 times) faster than "Tesla
 C2070".
- 20

21 Variant Calling Performance:

- Next, we considered SOAP3-dp's alignment quality for variant calling. The full sets of both PE100 and PE150 were aligned using SOAP3-dp. We used widely adopted BWA as benchmark. With the alignment results, variants were called using GATK's UnifiedGenotyper[20] and filtered by VariantRecalibrators, with and without GATK's local realignment (see Methods). Before we detail the results on variant detection, it is worth mentioning that BWA, even running in a slower mode to allow a longer gap (one gap up to 50bp, without "-m" parameter
- to allow hit entries higher than 2 million due to out of memory error), still
- 30 cannot catch up the sensitivity of SOAP3-dp in default setting (for a 100bp read,
- 31 one gap up to 68bp, to 23 one-bp gaps, Table 3).
- 32
- SOAP3-dp's better sensitivity is due to its ability of extensive gapped alignment;
 the extra reads aligned are crucial for variant detection (in particular, Indels).
- 35 SOAP3-dp allowed 2.4% and 4.0% more SNPs than BWA, and 6.1% and 9.8%
- 36 more Indels for the two datasets PE100 and PE150, respectively. Intuitively,
- 37 longer reads are more favorable for variant detection; this is indeed reflected in
- 38 SOAP3-dp's performance, but not for BWA (Table 3). We further checked the
- 39 SNPs detected against dbSNP v135 (an archive of SNPs validated by previous
- 40 studies); SOAP3-dp has notably 2.1% and 3.6% more SNPs found in dbSNP,
- 41 confirming a higher sensitivity.
- 42

SOAP3-dp allowed more Indels to be detected than BWA, especially more long
Indels (Figures 5a, b, Supplementary Figure 1). To validate the novel Indels
detected, we randomly selected 50 SOAP3-dp exclusively detected deletions not
yet archived in dbSNP v135 and verified them using Fosmid sequencing (see
Methods). The Fosmids were sequenced, assembled and then aligned to the
reference genome. The 50 deletions were covered by 460 Fosmid sequences. As

49 a result, 6 deletions were inconclusive due to insufficient coverage of Fosmid

- 1 sequences, 40 deletions were validated, and 4 rejected, revealing a FDR of 9.1%
- 2 (Supplementary Tables 14, 15). SOAP3-dp's ability to allow long gaps without
- 3 speed penalty provides an unprecedented opportunity to come up with a more
- 4 comprehensive Indel identification in large-scale genome studies.
- 5
- 6 With SOAP3-dp's ability to authentically align more reads, more multi-nucleotide
- 7 polymorphisms (MNP) among the whole genome were identified
- 8 (Supplementary Figure 2). Notably, GATK's local realignment can eliminate
- 9 inauthentic alignments and rescue true variants. For SOAP3-dp, the number of
- 10 MNP increased by 4.1% after realignment; yet for BWA, the number decreased
- 11 by 6.3% (Supplementary Table 16), indicating that SOAP3-dp initially provided
- 12 much more reliable alignments and led to more accurate variant calling.
- 13

14 **Discussion**:

- 15 SOAP3-dp has been successfully deployed on Amazon EC2, NIH BioWulf and
- 16 Tianhe-1A computing-cloud. On Amazon EC2, users can access SOAP3-dp's
- 17 program and a testing dataset by mounting EBS snapshot "snap-154f1c54"
- 18 named "SOAP3-dp" while creating a GPU instance (Supplementary Note). To test
- 19 out SOAP3-dp's performance on Amazon EC2, we selected 10 Illumina PE
- 20 datasets from 1000 genomes project, comprising 131.44Gbp of raw reads (43.8-
- fold). The 10 datasets was distributed to the two available Tesla M2050 GPU
- 22 cards (see Methods and Supplementary Note. Notably, Tesla M2050 is slower
- than the GTX680 and Tesla C2070 we have used for real and simulated data
- evaluation) with one copy of the index shared in host memory. Using default
- parameters and BAM output, the alignment finished in 3.8 hours, yielding a total
 cost of \$7.98, or \$0.061 per Gbp reads aligned.
- 27
- 28For users' convenience, SOAP3-dp separates the output of reads into three
- 29 categories: 1) alignments involve gaps and extensive mismatches; 2) few
- 30 mismatches only and 3) improperly paired or unaligned (file suffix "dpout",
- "gout", and "unpair" respectively). The separate file scheme fits well with the
 production environments, where files could be sorted separately in parallel and
 then merged together, which saves time than sorting a single SAM file. The files
- 34 could also be concatenated by SAMTOOLS easily.
- 35

SOAP3-dp does not enforce a maximum read length. However, read length longer
than 500bp is not recommended while the current version of SOAP3-dp is tailormade for Illumina reads. A version for longer 454 reads and Ion Torrent reads
without performance degradation is out next task.

- 40
- 41 Overall, SOAP3-dp is an efficient alignment tool that targets the future of genome
 42 analysis where reads are longer and the volume is larger. SOAP3-dp is much
- 43 faster than existing tools while retaining the ability to align more reads correctly.
- 44 To be flexible, SOAP3-dp outputs both SAM and BAM formats that are compatible
- 45 with most downstream analysis tools. SOAP3-dp is free, open-source software
- 46 available at http://www.cs.hku.hk/2bwt-tools/soap3-dp/.

Acknowledgements 1

- 2 Tak-Wah Lam was partially supported by RGC General Research Fund 10612042.
- 3 Ruibang Luo was supported by Hong Kong ITF Grant GHP/011/12.

Author Contributions 4

- 5 Ruib. L., T. W., J. Z. and C. L. contributed equally to this work. T. L. led the project
- and contributed the main ideas. Ruig. L., Y. L., D. C., H. T. and S. Y. co-managed the 6
- 7 project. J. Z., T. W., C. L., Ruib. L., E. W., L. L. wrote the program. Ruib. L. and T. W.
- 8 designed the analysis. Ruib. L. and T. W. conducted the simulation analysis. Ruib.
- L., H. L. and W. Z. conducted the real data analysis. C. Y. assisted in deploying the 9
- 10 program. T. L., Ruib. L. and T. W. wrote the paper.

Competing Interests 11

12 The authors have declared that no competing interests exist.

Tables 13

- 14
- 15 Table 1. Benchmarking using real reads. The percentage of reads aligned and
- time consumption aligners other than SOAP3-dp are recorded as the difference 16
- and ratio based on SOAP3-dp's figures. The '%' column represents 'Properly 17
- 18 paired' for PE reads.
- 19

Dataset	Volume (Gbp)	SOAP3-dp		BWA or BWASW		Bowtie2		SeqAlto		GEM		BarraCUDA		CUSHAW		SOAP3	
		%	Time (s)	94	Time	%	Time	%	Time	%	Time (fold)	%	Time	%	Time	%	Time
				70	(fold)		(fold)		(fold)	70			(fold)		(fold)	70	(fold)
realYHPE100	12.24	98.12%	1,079	-4.14%	16.03	-3.87%	12.04	-0.48%	14.25	-2.86%	3.54	-6.39%	14.03	-3.13%	46.79	-10.60%	2.63
realYHPE150	56.23	97.16%	6,835	-8.05%	15.26	-7.45%	7.82	-3.65%	17.69	-5.85%	3.76	-10.99%	7.22	-10.11%	24.57	-28.33%	0.64
SRR211279	5.07	97.21%	439	-2.81%	16.27	-0.69%	12.00	-0.48%	17.29	-2.44%	4.11	-5.08%	13.42	-2.22%	54.22	-6.83%	1.05

20 21

22 Table 2. Comparison on 15 sets of programs and parameters using 100bp paired-23 end simulated reads.

24

					(SPU			CPU								
6M 100bp Paired-end reads, 1.2Gbp bases. 500bp insert size, 25bp standard deviation.			SOAP3-dp ³ (1/4 SA)	SOAP3-dp ¹ (1/2 SA)	SOAP3-dp ¹ (Full SA)	SOAP3-dp ¹ (Full SA, Binary)	SOAP3-dp ¹ (Full SA, Succinct)	SOAP3	Bowtie2 (Sensititve)	Bowtie2 (Very- Sensititve)	Bowtie2 (Very-fast)	BWA ²	SeqAlto	SeqAlto (Fast alignement)	GEM ³	GEM ³ (Fast Mapping: adaptive)	GEM ³ (Fast Mapping: 0)
Configuration CPU (thread: core i7-3930k)		4			4	4	4	4	4	4	4	4	4	4	4	4	
comguration	GPU (device: GTX680)		1			1	1	1	0	0	0	0	0	0	0	0	0
Computational Resources	Total Elapsed	sec.	132	137	162	111	112	132	966	1974	672	1154	495	379	416	446	298
		Fold	-	1.04	1.23	0.84	0.85	1.00	7.32	14.95	5.09	8.74	3.75	2.87	3.15	3.38	2.26
	Loading Index ⁴	sec.	32	46	74	74	74	74	38	38	38	53+1+1	96	96	40+1	40+1	40+1
	Alignment ⁵	sec.	100	91	88	37	38	58	928	1936	634	370+369+360	399	283	199+176	238+167	90+167
		Fold	-	0.91	0.88	0.37	0.38	0.58	9.28	19.36	6.34	10.99	3.99	2.83	3.75	4.05	2.57
	Avg. Memory	GB	9.3	11.9	17.2	17.2	17.2	17.3	3.3	3.3	3.3	3.5	7	6.9	4.3	4.3	4.3
	Peak Memory	GB	9.7	12.5	18.1	18.1	18.1	19.2	3.5	3.5	3.5	4.8	7.2	7.2	4.3	4.3	4.3
	Aligned	#	11,999,827					11,870,740	11,999,763	11,999,936	11,998,226	11,998,804	12,000,000	11,995,872	11,999,763	11,999,484	11,995,422
		Diff.			-			-129,087	-64	109	-1,601	-1,023	173	-3,955	-64	-343	-4,405
	Properly Paired	#			11,999,460			11,742,902	11,998,912	11,999,344	11,996,528	11,997,254	11,999,976	11,995,410	11,998,994	11,997,702	11,991,992
		Diff.			-			-256,558	-548	-116	-2,932	-2,206	516	-4,050	-466	-1,758	-7,468
Alignment Metrics	Incorrectly Aligned	#			40,561			138,655	143,012	141,373	147,764	85,297	95,672	99,243	56,514	61,642	61,887
		Diff.			-			98,094	102,451	100,812	107,203	44,736	55,111	58,682	15,953	21,081	21,326
	Sensitivity ⁶	%			99.66%			97.77%	98.81%	98.82%	98.75%	99.28%	99.20%	99.14%	99.53%	99.48%	99.45%
	Sensitivity	Diff.			-			-1.89%	-0.85%	-0.84%	-0.91%	-0.38%	-0.46%	-0.52%	-0.13%	-0.18%	-0.21%
	EDR7	%			0.34%			1.17%	1.19%	1.18%	1.23%	0.71%	0.80%	0.83%	0.47%	0.51%	0.52%
	- DR	Diff.			-			0.83%	0.85%	0.84%	0.89%	0.37%	0.46%	0.49%	0.13%	0.18%	0.18%

¹ Alignment results by the three entries of SOAP3-dp (1/4 SA, 1/2 SA, Full SA) are identical. ² The time consumption of BWA is calculated as "align left reads"+"align right reads"+"sampe". The index loading times of "align right reads" and "sampe" modules are 1 second due to the reason that, index files were cached during "align left reads". However, datasets larger than the host memory will flush the cache during alignment.

³ The alignment time consumption of GEM is calculated as "alignment"+" convert to SAM format ". The

conversion module was run with 4 threads in consistent with the alignment module.

⁴ SOAP3-dp, SOAP3, SeqAlto and GEM aligners explicitly provide index loading time consumption. The index loading time for Bowtie2 and BWA are calculated by the total size of index, divided by 100MB/s, which is

the average network file system speed of the testing environment. The index loading time maybe

underestimated while the time processing the index was not calculated.

1 2 3 4 ⁵ The alignment times were explicitly provided by the aligners (include results processing and input/output time) or calculated by total elapsed time minus estimated index loading time.

⁶ Sensitivity is calculated as "Correctly aligned reads"/"All simulated reads". The higher the better.

⁷ FDR is calculated as "Incorrectly aligned reads"/"All aligned reads". The lower the better.

5

6 Table 3. Summary of alignment and variation calling using SOAP3-dp and BWA

7 with different parameters and datasets. PE100 and PE150 represent the 100bp

8 and 150bp paired-end reads of YH sample. 'w/' and 'w/o' indicates whether the

9 alignments were processed with and without GATK's local realignment,

- 10 respectively. 'dbSNP' is the number of SNPs already archived in dbSNP v135. The
- 11 running times were rounded to hour.
- 12

	Running time (h)		Poads	Properly		Indole						
Program / Type			aligned			w/o			w/	indels		
	Aln.	GATK	aligheu	Falleu	Total	dbSNP	Homozygous	Total	dbSNP	Homozygous	w/o	w/
BWA (Default, PE100)	99	32	93.30%	92.10%	3,996,730	3,926,590	1,620,269	3,922,841	3,876,262	1,609,899	370,231	383,347
BWA (-o 1 -e 50 -L, PE100)	3363	36	96.74%	94.80%	4,040,813	3,955,380	1,609,729	3,946,627	3,896,737	1,601,049	385,216	417,832
SOAP3-dp (Default, PE100)	6	35	98.03%	96.98%	4,120,165	4,028,414	1,615,376	4,015,725	3,956,371	1,604,106	407,533	406,845
BWA (Default, PE150)	266	46	90.30%	89.17%	3,963,085	3,897,866	1,621,742	3,911,386	3,862,977	1,616,767	350,201	384,684
SOAP3-dp (Default, PE150)	14	92	97.66%	96.54%	4,139,625	4,055,314	1,612,213	4,067,552	4,002,750	1,612,351	409,925	422,427

13 14

1 Figure legends

2

3 Figure 1, Alignment workflow.

- 4 For each read (paired-end specifically, single-end is only with step 1 and step 3),
- 5 the alignment would be decided in at most three steps. In step 1, SOAP3-dp
- 6 aligns both ends of a read-pair to the reference genome by using GPU version
- 7 2way-BWT algorithm (Supplementary Note). Pairs with only one end aligned
- 8 proceed to step 2 for a GPU accelerated dynamic programming (Supplementary
- 9 Note) alignment at candidate regions inferred from the aligned end. Pairs with
- both ends unaligned in step 1 and those ends failed in step 2 proceed to step 3 to
- 11 perform a more comprehensive alignment across the whole genome until all 12 acad bits (substrings from the read) are supplied as with a sufficient such as a supplied of the supplier of the supplice of t
- 12 seed hits (substrings from the read) are examined or until a sufficient number of
- 13 alignments are examined.
- 14
- 15 Figure 2,
- 16 Speed and sensitivity of alignment using simulated paired-end reads. We
- 17 recorded the number of correct and incorrect alignments stratified by reported
- 18 mapping quality for each dataset. We then calculated the cumulative number of
- 19 correct and incorrect alignments from high to low mapping quality. We
- 20 considered an alignment correct only if the leftmost position was within 50bp of
- 21 the position assigned by the simulator on the same strand according to the
- 22 previous study of Bowtie2 to avoid soft-clipping artifacts.
- 23
- 24 Figure 3,
- The accumulated number of incorrectly aligned reads categorized at differentmapping quality scores by the five aligners.
- 27
- 28 Figure 4,
- Alignment time consumption of using GPU card "GTX680" and previous
- 30 generation GPU card "Tesla C2070" respectively.
- 31
- 32 Figure 5,
- **a.** Indel length distribution of indels smaller than or equal to 20bp, **b.** larger than
- 20bp identified by SOAP3-dp and BWA respectively using full set of 100bp
- 35 paired-end YH sample reads.

36 Method:

37 Simulation of single-end and paired-end reads:

- 38 Mason 0.1[19] was used to simulate reads using the GRCh37 major build human
- 39 reference genome, including 22 pairs of autosomes, 2 sex chromosomes and a
- 40 mitochondrial chromosome . For the paired Illumina-style datasets with read
- 41 length ranging from 50bp to 250bp, Mason was run in 'Illumina' read mode with
- 42 options -N 6000000 --source-no-N -mp -sq -ll 500 -le 25 -rn 2 -hn 2 --haplotype-
- 43 snp-rate 0.001 --haplotype-indel-rate 0.0001 --haplotype-no-N -n 100 -pi 0 -pd 0

- 1 --no-N'. Each PE set was simulated to contain exactly 12 million reads (6M pairs),
- 2 which is about 2 times the default batch size of SOAP3-dp.
- 3
- 4 All simulated datasets are available at <u>http://bio8.cs.hku.hk/dataset/</u>.

5 Simulation comparison between tools:

6 Executable files for SOAP3-dp v2.3, Bowtie2 v2.0.0-beta4, BWA 0.6.2, SeqAlto 7 basic 0.5-r123, BarraCUDA r232, CUSHAW-1.0.40, SOAP3 version146 and GEM-8 core_i3-20121106-022124 were obtained via standard build procedures with 9 default arguments. We indexed the reference genome with each tool's default 10 indexing parameters. SeqAlto uses 22bp seed length and sub-sampled mode. "Running time" was measured from initial call of the aligner to the completion of 11 12 SAM-format output. 'Reads aligned' was measured as the number of reads for 13 which the tool found at least one alignment regardless of mapping score. 14 'Properly paired' was measured as the number of read-pairs aligned with proper read orientation and insert-size range (mean insert-size ± 3*standard deviation). 15 16 'Peak memory' and 'Average memory' usage was measured by tracking the 17 Linux's proc file-system with respective process id. 18 For BWA and BarraCUDA, separate calls of the software modules were required 19 for aligning each end and for processing intermediate alignment results into a 20 final SAM file. 'Running time' was measured for separate modules while 'Peak 21 memory' and 'Average memory' were measured across all modules. For GEM, to 22 be consistent with other tools, a conversion is necessary after alignment to 23 obtain SAM format results with mapping quality. Time consumptions were 24 measured separately and then summed for comparison to other tools. All tools 25 or components were run with 4 threads (except for the alignment module of 26 BarraCUDA, where the CPU thread is constantly 1). Parameters were listed in 27 Supplementary Note as receipts. 28 The experiments used a single computing node running CentOS v6.3 with an 29 Intel i7-3930k 3.2Ghz quad-core processor, an Nvidia GTX 680 GPU card with 30 4GB non-ECC (Error-correcting code) graphic memory and 64GB non-ECC

- 31 memory.
- 32

33 Scripts and command lines to evaluate the authenticity of aligned reads and

34 generate the ROC curves are available at <u>http://bio8.cs.hku.hk/dataset/.</u>

35

36 YH data production:

37 Genomic DNA was isolated using standard molecular biology techniques. For

38 each short insert library, 5 µg of DNA was fragmented, end-repaired, A-tailed and

- ligated to Illumina paired-end adapters. The ligated fragments of 100bp pairedend reads (PE100) were size-selected at 170bp and 500bp on agarose gels, while
- 41 PE150 were size-selected at 240bp. All libraries are amplified by LM-PCR to yield
- 42 the corresponding short insert libraries. PE100 were sequenced using TruSeq v2
- 43 while PE150 were sequenced using TruSeq v3 reagent on the Illumina
- 44 sequencing platform.
- 45

1 Real data comparison:

The 100bp and 150bp paired-end Illumina HiSeq 2000 reads of YH sample were
sequenced and deposited to EBI with study accession number ERP001652. The
data are also available at http://yh.genomics.org.cn.

5

6 SOAP3-dp uses default parameters. BWA uses both default parameters and "-o 1 7 -e 50", which allows at most a gap not longer than 50bp (-m option to elevate the 8 2M hits limit for each read was not applied due to out of memory error, the 9 option allows more reads to be aligned but consumes much more memory and 10 longer alignment time). The latter option allows more reads to be aligned and more indel signals to be discovered, but would enormously decrease the running 11 12 speed. Alignments were post-processed by following procedures: 1) local 13 realignment by GATK v2.1, 2) duplication removal by Picard v1.74, 3) base 14 quality score recalibration, 4) variants calling by UnifiedGenotyper and 5) 15 variants quality score recalibration by GATK v1.6. Step 1 is optional according to the experiment while steps 2 to 5 are mandatory. Parameters and known variant 16

- 17 sets were set according to the GATK's Best Practice v4 on GATK's website.
- 18

19 SOAP3-dp used a single node with a quad-core Intel Xeon E5570 2.93Ghz CPU

 $20 \qquad \text{and a GPU while BWA used 10 nodes with the same CPU. To imitate the real} \\$

21 production environment, we used Nvidia Tesla C2070 GPU device with 6G

22 graphic memory and with ECC enabled to perform full YH dataset alignment.

23 Fosmid sequencing:

24 Fosmid libraries (averagely 40kbp in size) were constructed according to Kim et 25 al. [21]. In total, ~100k Fosmid clones were created and every 30 Fosmids were pooled together sharing a barcode for Illumina HiSeq 2000 sequencing. For each 26 pool, one 200bp and one 500 bp insert size libraries were constructed and 27 28 sequenced at 20x respectively. If a library had problem of abnormal base content 29 bias or a relatively high base error rate reported by base-calling software, we 30 took it as a non-qualified library and performed sequencing again. Each pool 31 was assembled with SOAPdenovo[22] using 63-mer and other parameters as 32 default. Sequences solved by SOAPdenovo's repeat solving module were 33 remembered and soft-masked in final sequences in order not to obscure the 34 following alignment. The assembled sequences were aligned to the human 35 reference genome using BWASW with default parameters. Most of the Fosmid 36 clones could be assembled to full length. For fragmented Fosmid clone sequences, 37 we further assembled the fragments according to the in-pool linkage information 38 during alignment. We require over 90% of a Fosmid clone sequence to be 39 linearly aligned to only one location in the reference genome. Finally, 460 40 Fosmid clones were found covering the 50 randomly selected SOAP3-dp specific 41 deletion calls that are not yet archived in dbSNP v135. We define "a Fosmid 42 sequence supporting a deletion" as over 80% of the deleted bases in reference genome cannot be covered by the aligned Fosmid sequence (excluding soft-43 44 masked bases), and the identity of the 200bp alignments flanking the deletion 45 should exceed 90%. While a Fosmid clone can only come from a haploid, we 46 require a heterozygous deletion has at least a Fosmid supporting the deletion, 47 while a homozygous deletion should only have Fosmids supporting the deletion. Heterozygous deletions with lower than 5 spanning Fosmids and without a 48

- 1 Fosmid supporting the deletion will be classified as "not clear". Homozygous
- 2 deletions without a spanning Fosmid will also be classified as "not clear"
- 3 (Supplementary Tables 14, 15).
- 4
- 5 The assembled Fosmid sequences are available as Supplementary Data in BAM
- 6 file format. Raw reads are available upon request.

7 Experiment on Amazon EC2:

- 8 The experiment used a single GPU Quadruple Extra Large Instance (cg1.4xlarge,
- 9 \$2.1 per hour) rented from the Amazon Elastic Compute Cloud (EC2) service
- 10 (http://aws.amazon.com/ec2). The instance has 2 quad-core Intel Xeon X5570 at
- 11 2.93GHz with hyper-threading, 2 Nvidia Tesla M2050 GPU cards with 3GB ECC
- 12 graphic-memory per card, 22 gigabytes of physical memory and runs Amazon
- 13 Linux AMI v2012.09 operating system. Alignments were distributed onto the two
- 14 GPU cards with two SOAP3-dp processes sharing the same copy of index in host
- 15 memory. Each process occupies at most 7 threads. 10 sets of Illumina HiSeq
- 16 2000 reads generated in 1000 genomes project (Supplementary Note) were
- 17 downloaded from the Amazon Simple Storage Service (S3). Additional tests have
- 18 been carried out by NIH biowulf laboratory (<u>http://biowulf.nih.gov/apps/bioinf-</u>
- 19 gpu.html) and Tianhe-1A super-computing center
- 20 (<u>http://en.wikipedia.org/wiki/Tianhe-I</u>). While these clouds host CPU and GPU
- 21 computing nodes and centralized storage system, the ultra-fast SOAP3-dp could
- 22 be easily integrated into existing pipelines.
- 23

Reference: 1

- 2 1. (2010) A map of human genome variation from population-scale sequencing.
- 3 Nature 467: 1061-1073.
- 4 2. Li H, Ruan J, Durbin R (2008) Mapping short DNA sequencing reads and calling
- 5 variants using mapping quality scores. Genome Res 18: 1851-1858.
- 3. Langmead B, Trapnell C, Pop M, Salzberg SL (2009) Ultrafast and memory-6
- 7 efficient alignment of short DNA sequences to the human genome. Genome 8 biology 10: R25.
- 9 4. Li H, Durbin R (2009) Fast and accurate short read alignment with Burrows-
- 10 Wheeler transform. Bioinformatics 25: 1754-1760.
- 11 5. Li R, Yu C, Li Y, Lam TW, Yiu SM, et al. (2009) SOAP2: an improved ultrafast 12 tool for short read alignment. Bioinformatics 25: 1966-1967.
- 13 6. Li Y, Zheng H, Luo R, Wu H, Zhu H, et al. (2011) Structural variation in two
- 14 human genomes mapped at single-nucleotide resolution by whole genome de 15 novo assembly. Nat Biotechnol 29: 723-730.
- 7. Mullaney JM, Mills RE, Pittard WS, Devine SE (2010) Small insertions and 16
- 17 deletions (INDELs) in human genomes. Hum Mol Genet 19: R131-136.
- 18 8. Mu JC, Jiang H, Kiani A, Mohiyuddin M, Bani Asadi N, et al. (2012) Fast and
- 19 accurate read alignment for resequencing. Bioinformatics 28: 2366-2373.
- 20 9. Marco-Sola S, Sammeth M, Guigo R, Ribeca P (2012) The GEM mapper: fast,
- 21 accurate and versatile alignment by filtration. Nat Methods 9: 1185-1188.
- 22 10. Zaharia M, Bolosky WJ, Curtis K, Fox A, Patterson D, et al. (2011) Faster and
- 23 More Accurate Sequence Alignment with SNAP. arXiv 1111.5572v1.
- 24 11. Ewing B, Green P (1998) Base-calling of automated sequencer traces using
- 25 phred. II. Error probabilities. Genome Res 8: 186-194.
- 26 12. Nvidia (2012) Bioinformatics and Life Sciences.
- 27 13. Klus P, Lam S, Lyberg D, Cheung MS, Pullan G, et al. (2012) BarraCUDA - a fast
- 28 short read sequence aligner using graphics processing units. BMC Res Notes 5: 27.
- 29
- 30 14. Liu CM, Wong T, Wu E, Luo R, Yiu SM, et al. (2012) SOAP3: ultra-fast GPU-
- 31 based parallel alignment tool for short reads. Bioinformatics 28: 878-879.
- 32 15. Wang J, Wang W, Li R, Li Y, Tian G, et al. (2008) The diploid genome sequence
- 33 of an Asian individual. Nature 456: 60-65.
- 34 16. Li H, Durbin R (2010) Fast and accurate long-read alignment with Burrows-
- 35 Wheeler transform. Bioinformatics 26: 589-595.
- 17. McKenna A, Hanna M, Banks E, Sivachenko A, Cibulskis K, et al. (2010) The 36
- 37 Genome Analysis Toolkit: a MapReduce framework for analyzing next-
- 38 generation DNA sequencing data. Genome Res 20: 1297-1303.
- 39 18. Li H, Handsaker B, Wysoker A, Fennell T, Ruan J, et al. (2009) The Sequence
- 40 Alignment/Map format and SAMtools. Bioinformatics 25: 2078-2079.
- 41 19. Doring A, Weese D, Rausch T, Reinert K (2008) SeqAn an efficient, generic
- 42 C++ library for sequence analysis. BMC Bioinformatics 9: 11.
- 43 20. DePristo MA, Banks E, Poplin R, Garimella KV, Maguire JR, et al. (2011) A
- 44 framework for variation discovery and genotyping using next-generation DNA
- 45 sequencing data. Nat Genet 43: 491-498.

- 21. Kim UJ, Shizuya H, de Jong PJ, Birren B, Simon MI (1992) Stable propagation of cosmid sized human DNA inserts in an F factor based vector. Nucleic Acids Res 20:1083-1085.
- 22. Li R, Zhu H, Ruan J, Qian W, Fang X, et al. (2010) De novo assembly of human
- genomes with massively parallel short read sequencing. Genome Res 20: 265-
- 272.

Step 1: Use 2way-BWT to align pair-ended reads











simPE50

simPE75

simPE100



simPE150

simPE250



Indel length distribution (≤20bp)

Indel length distribution (>20bp)

