# CLUSTER ANALYSIS OF HIGH-DIMENSIONAL SCRNA SEQUENCING DATA

Jiawei Long UCLA Fielding School of Public Health University of California, Los Angeles Los Angeles, CA 90024 peterljw@g.ucla.edu Yu Xia UCLA Fielding School of Public Health University of California, Los Angeles Los Angeles, CA 90024 xiayu960112@g.ucla.edu

December 19, 2019

### ABSTRACT

With ongoing developments and innovations in single-cell RNA sequencing methods, advancements in sequencing performance could empower significant discoveries as well as new emerging possibilities to address biological and medical investigations. In the study, we will be using the dataset collected by the authors of *Systematic comparative analysis of single cell RNA-sequencing methods*. The dataset consists of single-cell and/or single nucleus profiling from three types of samples – cell lines, peripheral blood mononuclear cells, and brain tissue, which offers 36 libraries in six separate experiments in a single center. Our quantitative comparison aims to identify unique characteristics associated with different single-cell sequencing methods. Our procedures also incorporate evaluations of every method's capacity for recovering known biological information in the samples through clustering analysis.

## 1 Introduction

With the provided dataset of aggregated count matrices of the three samples, we will preprocess each aggregated count matrix into sub count matrices that correspond to separate replicates and sequencing methods. After preprocessing the scRNA sequencing data, we will quantitatively compare and evaluate various methods used on each sample replicate by assessing metrics such as dropout rate and sensitivity measures. Lastly, we will perform various dimensionality techniques along with a number of clustering methods on PBMC and cortex samples to investigate each utilized method's capacity for recovering known cell types in the samples.

# 2 Data Understanding and Preprocessing

The dataset provides three data files for each sample, which includes the aggregated read count matrix for both replicates and all methods, as well as the gene annotation and cell annotation for the count matrix. The first stage of the study involved dividing the aggregated count matrix into a number of sub count matrices which correspond to different replicates, and sequencing methods and annotating rows with cell names and columns with gene names.

Among the whole dataset, seven scRNA sequencing methods are present. There are five high-throughput methods, 10x-Chromium-v2, 10x-Chromium-v3, inDrops, Drop-seq, and sci-RNA-seq, alongside two low-throughput methods, Smart-seq2 and CEL-Seq2. By observing the dimensions of the count matrices, low-throughput methods such as Smart-seq2 and CEL-Seq2 tend to be associated with a relatively low number of cells. Figure 1 displays a summary of the count matrices generated after the division of the count matrix for each of the three samples.

Sample	Replicate	Sequencing Method	Pre-filtering Dimension	Post-filtering Dimension
Cortex	Cortex 1	Smart-seq2	28692 x 295	4691 x 1480
		Drop-seq	28692 x 2195	2049 x 2195
		sci-RNA-seq	28692 x 1886	2086 x 1886
		10x-Chromium-v2	28692 x 1480	8485 x 295
	Cortex 2	Smart-seq2	28692 x 349	2484 x 4091
		Drop-seq	28692 x 935	2959 x 935
		sci-RNA-seq	28692 x 3944	1513 x 3944
		10x-Chromium-v2	28692 x 4091	7646 x 349
PBMC	PBMC 1	Smart-seq2	33694 x 311	716 x 5172
		CEL-Seq2	33694 x 257	1069 x 3057
		10x-Chromium-v2-A	33694 x 5172	1996 x 4033
		10x-Chromium-v2-B	33694 x 3057	3909 x 257
		10x-Chromium-v3	33694 x 4033	595 x 4683
		Drop-seq	33694 x 4683	223 x 6184
		Seq-Well	33694 x 5125	357 x 5125
		inDrops	33694 x 6184	3299 x 311
	PBMC 2	Smart-seq2	33694 x 273	1034 x 3362
		CEL-Seq2	33694 x 307	3722 x 307
		10x-Chromium-v2	33694 x 3362	876 x 6412
		Drop-seq	33694 x 6412	480 x 4526
		Seq-Well	33694 x 913	197 x 913
		inDrops	33694 x 4526	4017 x 273
Mixture	Mixture 1	Smart-seq2	62046 x 342	-
		Drop-seq	62046 x 2594	-
		sci-RNA-seq	62046 x 299	-
		10x-Chromium-v2	62047 x 3159	-
		inDrops	62048 x 3081	-
		Seq-Well	62049 x 1627	-
		CEL-Seq2	62050 x 359	-
	Mixture 2	Smart-seq2	62051 x 343	-
		Drop-seq	62052 x 3812	-
		sci-RNA-seq	62053 x 5023	-
		10x-Chromium-v2	62054 x 3187	-
		inDrops	62055 x 2529	-
		Seq-Well	62056 x 1012	-
		CEL-Seq2	62057 x 347	-

Figure 1: Count Matrices

# **3** Comparative Evaluations

#### 3.1 Dropout Rate

A dropout event refers to the occurrence when a transcript is expressed in a cell but is undetected in its mRNA profile. The reason behind the occurrence of dropout events is the low amounts of mRNA in individual cells. The frequency of dropout events depends on scRNA-seq protocols. There exists a general trade-off between the number of cells and the frequency of dropout events, whereas scRNA-seq protocols that generate more cells tend to have a higher frequency of dropout events. Figure 2 demonstrates a graphic illustration of dropout events.



Figure 2: Dimensions of Count Matrices

The dropout rate of a sequencing method indicates the level of the sparsity of a scRNA sequencing method. For each of the count matrix that we have obtained from data preprocessing, we calculated the dropout rate by dividing the number of cells with zero entries by the total number of cells. The outcomes are demonstrated in figure 3 to figure 5 below.

Method	Dropout Rate	Method	Dropout Rate
DroNc-seq	94.03%	sci-RNA-seq	95.31%
sci-RNA-seq	94.01%	10x-Chromium-v2	93.39%
10x-Chromium-v2	89.11%	DroNc-seq	92.50%
Smart-seq2	80.72%	Smart-seq2	82.88%

Figure 3: Dropout rates of Cortex 1 (left) and Cortex 2 (right)

Method	Dropout Rate		
inDrops	98.92%		
Seq-Well	98.32%	Method	Dropout Rate
Drop-seq	97.87%	Seq-Well	98.89%
10x-Chromium-v2-A	97.48%	inDrops	97.93%
10x-Chromium-v2-B	96.61%	Drop-seq	97.11%
10x-Chromium-v3	95.09%	10x-Chromium-v2	96.64%
Smart-seq2	92.77%	CEL-Seq2	92.38%
CEL-Seq2	92.13%	Smart-seq2	91.70%

Figure 4: Dropout rates of PBMC 1 (left) and PBMC 2 (right)

Method	Dropout Rate	Method	Dropout Rate
inDrops	97.48%	inDrops	96.53%
Drop-seq	95.85%	Drop-seq	95.95%
sci-RNA-seq	93.74%	sci-RNA-seq	92.99%
Seq-Well	93.35%	10x-Chromium-v2	92.56%
10x-Chromium-v2	92.93%	Seq-Well	92.37%
Smart-seq2	89.85%	CEL-Seq2	86.39%
CEL-Seq2	86.99%	Smart-seq2	85.97%

Figure 5: Dropout rates of Mixture 1 (left) and Mixture 2 (right)

According to the dropout rates table,

- low-throughput plate-based methods such as Smart-seq and Cel-seq are associated with the lowest dropout rates among all methods.
- high-throughput methods such as Drop-seq, Seq-well, and inDrops have significantly higher dropout rates.

#### 3.2 Sensitivity

The sensitivity of a scRNA-Seq method refers to the likelihood to capture and convert a particular mRNA transcript present in a single cell into a cDNA molecule present in the library.

To assess the sensitivity of each scRNA sequencing method, we have generated visualizations of two metrics, gene detection per cell and cumulative gene detection. Please refer to figure 6 to figure 9 for the visualisations.

- Gene detection per cell (left): The number of genes detected (counts ≥ 1) per cell. Each dot represents a cell and each box represents the median and first and third quartiles per replicate and method.
- Cumulative gene detection (right): The cumulative number of genes detected as more cells are added. The x-axis represents the number of cells accounted while the y-axis represents the number of genes detected.



Figure 6: Sensitivity Visualizations for Cortex 1



Figure 7: Sensitivity Visualizations for Cortex 2



Figure 8: Sensitivity Visualizations for PBMC 1



Figure 9: Sensitivity Visualizations for PBMC 2

It is noticeable that low-throughput methods appear to have high sensitivity, in which case they tend to have a higher number of genes detected per cell and they detected the most genes in an efficient manner. In contrast, Drop-seq is consistently associated with low sensitivity across all samples, displaying relatively low gene detection efficiency. Other high-throughput also displayed relatively lackluster performance in terms of sensitivity.

Based on the graphs which display the cumulative number of genes detected with additional cells added, high-throughput methods tend to converge to a relatively low number of detected genes despite increments in the number of cells, implying that there is a group of genes that are potentially difficult to capture with high-throughput methods.

# 4 Cluster Analysis

Sub-population identification, usually via some form of unsupervised clustering, is a fundamental step in the analysis of many single-cell RNA-seq data sets as it helps with identifying underlying gene types. In this section, we have conducted different clustering methods on the data after applying filtering, normalization, and dimensionality reductions.

## 4.1 Filtering

Although RNA-seq technology has improved the dynamic range of gene expression quantification, low-expression genes may be indistinguishable from sampling noise. The presence of noisy, low-expression genes can lead to poor analysis results. Thus, it is beneficial to remove the low-expression genes in the data before applying clustering techniques. Please refer to figure 1 for dimensions of count matrices after filtering.

We have taken the following two procedures to filter the data.

- Row sparsity filtering: Remove genes with more than 80% zero counts.
- CV filtering: Remove out the bottom 15% of the remaining genes which have the lowest coefficient of variation.

#### 4.2 Normalization

Normalization is an essential step in an RNA-Seq analysis, in which the read count matrix is transformed to allow for meaningful comparisons of counts across samples. In this case, we have applied Quantile normalization. Quantile

normalization is a technique for making distributions identical in statistical properties. Each column is converted to rank values before sorting and averaging the ranks to replace the previous ranks.

## 4.3 Dimensionality Reduction

To achieve better and more reasonable clustering visualization results on a two-dimensional space, we needed to apply dimensionality techniques to implement feature engineering.

**PCA** Principal component analysis is a linear dimensionality reduction technique, and it performs a linear mapping of the data to a lower-dimensional space in such a way that the variance of the data in the low-dimensional representation is maximized. The major procedure involves calculating the eigenvectors from the covariance matrix, where the eigenvectors that correspond to the largest eigenvalues are used to reconstruct a significant fraction of the variance of the original data. However, in the case of high dimensions, eigenvalues and eigen-vectors of the sample covariance matrix are not consistent, which may affect the performance of PCA.

**t-SNE** t-SNE is a non-linear technique for dimensionality reduction technique. It calculates the probability of similarity of points in high-dimensional space and calculating the probability of similarity of points in the corresponding low-dimensional space, and it minimizes the difference between these conditional probabilities (or similarities) in higher-dimensional and lower-dimensional space for a representation of data points in lower-dimensional space. However, t-SNE is computationally complex and non-deterministic, which makes the technique less robust.

**\*Projection Pursuit** Projection Pursuit is a component transform technique that looks for a component whose projection vector points to a direction of the quality of interest in data space which can be determined by a Projection Index (PI). In general, projections that go astray more from a normal distribution are viewed as more interesting. As each projection is found, the data are reduced by removing the component along that projection, and the process is repeated to find new projections. The idea of projection pursuit is to locate the projection or projections from high-dimensional space to low-dimensional space that reveal the most insights regarding the structure of the data set. However, one potential issues associated with the technique is that it generally uses random initial conditions to produce projection index components. As a result, when the same projection pursuit is performed on different occasions or by different users, the resulting projection index components are generally not the equivalent. This approach has yet to be implemented in our study.

## 4.4 Clustering Methods

**K-means Clustering** K-Means clustering intends to partition n objects into k clusters in which each object belongs to the cluster with the nearest mean. This method produces exactly k different clusters of greatest possible distinction. K-Means clustering minimizes the squared error function with a given k.

**Hierarchical Clustering** Hierarchical clustering works by grouping the data one by one on the basis of the nearest distance measure of all the pairwise distance between the data point. Distance could be defined differently and we have used euclidean distance and 1 - correlation to experiment with different distance measurements. Distance calculation methods also vary, which includes single linkage, complete linkage, average linkage, and ward's method (sum of squared euclidean distance is minimized), which is our choice of method in this case. The results from hierarchical clustering will not be displayed due to its relatively poor performance on the data.

From the Single Cell Comparison studies on Single Cell Portal, we acquired the true number of cell types for Cortex and PBMC, respectively.

- Cortex: 9 cell types
- PBMC: 9 cell types along with an unidentified cell type

Lastly, We used **silhouette** as a method of interpretation and validation of consistency within clusters of data. The silhouette value is a measure of how similar an object is to its own cluster (cohesion) compared to other clusters (separation). The silhouette ranges from -1 to +1, where a high value indicates that the object is well matched to its own cluster and poorly matched to neighboring clusters. Below are the visualizations of the clustering results.



Figure 10: Cortex 1 with PCA and t-SNE



Figure 11: Cortex 2 with PCA and t-SNE



Figure 12: PBMC 1 with PCA



Figure 13: PBMC 1 with t-SNE



Figure 14: PBMC 2 with PCA



Figure 15: PBMC 2 with t-SNE

Sample_Method	Silhouette(PCA)	Silhouette(TSNE)
Cortex1_10x-Chromium-v2	0.186910197	0.518623829
Cortex1_DroNc-seq	0.105032225	0.455349474
Cortex1_sci-RNA-seq	0.141394425	0.454163008
Cortex1_Smart-seq2	0.274517545	0.43116908
Cortex2_10x-Chromium-v2	0.205544537	0.494529206
Cortex2_DroNc-seq	0.175472146	0.478726973
Cortex2_sci-RNA-seq	0.158629506	0.388148937
Cortex2_Smart-seq2	0.216699042	0.532543572

Figure 16: Silhouette Scores for Cortex

Sample_Method	Silhouette(PCA)	Silhouette(TSNE)
PBMC1_10x-Chromium-v2-A	0.171863515	0.415434872
PBMC1_10x-Chromium-v2-B	0.161875009	0.45348634
PBMC1_10x-Chromium-v3	0.223967103	0.447270315
PBMC1_CEL-Seq2	0.156729312	0.37057283
PBMC1_Drop-seq	0.081447741	0.441737917
PBMC1_inDrops	0.077024669	0.411455969
PBMC1_Seq-Well	0.093972495	0.41759837
PBMC1_Smart-seq2	0.167225721	0.434217715
PBMC2_10x-Chromium-v2	0.201661206	0.451267706
PBMC2_CEL-Seq2	0.109379336	0.458700923
PBMC2_Drop-seq	0.102868241	0.381704566
PBMC2_inDrops	0.166264397	0.419217714
PBMC2_Seq-Well	0.045568644	0.37349809
PBMC2_Smart-seq2	0.145005358	0.386656882

Figure 17: Silhouette Scores for PHBMC

# References

- [1] Jiarui Ding, Xian Adiconis, Sean K. Simmons. Systematic comparative analysis of single cell RNA-sequencing methods. https://www.biorxiv.org/content/10.1101/632216v1.full.
- [2] Wei Vivian Li and Jingyi Jessica Li. Accurate and Robust Imputation for Single-cell RNA-seq data. http://jsb.ucla.edu/sites/default/files/Dahshu\_scImpute\_043018.pdf.
- [3] F. William Townes and Rafael A. Irizarryi. Quantile normalization of single-cell RNA-seq read counts without unique molecular. https://www.biorxiv.org/content/biorxiv/early/2019/10/24/817031.full.pdf.
- [4] Christoph Ziegenhain, Beate Vieth, Swati Parekh. Comparative Analysis of Single-Cell RNA Sequencing Methods. https://www.sciencedirect.com/science/article/pii/S1097276517300497
- [5] Study: Single Cell Comparison: PBMC data. https://singlecell.broadinstitute.org/single\_cell/study/SCP424/single-cell-comparison-pbmc-data-study-visualize.
- [6] Study: Single Cell Comparison: Cortex data. https://singlecell.broadinstitute.org/single\_cell/study/SCP425/singlecell-comparison-cortex-data-study-visualize
- [7] Soledad Espezuaa, Edwin Villanuevaa. Carlos D. Macielb, André Carvalhoa. A Projection Pursuit Algorithm for Exploratory Data Analysis. IEEE Transactions on Computers.
  J. H. Friedman and J. W. Tukey