XLMOD: Cross-linking and chromatography derivatization reagents ontology

Authors: Gerhard Mayer¹

1 Ruhr University Bochum, Faculty of Medicine, Medizinisches Proteom-Center, D-44801 Bochum, Germany

Abbreviations:

CID Collision Induced Dissociation

CV Controlled Vocabulary

CX-MS Cross-Linking-Mass Spectrometry

Da Dalton (atomic mass unit)
ESI ElectroSpray Ionization

HUPO-PSI HUman Proteome Organization – Proteomics Standards Initiative

IUPAC International Union of Pure and Applied Chemistry

InChI International Chemical Identifier

MS Mass Spectrometry

NMR Nuclear Magnetic Resonance

SMILES Simplified Molecular Input Line Entry Specification

XRC X-Ray Crystallography

Summary/Abstract:

Mass spectrometry has experienced a rapid development since its first application for protein analysis in the 1980s. While the most common use of mass spectrometry for protein analysis is identification and quantification workflows on peptides (digested from their parent protein), there is also a rapidly growing use of mass spectrometry for "structural proteomics". One example is the analysis of cross-linked proteins that can give valuable structural information, complementing the information gained by classical protein structure determination methods, useful for integrated methods of structure determination and modeling. For a broad and reproducible application of cross-linking mass spectrometry a standardized representation of cross-linking experimental results is necessary. This paper describes the developing and release of the xlmod ontology from the HUPO-PSI. xlmod contains terms for the description of reagents used in cross-linking experiments and of cross-linker related chemical modifications together with their main properties relevant for planning and performing cross-linking experiments. We also describe how xlmod is used within the new release of HUPO-PSI's mzldentML data standard, for reporting the used cross-linking reagents and results in a consistent manner. In addition xlmod contains terms for GC-MS and LC-MS derivatization reagents for specifying them in the upcoming mzTab-M and mzTab-L formats.

Keywords: Cross-linking reagents, cross-linker related chemical modifications, distance constraint, derivatization, xlmod ontology, mzldentML 1.2, mzTab

Download link: https://raw.githubusercontent.com/HUPO-PSI/mzldentML/master/cv/XLMOD.obo

1. Introduction

Traditionally proteins investigations are separated into two domains: At first, there are plenty of experimental methods for isolation, purification, identification and quantification of proteins (e.g. Difference Gel Electrophoresis (DIGE), SDS-PAGE (Sodium Dodecyl Sulfate – PolyacrylAmide Gel Electrophoresis), mass spectrometry (MS)). These aim at characterizing the proteins by their properties like e.g. their mass and sequence, their modifications, their physicochemical properties, interaction partners and so on. Second, there is the area of protein structural biology, aiming at the determination of tertiary and quaternary three-dimensional structure of the proteins, which is essential for their functional and dynamical behavior. With the technological progress of protein mass spectrometry (MS) in the last 3 decades MS is becoming more and more a useful tool for structural biology by giving information like dynamical, spatial and accessibility constraints, which can either be combined with other existing experimental protein structure determination or with computer modelling methods [1, 2]. Especially cross-linking mass spectrometry (CX-MS) [3, 4] is able to complement traditional methods like X-ray crystallography (XRC), nuclear magnetic resonance (NMR) and cryo-electron microscopy (cryo-EM), since each of these methods has some limitations: X-ray crystallography and NMR are both time-consuming and costly methods needing a relatively large amount of pure proteins. Furthermore they are either difficult to apply for membrane proteins and analyze the protein not in their native conformation (XRC), can be used only for smaller soluble proteins (NMR) or give only medium resolution information (cryo-EM), whereas CX-MS is a relatively fast high-throughput (HT) method by which also small amounts of proteins from protein mixtures can be analyzed.

The cross-linking reagents used in CX-MS contain reactive groups having specific reactivity's for reactions with the targeted functional groups of the amino acids and possessing spacer arms of varying lengths. Supplementary Table S5 shows some typical examples of such spacer arms. If such a cross-linking reagent reacts with residues from the same protein, one can derive distance constraints between these residues. These constraints are used as input for structural modeling. If the reagent links two different proteins, one can get information about the protein interaction partners and about their quaternary structure. Since the cross-links are covalent bonds, this allows also to detect transiently and weak binding interacting proteins. If the cross-linking reaction is terminated by a quencher, then the proteins contain also cross-linker specific, so-called dead end modifications. In general one distinguishes three types of cross-links [5]: Type 0 cross-links are cross-linker related dead end modifications, also called mono-links. They can give information about the solvent accessibility of the modified residues. Type 1 cross-links are intra-linked (intra-peptide or intra-protein) cross-links (also called loop-links). From intra-protein cross-links one can derive distance constraints for the maximum distance between the two cross-linked residues inside the protein. Type 2 cross-links are inter-linked (inter-peptide or inter-protein) cross-links. Inter-protein cross-links can give information about interacting protein partners and deliver distance constraints between subunits of protein complexes.

Since the efficiency of the cross-linking reactions is low, the cross-linked peptides occur in low abundance. Therefore, an enrichment step after the enzymatic digestion, e.g. a size exclusion chromatography (SCX), is

done before the cross-linked peptides are identified [6]. Another possibility is to use trifunctional cross-linkers, which contain a handle, e.g. biotin, which can later be used to separate the cross-linked peptides by affinity chromatography [7].

The identification of the cross-linked peptides by a search engine is challenging, because the identification of the cross-links has in general quadratic computational complexity compared with the linear complexity of a database search in normal MS [8]. This is due to the fact that one has to search for all possible binary combinations of peptides. This can either be done in an exhaustive search [9] or one explicitly encodes all the cross-linked peptides in the search database [10]. Therefore cross-linkers with certain properties, which can be utilized to develop specialized and faster search strategies, were developed. Examples are CID-cleavable cross-linkers [11, 12], isotope-labelled cross-linkers [13, 14], Protein Interaction Reporters (PIR) [15], photo-cleavable cross-linkers or cross-linker with MS2 labile bonds, allowing to identify the cross-linked peptides in the following MS3 step [16]. Another category is photoactivatable cross-linkers [17-19], which are mostly not selective regarding the targeted amino acid, but allow controlling the cross-linking reaction using a UV light source. Aldehydes such as glutardialdehyde or formaldehyde are broadly specific, i.e. they react with many amino acid residues and also with DNA and can also be used for the detection of protein-DNA interactions [20].

2. Cross-linking reagents

Besides deriving distance constraints for structure determination, cross-linking is useful to study receptor-ligand interactions, for the location determination of membrane proteins and cell surface receptors, for the immobilization on solid surfaces, for the production of immunotoxins and for labelling of proteins with fluorescent dyes (Table 1). For instance, the photo-affinity labeling method gets three-dimensional information about the binding interfaces, either between two proteins or between a protein and a drug in the binding pocket of that protein receptor for a small molecule drug [21]. This method uses cross-linkers with the photo-reactive groups (azido, benzophenone and diazirines) mentioned above. These cross-linkers are able to react with C-H bonds, which is the reason for the broad non-selective reactivity of these photoactivatable cross-linkers [21, 22]. In addition, there are cross-linkers available, which can link proteins with nucleic acids, lipids or carbohydrates. For instance cross-linkers with the photo-reactive psoralen group react with pyrimidines, especially the thymines of nucleic acids [23] and diazirine reactive groups are able to react with aldehydes of sugar chains [24].

Another possibility is the site-specific incorporation of cross-linkers into a cell's proteome by using photo-reactive amino acids analogues in order to study protein-protein interactions (PPIs) in living cells (*in vivo* crosslinking) [25]. One possibility is to genetically encode - by using an expanded genetic code - for the incorporation of a photo-reactive amino acid analogue like e.g. p-benzoyl-L-phenylalanine (pBpa) into the targeted protein [26-28]. For easy identification by doublets at M and M + 11, a deuterium-labelled pBpa-d11 version is used [29].

Other photo-reactive amino acids analogues that can be used are photo-L-Leu, photo-L-Ile and photo-L-Met [30] containing a photo-reactive diazirine group by which they can insert into C-H, N-H or S-H bonds, where no structural changes of the peptides due to the insertion of the diazirine ring were observed [31]. For all of them it was shown that their incorporation retains the protein structure and functions [28].

Table 1: Example use cases for cross-linking reagents usage.

Application	Comment
Deriving distance constraints for	Uses cross-linkers with different chemical reactivity and with
investigation of protein structure and	different spacer lengths to derive a set of distance constraints
guiding a structure-based modelling	which can restrict the search space in a subsequent structure
process.	modelling step. Often isotopically labelled heavy/light
	cross-linker pairs are used, since then the cross-linked
	peptides can be easily identified by the search software due to
	their characteristic doublet signals in the MS spectrum. Low
	concentration of protein and high cross-linker concentration
	favors intramolecular crosslinking. For the identification of
	subunits, cleavable cross-linkers are useful, so that by
	2D-gelelectrphoresis the subunits of the complex can be
	detected.
Study of protein-protein or	Uses protein cross-linkers to establish an association between
protein-DNA and receptor-ligand	two interacting proteins resp. a receptor protein and its ligand
interactions for the investigation of	for the derivation of the quaternary structure of protein
signal transduction.	complexes [32].
Determination of the location of	Integral membrane proteins react with water-insoluble
membrane proteins and	cross-linkers, but not with hydrophilic cross-linkers. Membrane
identification of cell surface	impermeable cross-linkers like Sulfo-NHS esters and
receptors.	pegylated cross-linkers, due to their hydrophilic nature, react
	only with proteins located on the cell surface.
Immobilization of proteins on solid	Typical supports are e.g. nitrocellulose membranes, glass,
supports.	agarose, polystyrene plates or beads. Allows e.g. to produce
	protein nanosurfaces for the production of biosensors or
	analytical devices like e.g. protein chips [33].
Use as immunological tools, e.g. for	Linking toxins to antibodies for getting molecules that can be
the production of immunotoxins or	used for targeted killing of cancer cells. Often cross-linkers
the coupling of haptens to carrier	with a cleavable disulfide S-S bond are used, because most
proteins for getting antigens for	cells are able to cleave these and thereby the bound toxin is
antibody production.	released.
Transfer of labels or reactive groups	Cross-linkers can also be used to modify target groups on a

	protein or to add spacers for later coupling reactions.	
Labeling of proteins.	Enable identification, enrichment, purification or detection of proteins by transferring labels like e.g. biotin or fluorescent	
	dyes from a known bait protein to an unknown interacting prey	
	protein.	
In vivo crosslinking [34].	Stabilizing transient endogenous protein-protein or	
	receptor-ligand interactions. Photo-reactive amino acids	
	analogues are used to incorporate the cross-linkers into the	
	cell's proteome.	

3. Properties and selection criteria for choosing cross-linking reagents

Criteria for the selection of a cross-linker were defined by Jin Lee [35] and Back et.al. [36]. The most important criteria are the specificity, which is given by the reactive groups, and the spacer lengths, especially with respect to getting distance constraints for later protein structure modeling.

Further possible criteria, e.g. the ability to permeate membranes, determine which kinds of proteins can be cross-linked. For an easier detectability of the cross-linked peptides by the used search machines, the availability of isotopically coded versions of the cross-linker or the cleavability of the cross-linker can be important, since then such cross-linked peptides are more easily detectable by the used search engines. When such cross-linkers are not available, one can make use of trifunctional cross-linkers with an additional handle, which allow for an additional enrichment step before the mass spectrometry analysis. Another possibility is the enrichment by using size exclusion chromatography [37]. Table 2 gives an overview about cross-linker properties, documented in the xlmod ontology and useful for cross-linker selection. In order to expand the applicability of the xlmod ontology, we included also some properties, which are not specific for XL-MS. Some vendors also offer cross-linker selection guides and/or interactive selection tools for choosing proper cross-linkers (see Supplementary Table S3).

Table 2: Cross-linker properties

Property	Comment / Reason / Example	
Chemical specificity	The specificity is determined by the reactive groups, targeting	
	specific groups of the different amino acids.	
Reaction conditions	e.g. buffer system, pH, concentrations,	
Spacer length	Distance between the cross-linked proteins acting as molecular	
	rulers for deriving spatial constraints for protein structure modelling	
	[38].	
Spacer composition	Determines hydrophobicity. For instance, a spacer containing PEG	
	(polyethylenglycol) makes the cross-linker hydrophobic and	
	therefore prevents hydrophobic aggregation.	

Spacer cleavability	If the cross-linker can be cleaved, either chemically, photo-activated
	or by fragmentation, e.g. CID (Collision Induced Dissociation). This
	cleavability e.g. allows the release of labeled cross-linked peptides
	from the supports used in an affinity based purification step [39].
Water solubility /	Determines if the cross-linker can detect cell surface proteins or can
membrane permeability	cross-link integral membrane proteins (which are hydrophobic) resp.
	if a cross-linker can permeate into cells (e.g. for in vivo cross-linking).
Isotope-coded versions	Makes the cross-link identification process by the search program easier, e.g. by searching for the characteristic doublets [29] of cross-links containing the heavy resp. light version of the cross-linker.
Selectivity of the reactive group	If spontaneously reactive with the targeted functional groups of residues or photo-reactive (non-selective), e.g. for two-step cross-linking reactions.
Handle for detection / identification or purification / enrichment	Using cross-linkers with a handle (e.g. biotin or a fluorophore) for easy enrichment of the cross-linked peptides.

4. Derivatization reagents

Because derivatization reagents like cross-linker have functional groups specifically reacting with chemical target groups, their description is very similar to that of the cross-linker reagents. Therefore, we decided to include them into the XLMOD ontology, so that beside the cross-linker reagents used in proteomics we also included derivatization reagents, useful for annotation in metabolomics and lipidomics. These derivatization reagents are either used in GC/MS for enhancing the volatility of polar metabolites, allowing them to pass easily into the gas phase [40, 41]. Similarly in LC-MS approaches the derivatization is done in order to enhance the LC signal intensities, the ionization in the ESI, (ElectroSpray Ionization) and/or the MS-detectability of the metabolite by a more efficient fragmentation [42]. The used derivatization reagents should specifically react with the targeted metabolites producing a high yield of the derivative [42]. The main groups of derivatization reagents for GC-MS are acylation, alkylation and silylation reagents, which replace active hydrogens (e.g. from -SH, -OH, -NH, -COOH groups) by an acyl, alkyl resp. silyl group, which confers the derivatized target molecule a more hydrophobic character, which again enhances their volatility.

Derivatization reagents for LC-MS at one hand have a functional group for the formation of the derivative, and on the other hand they introduce a group, which allows an easier detection of the derivative. This improved detectability is accomplished either in the LC part by introducing a chromophore, a fluorophore or an electrophore group for detection by UV/VIS light, by fluorescence or by electrochemical detectors [43]. Another possibility is to enhance the detectability in the MS part by improved ion formation in the ESI or by making the fragmentation more efficient. In mzTab-M [44] the used derivatization reagents should be

reported in the metadata section under derivatization_agent[1-n] as specified by the mzTab-M specification document available at https://github.com/HUPO-PSI/mzTab.

5. xlmod ontology for cross-linking and derivatization reagents

Although there are several articles listing cross-linkers for use in structural proteomics [45-47], there was still an urgent need for a structured ontology, which can serve both as a database for choosing proper cross-linkers, as well as source for uniquely describing identified cross-linked peptides and peptides with post-translational modifications derived from cross-linkers.

For building up the xlmod ontology we scanned vendor websites, publications and catalogues (see Supplementary Tables S1 - S3), relevant papers from the PubMed literature database [48] and some relevant books [49, 50].

The most important feature for the selection of a cross-linker reagent is its specificity for amino acids [51] and then it's spacer length. The amino acid residues, which are currently targetable by cross-linkers, are summarized in Supplementary Table S4.

It was shown, that the availability of cross-linking reagents targeting the acidic carboxyl groups of Asp and Glu allows deriving distance constraints that orthogonally complement constraints obtained from Lys cross-linking [51, 52]. It can be expected that cross-linkers specifically targeting other amino acids, e.g. the indole amine of tryptophan, can further enhance the power of the protein cross-linking method as a 3D structure determination tool.

The reactive groups occurring in cross-linkers and their amino acid specificities are reviewed in [25, 46, 47, 53], and summarized in Table 3.

Table 3: Overview about the reactive groups of cross-linking reagents

Reactive group	Target functional	Specificity for	Photo-re	Reference
	group	amino acids	active	(PMID)
formaldehyde	primarily amines,	K,N,R,Q,Y, N-term		25979347
	but broadly specific	secondary:		
		C,F,H,S,T,W		
1-hydroxy-7-azabenzotriaz	amines	K,N,Q,R,Protein		19994840
ole		N-term		
amine	carboxyl groups	D,E		
aryl halide	primary amines	K,N,Q,R,Protein		
		N-term		
aryl azides (e.g. phenyl-,	insertion into C-H	nonselective	yes	22641729
nitrophenyl-,	and N-H bonds;			

hydroxyphenyl-,	primary amines			
tetrafluorophenyl,-				
perfluoroaryl azide)				
azido-methylcoumarin		nonselective	yes	
benzophenone		nonselective	yes	
carbodiimide	carboxyl-to-amine	(K,N,Q,R,Protein		
	reactive	N-term)&(E,D,Prote		
		in C-term)		
carbonyl	hydrazines			
carboxyl	amines	K,N,Q,R,Protein		
		N-term		
diazirine		nonselective	yes	21778064
diazo group		nonselective	yes	
dimethylether	carboxyls	D,E		
glyoxal	arginine-specific	R		
haloacetyl (chloro-,	sulfhydryls	С		
bromo-, iodoacetyl)				
hydrazide	carbonyls	D,E,Protein C-term		
hydrazone	carbonyls	D,E,Protein C-term		
hydroxybenzotriazole	amines	K,N,Q,R,Protein		19994840
		N-term		
hydroxymethyl phosphine	amines	K,N,Q,R,Protein		
		N-term		
imidoester	amines	K,N,Q,R,Protein		
		N-term		
isocyanate	hydroxyl	S,T,Y		
maleimide	sulfhydryl			
methanethiosulfonate	sulfhydryl			
N-hydroxyphthalimide	amines	K,N,Q,R,Protein		19994840
		N-term		
NHS	lysines and	K,Protein N-term		
(N-hydroxysuccinimide)	N-termini, but also	secondary: S,T,Y		
ester	hydroxyls			
pentafluorophenyl	amines	K,N,Q,R,Protein		22067100
		N-term		
platinum (II)	thioether,	C,H,M		21591778
	sulfhydryl, and			

	imidazole		
psoralen	pyrimidines, esp.	Thy	
	thymins		
pyridinyldisulfide	sulfhydryls	С	
(pyridyldithiol)			
S-acetyl	incorporates		9177838
	sulfhydryls		
Sulfo NHS	lysines and	K,Protein N-term	
(N-hydroxysulfosuccinimid	N-termini, but also	secondary: S,T,Y	
e) ester	hydroxyls		
thioimidate	amines	K,N,Q,R,Protein	20795639
		N-term	
vinyl sulfone	amines, sulfhydryls,	C,K,N,Q,R,S,T,Y	
	hydroxyls		

It must be mentioned here that the "specificities" for amino acids are not absolute, but rather describe a tendency for preferred reaction targets, because the specificities can vary with the reaction conditions [54]. For instance NHS esters and Sulfo NHS esters react primarily with the amine group of lysines and with protein N-terminals, but to a lower degree they can also react with the hydroxyl groups of serine, threonine and tyrosine [25, 55], as specified by the property "secondarySpecificities".

In general, each cross-linking reagent term has a unique and fixed id, a short name, a definition, an optional list of synonyms and a list of properties with their values, the specification of its parent terms and then a list of relationships, describing e.g. the chemical reactivity.

As a concrete example, the entry for the reactive group "NHS ester" is defined as follows:

[Term]

id: xlmod:00101

name: NHS ester

def: "A reactive group (N-hydroxysuccinimide) that reacts with lysines and N-termini but also

serines, threonines and tyrosines." [PSI:XL]

synonym: "N-hydroxysuccinimide ester" EXACT []

property_value: specificities: "(K,Protein N-term)" xsd:string property_value: secondarySpecificities: "(S,T,Y)" xsd:string

is a: xlmod:00100 ! reactive group

relationship: has_property xlmod:00013 ! membrane permeable relationship: is_reactive_with xlmod:00028 ! primary amine reactive

6. General structure of the ximod ontology

The xlmod ontology consists of four main branches of cross-linker or derivatization reagent terms and defines several relations, which can be used to describe such terms in general (see Figure 1). The three main branches are:

- cross-linker related chemical modification
- cross-linking entity
- derivatization entity
- label transfer reagent

```
□ Classes
□ Proteomics Standards Initiative cross-linking and derivatization controlled vocabulary

☐ ← ② cross-linker related chemical modification

  ← III BDP-NHP-stump
  □ ← ② cross-linking entity
  □ ← ② derivatization entity
   ← II chromophore
   ← II electrophore
  □ Relations

    has_chromophore

 has_click_enabled_linker
 has_electrophore
 — has_fluorophore
 — has_handle
 — has_neutral_loss_reporter
 — has_property

    has_reactive_group

  — is_activatable
 is_analog_of
 is_cleavable
 — is_labelled
 — is_partially_reacted
 - is_reactive_with
 — is_side_product_of
____ part_of
```

The *cross-linker related chemical modification* branch contains partially reacted products of cross-linkers, e.g. modifications resulting from dead-ends (mono-links), which are no real cross-links, since they do not link two chains of proteins. In principle one can quench the cross-linking reaction with any amine, but the mostly ammonium salts or Tris are used. Therefore mainly amidated or hydrolyzed cross-linkers, but also thioesterized or alkenized forms are listed in xlmod. The amidated modifications are subdivided into ammonium amidated and Tris amidated forms, since these differ slightly in the resulting modifications and their masses, because ammonium salts convert to -C(=O)-NH2, whereas the use of Tris(hydroxymethyl)-aminomethane as a quenching substance leads to the formation of dead ends of the form -C(=O)-NHC(CH2OH)3. At the moment we have no instances of Tris amidated forms, but we distinguished these cases to allow them to be added later. The label transfer reagents are either trifunctional cross-linkers containing a handle, e.g. biotin, for enrichment or for transferring a label to a prey protein.

The main branches under cross-linking entity are:

- amino acid and their derivatives (containing photo-reactive amino acid analogs)
- the cross-linkers in the narrower sense (resp. the cross-linked products derived from them)
- cross-linking attributes for a more detailed characterization of the cross-linkers, e.g. their cleavability, and cross-linker derived chemical modifications
- handles used for detection, affinity enrichment or purification of cross-linked peptides (see Supplementary Table S6)
- reactive group attributes for a more detailed description of reactive groups

Important cross-linking attributes are the reactive groups, which have specific chemical reactivity towards functional target groups, and which are listed under the reactive group attributes.

The branch derivatization entity contains sub-branches for

- chromophores (for UV/VIS detection)
- fluorophores (for fluorescence detection)
- electrophores (for electrochemical detection)
- GCMS derivatization reagents
- LCMS derivatization reagents.

A typical example for a cross-linker term is:

```
[Term]
id: xlmod:02001
```

name: DSS

def: "Disuccinimidyl suberate." [CAS:68528-80-3, PubChem_Compound:100658, PMID:16944939,

ChemSpiderID:90944, ChemicalBookNo:CB1219157]

synonym: "Bis(succinimidyl) suberate" EXACT []

synonym: "Suberic acid bis(N-hydroxysuccinimide ester)" EXACT []

synonym: "Disuccinimidyl octanedioate" EXACT []

synonym: "DSS-d0" EXACT [] synonym: "DSS-H12" EXACT []

synonym: "1,1'-[(1,8-Dioxooctane-1,8-diyl)bis(oxy)]dipyrrolidine-2,5-dione" EXACT []

property_value: bridgeFormula: "C8 H10 O2" xsd:string

property_value: monolsotopicMass: "138.06807961" xsd:double

property_value: reactionSites: "2" xsd:nonNegativeInteger

property_value: spacerLength: "11.4" xsd:float is_a: xlmod:00005! homofunctional cross-linker

relationship: has_reactive_group xlmod:00101 ! NHS ester

The cross-linker entries possess a unique ID and (short) name. Then follows a definition line containing the long name as well as a list of database cross references, which can be of one of the types listed in Table 4. The PMID cross references and web links can occur several times for a given definition, whereas the other possible cross references must be unique.

Table 4: Used database cross references (dbxref's) in the xlmod.obo ontology

DB cross reference	DB name	Link
Beilstein	Beilstein Registry Number	https://www.elsevier.com/solutions/reaxys
CAS	Chemical Abstract Service	https://www.cas.org/content/chemical-substances
CBNumber	Chemical Book	http://www.chemicalbook.com/Search_EN.aspx
ChemSpiderID	ChemSpider	http://www.chemspider.com
MDL	former Molecular Design	http://accelrys.com
	Ltd., now BIOVIA	https://www.3ds.com/products-services/biovia/
PMID	Pubmed ID	http://www.ncbi.nlm.nih.gov/pubmed/advanced
PubChem_Compound	PubChem	https://www.ncbi.nlm.nih.gov/pccompound/advanced
http://	Web link	specific web link, where the cross-linker is defined
DOI:	Web link	web link, e.g. for publications not indexed in PubMed

The chemical substance registry references can be used, if one wants to get more detailed information about a cross-linker, e.g. the IUPAC name, a full list of synonyms, the SMILES string, InChI key, 3D formula, and so on. From the Chemical Book repository, one can get also MDL .mol files containing the exact coordinates of the 3D structure. In principle one can then calculate the spacer arm lengths from the coordinates in these mol files [56], if one is acquainted with the reaction mechanism, so that one exactly knows which atoms belong to the spacer chain of the resulting cross-linked peptides. Details of that format are described at http://accelrys.com/products/collaborative-science/biovia-draw/ctfile-no-fee.html (accessed 08/2018) and http://media.accelrys.com/downloads/ctfile-formats/ctfile-formats.zip. Other possibilities to get the spacer arm length are to calculate them from force filed simulations [38] or to calculate

them from 3D coordinates. If the coordinates are unknown, one can calculate them from the SMILES strings, which one can get from PubChem, by using cheminformatics programs like e.g. the 3D structure generator CORINA, available from http://www2.chemie.uni-erlangen.de/services/telespec/corina/ (accessed 08/2018) [57], which converts the SMILES string into 3D coordinates. Another possibility is to measure the spacer arm lengths with chemical drawing programs like e.g. BIOVIA Draw (http://accelrys.com/resource-center/downloads/freeware/, accessed 08/2018). A tool that can be used to calculate the m/z ratio of cross-linked peptides is the computer program GPMAW [58] (http://www.gpmaw.com/html/cross-linking.html and http://www.gpmaw.com/html/cross-linker.html,

accessed 08/2018).

Further available software tools for cross-linking analysis and identification are already summarized by Mayne and Patterton [59], Sinz [60], Leitner *et.al.* [16], Tran *et.al.* [32], and in the list of cross-linking software tools at https://omictools.com/cross-linked-peptide-identification-category (accessed 08/2018).

In the xlmod ontology we only listed synonyms, which are commonly used by the cross-linking community. The bridge formula (sum formula of the spacer arm) is given in Hill notation [61]. In addition, the monoisotopic mass calculated from the bridge formula, is specified. Then the number of reaction sides is given and the cross-linker is either classified as homofunctional (i.e. the cross-linker has only one sort of reactive group) or heterofunctional (i.e. containing more than one different reactive groups). In the DSS example above the number of reaction sites is 2, so that it's a homobifunctional cross-linker. At the end the type of reactive groups are listed.

In order to avoid redundancy, the specificities for the amino acids are not repeated in the term for the cross-linker, so that the specificities must be looked up at the terms for the listed reactive groups of that cross-linker. This means that the properties from reactive groups are inherited and apply to the cross-linking reagents containing the respective reactive group. For example, it follows that a cross-linker with a NHS ester reactive tends to be more membrane permeable, whereas one with a Sulfo-NHS ester reactive group is more hydrophilic. The overall hydrophobicity is also influenced by the hydrophilicity of the spacers and affinity arms, e.g. PEG spacers make a reagent the more hydrophilic the longer the PEG spacer group is. Likewise, it is clear, that a cross-linker with a pyridinyldisulfide reactive group possesses a cleavable S-S bond or that a cross-linker with a biotin handle is enrichable.

An example for a cross-linker related post-translational modification is the term for "hydrolyzed DSS", where the relationship "is_side-product_of" specifies from what cross-linker reagent the chemical modification is derived from:

[Term]

id: xlmod:01002

name: hydrolyzed DSS

def: "Hydrolyzed disuccinimidyl suberate." [PSI:XL]

property_value: deadEndFormula: "C8 H12 O3" xsd:string

property_value: monolsotopicMass: "156.07864431" xsd:double

is a: xlmod:00096! hydrolization cross-linker related chemical modification

relationship: is_partially_reacted xlmod:00011 ! hydrolyzed

relationship: is_side_product_of xlmod:02001 ! DSS

The following is an example for a labelled cross-linker with the relationship "is_labelled":

[Term]

id: xlmod:02002 name: DSS-d4

def: "Deuterium labelled disuccinimidyl 2,2,7,7-suberate." [PubChem_Compound:91757798,

PMID:11354472]

property_value: bridgeFormula: "C8 D4 H6 O2" xsd:string property_value: doubletDeltaMass: "4.02508" xsd:double

property_value: monolsotopicMass: "142.093186586" xsd:double

property_value: reactionSites: "2" xsd:nonNegativeInteger

property_value: spacerLength: "11.4" xsd:float is_a: xlmod:00005! homofunctional cross-linker

relationship: has_reactive_group xlmod:00101 ! NHS ester relationship: is_labelled xlmod:00010 ! deuterium labelled

The property "doubletDeltaMass" gives the m/z difference between the labelled and the unlabeled versions of the cross-linker and can be used by cross-linking search-engines for easy detection of the corresponding cross-linked peptides.

As an example for an entry of a photo-reactive amino acid analogue the entry for pBpa-d11 is specified as:

[Term]

id: xlmod:01904 name: pBpa-d11

def: "Deuterium labelled p-benzoyl-L-phenylalanine." [PMID:18704231]

comment: A photoreactive deuterium-labelled amino acid analog of L-Phenylalanine for incorporation during protein synthesis that can be used for in vivo labelling, cross-linking and protein-protein interaction studies in live cells.

property_value: reactionSites: "1" xsd:nonNegativeInteger

is a: xlmod:00049 ! amino acid derivative

relationship: has_reactive_group xlmod:00119! benzophenone

relationship: is_activatable xlmod:00149! photoactivatable relationship: is_analog_of xlmod:00070! L-Phenylalanine relationship: is_labelled xlmod:00010! deuterium labelled

Since the PSI-MOD ontology is not maintained anymore, we also added a special branch to the xlmod ontology that contains label transfer reagents, which are primarily used for transferring labels and not for cross-linking. An example would be

[Term]

id: xlmod:01800

name: Psoralen-PEG3-Biotin

def: "Psoralen-PEG3-Biotin biotinylation reagent." [PSI:XL]

property_value: hydrophilicPEGchain: "3" xsd:nonNegativeInteger

property_value: reactionSites: "1" xsd:nonNegativeInteger

property_value: spacerLength: "36.86" xsd:float

is_a: xlmod:00003 ! label transfer reagent relationship: has handle xlmod:00051 ! biotin

relationship: has_property xlmod:00014! hydrophilic relationship: has_reactive_group xlmod:00134! psoralen

At the moment the xlmod ontology contains 16 type definitions and 1083 terms, of which 248 are bifunctional cross-linker terms, 103 terms belong to the class of cross-linker related chemical modifications and 84 terms describe chemical reactivities, reactive groups and reactive group attributes. In addition there are 174 terms for GC/MS derivatization reagents and 324 terms for LC/MS derivatization reagents. The other terms describe zero-length and trifunctional cross-linkers, label transfer reagents, amino acids and their derivatives, handles, chemical reactivities and other attributes and properties of the cross-linkers. For details see the Supplementary Table S7 and S8.

7. Reporting of cross-linking results

Controlled vocabulary (CV) terms from the xlmod ontology and from the psi-ms ontology [62] are used together for documenting cross-linking searches in version 1.2 of the XML-based mzldentML standard format [63], see Figure 2. For that one has to specify the CV term "crosslinking search" (MS:1002494) in the <SpectrumIdentificationProtocol> element of the mzldentML file [63]. The two cross-linked peptides must have <Modification> elements, which are marked as cross-linked peptides by using the CV terms "cross-link donor" and "cross-link acceptor", which are defined in the psi-ms ontology [62]. These two corresponding peptides are matched together by having the same value. The mzldentML validator [64] checks that these

two CV terms, grouped together by the same value, are located under different <Modification> elements. By convention, the donor peptide is the one with the longer peptide sequence. In case that both peptides have the same length, the donor is defined as the one, which comes first in the alphabetical order of the peptide sequence.

```
<Peptide·id="PEP_1">¶
        <PeptideSequence>STMLEKIK</PeptideSequence>¶
                                                                                                           XL·Donor·specifies·the·cross-
        ...¶
        <Modifi
                ation·location="6"·residues="K"·monoisotopicMassDelta="138.06807961">¶
                                                                                                           linking-reagent¶
                 cvParam-accession="xlmod:02001"-cvRef="PSI-MS"-name="DSS"/>¶
                   vParam-accession="MS:1002509"· cvRef="PSI-MS"· name="cross-link·donor"·value="CLM_1"/>¶
        </Modificat
</Peptide>¶
<Peptide·id="PEP_2">¶
        <PeptideSequence>LCHKNK</PeptideSequence>¶
                                                                                                           XL-Acceptor-reported-with-
               ation·location="5"·residues="K"·monoisotopicMassDelta="0">¶
        <Modifi
                                                                                                           monoisotopicMassDelta of 0 9
                  cvParam·accession="MS:1002510"·c<del>vRef</del>="PSI-MS"·name="cross-link·acceptor"·value="CLM_1"/>¶
        </Modifica
</Peptide>¶
<SpectrumIdentificationPro col·id="SIP_1"-analysisSoftware_ref="AS_1">¶
        <AdditionalSearchParams
                                    "MS:1002494"· cvRef="PSI-MS"· name="cros
                <cvParam·accession=
                                                                           -linking·search"/>¶
        ...¶
</SpectrumIdentificationProtocol->¶
<SpectrumIdentificationResult-spectraData_ref="SDATA_1"·spectrumID=\(\scan=1\),scan=\(\frac{1}{2}\)\\.id="SIR 1">\(\frac{1}{2}\)\\.
        <SpectrumIdentificationItem·passThreshold="1"·rank="1"·peptide</p>
                                                                   ref="PEP_1<sup>1</sup>"-experimentalMassToCharge="672.374450"-chargeState="3"-id="SII_1">¶
                 <cvParam·accession="MS:1002511"· cvRef="PSI-MS"· name
                                                                      cross-link-spectrum-identification-item"-value="CL_SII_1"/>"
                <cvParam·accession="MS:1002545"·cvRef="PSI-MS"·name=\vi:score"·value="1.23"/>¶
        </SpectrumIdentificationItem>¶
        ...¶
                 <cvParam·accession="MS:1002511"·cvRef="PSI-MS"·name="cross-link·spectrum·identification·item"·value="CL_SII_1"/>¶
                <cvParam·accession="MS:1002545"·cvRef="PSI-MS"·name="xi:score"·value="1.23"/>¶
        </SpectrumIdentificationItem>¶
</SpectrumIdentificationResult>¶
```

By convention, the cross-link acceptors are reported by a <Modification> element with a monoisotopic mass delta of zero, and the cross-link donor with the actual mass of the cross-link. In addition the peptide, which is the "cross-link donor" must have another CV term annotating its <Modification> element [63]. This term must be a term located under the branch "cross-linker" (xlmod:00004) from the xlmod ontology, and describes the cross-linking reagent. In contrast, a dead-end modification is annotated by only one CV term, referencing a child term of "cross-linker related chemical modification" (xlmod:00002) from the xlmod ontology.

Since there are always pairs of donor and receiver peptides, there must be two <PeptideEvidenceRef> elements (referencing the cross-linked peptides), that occur paired with the same charge state in two <SpectrumIdentificationItem> elements under the same <SpectrumIdentificationResult> element with the same value for the CV term "cross-link spectrum identification item" (MS:1002511) annotating these <SpectrumIdentificationItem> elements. If one uses isotope-labelled cross-linker reagents, there are even four <SpectrumIdentificationItem> elements for the same cross-link and charge state, since donor and acceptor can occur now cross-linked with the heavy and light version of the used cross-linker reagent. This rule is also checked by the mzldentML validator [64].

We want to mention here, that beside mzldentML 1.2 there is also an alternative open data format for reporting cross-linking results, developed by Hoopmann et.al. [65] and is an extension of the pepXML format [66].

8. Outlook to the further development of xlmod

Then HUPO-PSI intends to steadily further develop the xlmod ontology. Currently, on the one hand there are still some cross-linkers missing and on the other hand for some cross-linkers there are still some data missing. Our hope is that the cross-linking community will contribute such missing data and request terms for still missing respective newly developed cross-linkers via the 'psidev-ms-vocab' mailing list at https://lists.sourceforge.net/lists/listinfo/psidev-ms-vocab and thereby contribute to make this ontology a standard source for cross-linkers, their related chemical modifications, and derivatization reagents. Then this ontology can be utilized by specialized cross-linking identification and/or analysis software like e.g. xiFDR [67], SIM-XL [68] or OpenMS [69], which allow to report their results in the mzldentML 1.2 [63] standard format. The xlmod ontology is also open for bio-conjugative reagents used in cross-linking of either RNA-RNA [70], or proteins to other biomolecules like e.g. DNA / RNA [23, 71-73], lipids [74, 75] or sugars [76-78]. For instance there are lipid analogues defined, which are useful for describing the interaction of lipids with membrane proteins [75], the terms for nucleobases and their derivatives can be used for the description of protein-RNA interactions [72] and the sialic acid analogue terms are intended for the description of glycan-protein crosslinking [78] experiments.

Conceivable is the use of xlmod as a database for the development of an open-source tool for the planning of cross-linking experiments, allowing the selection of the most promising cross-linking reagents dependent of the detailed aims of the planned experiment. For the easy access to xlmod the development of a general Java access library, resolving all the inheritance dependencies between the different terms, e.g. between the reactive groups and the cross-linking reagents, would also be desirable.

Desirable for the reproducible application of CX-MS in integrated structure determination and prediction approaches are also standardized experimental protocols, since the specificities of the cross-linking reagents depends on the reaction conditions [79], as well as robust, standardized and automatable data analysis workflows [80, 81] for such integrated structure determination and modeling methods. The xlmod ontology is one building block for such a future highly automated high-throughput application of CX-MS in structural biology.

Acknowledgements:

Gerhard Mayer was funded by the BMBF grant de.NBI - German Network for Bioinformatics Infrastructure (FKZ 031 A 534A).

References:

- 1. Benesch JL, Ruotolo BT. Mass spectrometry: come of age for structural and dynamical biology, Curr Opin Struct Biol 2011;21:641-649.
- 2. Fabris D, Yu ET. Elucidating the higher-order structure of biopolymers by structural probing and mass spectrometry: MS3D, Journal of Mass Spectrometry 2010;45:841-860.
- 3. Holding AN. XL-MS: Protein cross-linking coupled with mass spectrometry, Methods 2015;89:54-63.
- 4. Leitner A, Faini M, Stengel F et al. Crosslinking and Mass Spectrometry: An Integrated Technology to Understand the Structure and Function of Molecular Machines, Trends in Biochemical Sciences 2016;41:20-32.
- 5. Schilling B, Row RH, Gibson BW et al. MS2Assign, automated assignment and nomenclature of tandem mass spectra of chemically crosslinked peptides, J Am Soc Mass Spectrom 2003;14:834-850.
- 6. Fritzsche R, Ihling CH, Gotze M et al. Optimizing the enrichment of cross-linked products for mass spectrometric protein analysis, Rapid Commun Mass Spectrom 2012;26:653-658.
- 7. Sinz A, Kalkhof S, Ihling C. Mapping protein interfaces by a trifunctional cross-linker combined with MALDI-TOF and ESI-FTICR mass spectrometry, J Am Soc Mass Spectrom 2005;16:1921-1931.
- 8. Yang L, Tang X, Weisbrod CR et al. A photocleavable and mass spectrometry identifiable cross-linker for protein interaction studies, Analytical Chemistry 2010;82:3556-3566.
- 9. Yu FC, Li N, Yu WC. ECL: an exhaustive search tool for the identification of cross-linked peptides using whole database, BMC Bioinformatics 2016;17.
- 10. Yilmaz S, Drepper F, Hulstaert N et al. Xilmass: A New Approach toward the Identification of Cross-Linked Peptides, Analytical Chemistry 2016;88:9949-9957.
- 11. Gotze M, Pettelkau J, Fritzsche R et al. Automated assignment of MS/MS cleavable cross-links in protein 3D-structure analysis, J Am Soc Mass Spectrom 2015;26:83-97.
- 12. Sinz A. Divide and conquer: cleavable cross-linkers to study protein conformation and protein-protein interactions, Anal Bioanal Chem 2016.
- 13. Petrotchenko EV, Borchers CH. ICC-CLASS: isotopically-coded cleavable crosslinking analysis software suite, BMC Bioinformatics 2010;11:64.
- 14. Rinner O, Seebacher J, Walzthoeni T et al. Identification of cross-linked peptides from large sequence databases, Nat Methods 2008;5:315-318.
- 15. Tang X, Bruce JE. A new cross-linking strategy: protein interaction reporter (PIR) technology for protein-protein interaction studies, Mol Biosyst 2010;6:939-947.
- 16. Leitner A, Walzthoeni T, Kahraman A et al. Probing native protein structures by chemical cross-linking, mass spectrometry, and bioinformatics, Mol Cell Proteomics 2010;9:1634-1649.
- 17. Preston GW, Wilson AJ. Photo-induced covalent cross-linking for the analysis of biomolecular interactions, Chem Soc Rev 2013;42:3289-3301.
- 18. Dubinsky L, Krom BP, Meijler MM. Diazirine based photoaffinity labeling, Bioorg Med Chem 2012;20:554-570.
- 19. Das J. Aliphatic Diazirines as Photoaffinity Probes for Proteins: Recent Developments, Chem Rev 2011;111:4405-4417.
- 20. Lu K, Ye WJ, Zhou L et al. Structural Characterization of Formaldehyde-Induced Cross-Links Between Amino Acids and Deoxynucleosides and Their Oligomers, J Am Chem Soc 2010;132:3388-3399.
- 21. Robinette D, Neamati N, Tomer KB et al. Photoaffinity labeling combined with mass spectrometric approaches as a tool for structural proteomics, Expert Rev Proteomics 2006;3:399-408.
- 22. Hatanaka Y. Development and leading-edge application of innovative photoaffinity labeling, Chem Pharm Bull (Tokyo) 2015;63:1-12.
- 23. Sastry SS, Ross BM, P'Arraga A. Cross-linking of DNA-binding proteins to DNA with psoralen and psoralen furan-side monoadducts. Comparison of action spectra with DNA-DNA cross-linking, J Biol Chem 1997;272:3715-3723.

- 24. Tomohiro T, Hashimoto M, Hatanaka Y. Cross-linking chemistry and biology: development of multifunctional photoaffinity probes, Chem Rec 2005;5:385-395.
- 25. Sinz A. Chemical cross-linking and mass spectrometry to map three-dimensional protein structures and protein-protein interactions, Mass Spectrom Rev 2006;25:663-682.
- 26. Hino N, Okazaki Y, Kobayashi T et al. Protein photo-cross-linking in mammalian cells by site-specific incorporation of a photoreactive amino acid, Nat Methods 2005;2:201-206.
- 27. Farrell IS, Toroney R, Hazen JL et al. Photo-cross-linking interacting proteins with a genetically encoded benzophenone, Nat Methods 2005;2:377-384.
- 28. Lee HS, Dimla RD, Schultz PG. Protein-DNA photo-crosslinking with a genetically encoded benzophenone-containing amino acid, Bioorganic & Medicinal Chemistry Letters 2009;19:5222-5224.
- 29. Wilkins BJ, Daggett KA, Cropp TA. Peptide mass fingerprinting using isotopically encoded photo-crosslinking amino acids, Molecular Biosystems 2008;4:934-936.
- 30. Lossl P, Sinz A. Combining Amine-Reactive Cross-Linkers and Photo-Reactive Amino Acids for 3D-Structure Analysis of Proteins and Protein Complexes, Methods Mol Biol 2016;1394:109-127.
- 31. Kolbel K, Ihling CH, Sinz A. Analysis of peptide secondary structures by photoactivatable amino acid analogues, Angew Chem Int Ed Engl 2012;51:12602-12605.
- 32. Tran BQ, Goodlett DR, Goo YA. Advances in protein complex analysis by chemical cross-linking coupled with mass spectrometry (CXMS) and bioinformatics, Biochim Biophys Acta 2016;1864:123-129.
- 33. Bhakta SA, Evans E, Benavidez TE et al. Protein adsorption onto nanomaterials for the development of biosensors and analytical devices: a review, Anal Chim Acta 2015;872:7-25.
- 34. Kaake RM, Wang X, Burke A et al. A new in vivo cross-linking mass spectrometry platform to define protein-protein interactions in living cells, Mol Cell Proteomics 2014;13:3533-3543.
- 35. Jin Lee Y. Mass spectrometric analysis of cross-linking sites for the structure of proteins and protein complexes, Mol Biosyst 2008;4:816-823.
- 36. Back JW, de Jong L, Muijsers AO et al. Chemical cross-linking and mass spectrometry for protein structural modeling, J Mol Biol 2003;331:303-313.
- 37. Leitner A, Reischl R, Walzthoeni T et al. Expanding the Chemical Cross-Linking Toolbox by the Use of Multiple Proteases and Enrichment by Size Exclusion Chromatography, Molecular & Cellular Proteomics 2012;11.
- 38. Green NS, Reisler E, Houk KN. Quantitative evaluation of the lengths of homobifunctional protein cross-linking reagents used as molecular rulers, Protein Science 2001;10:1293-1304.
- 39. Rudolf GC, Heydenreuter W, Sieber SA. Chemical proteomics: ligation and cleavage of protein modifications, Curr Opin Chem Biol 2013;17:110-117.
- 40. Danielson ND, Gallagher PA, Bao JJ. Chemical Reagents and Derivatization Procedures in Drug Analysis. Chichester: John Wiley & Sons Ltd, 2000.
- 41. Orata F. Derivatization Reactions and Reagents for Gas Chromatogryphy Analysis. IntechOpen, 2011.
- 42. Santa T. Derivatization in liquid chromatography for mass spectrometric detection, Drug Discoveries & Therapeutics 2013;7:9-17.
- 43. Poole CF. Derivatization in Liquid Chromatography. Elsevier, 2017.
- 44. Hoffmann N, Rein J, Sachsenberg T et al. mzTab-M: A Data Standard for Sharing Quantitative Results in Mass Spectrometry Metabolomics, Anal Chem 2019;91:3302-3310.
- 45. Mattson G, Conklin E, Desai S et al. A practical approach to crosslinking, Mol Biol Rep 1993;17:167-183.
- 46. Serpa JJ, Parker CE, Petrotchenko EV et al. Mass spectrometry-based structural proteomics, Eur J Mass Spectrom (Chichester, Eng) 2012;18:251-267.
- 47. Paramelle D, Miralles G, Subra G et al. Chemical cross-linkers for protein structure studies by mass spectrometry, Proteomics 2013;13:438-456.

- 48. Lu Z. PubMed and beyond: a survey of web tools for searching biomedical literature, Database (Oxford) 2011;2011:baq036.
- 49. Hermanson GT. Bioconjugate Techniques. Oxford: Elsevier, 2013.
- 50. Wong SS, Jameson, D.M. Chemistry of Protein and Nucleic Acid Cross-Linking and Conjugation. CRC Press Inc, 2011.
- 51. Hofmann T, Fischer AW, Meiler J et al. Protein structure prediction guided by crosslinking restraints A systematic evaluation of the impact of the crosslinking spacer length, Methods 2015;89:79-90.
- 52. Leitner A, Joachimiak LA, Unverdorben P et al. Chemical cross-linking/mass spectrometry targeting acidic residues in proteins and protein complexes, Proc Natl Acad Sci U S A 2014;111:9455-9460.
- 53. Zybailov BL, Glazko GV, Jaiswal M et al. Large Scale Chemical Cross-linking Mass Spectrometry Perspectives, J Proteomics Bioinform 2013;6:001.
- 54. Kalkhof S, Sinz A. Chances and pitfalls of chemical cross-linking with amine-reactive N-hydroxysuccinimide esters, Anal Bioanal Chem 2008;392:305-312.
- 55. Madler S, Bich C, Touboul D et al. Chemical cross-linking with NHS esters: a systematic study on amino acid reactivities, J Mass Spectrom 2009;44:694-706.
- 56. Dalby A, Nourse JG, Hounshell WD et al. Description of Several Chemical-Structure File Formats Used by Computer-Programs Developed at Molecular Design Limited, Journal of Chemical Information and Computer Sciences 1992;32:244-255.
- 57. Sadowski J, Gasteiger J, Klebe G. Comparison of Automatic Three-Dimensional Model Builders Using 639 X-ray Structures, Journal of Chemical Information and Modeling 1994;34:1000-1008.
- 58. Peri S, Steen H, Pandey A. GPMAW--a software tool for analyzing proteins and peptides, Trends in Biochemical Sciences 2001;26:687-689.
- 59. Mayne SL, Patterton HG. Bioinformatics tools for the structural elucidation of multi-subunit protein complexes by mass spectrometric analysis of protein-protein cross-links, Brief Bioinform 2011;12:660-671.
- 60. Sinz A, Arlt C, Chorev D et al. Chemical cross-linking and native mass spectrometry: A fruitful combination for structural biology, Protein Sci 2015;24:1193-1209.
- 61. Hill EA. ON A SYSTEM OF INDEXING CHEMICAL LITERATURE; ADOPTED BY THE CLASSIFICATION DIVISION OF THE U. S. PATENT OFFICE. J. Am. Chem. Soc., 1900, 478–494.
- 62. Mayer G, Montecchi-Palazzi L, Ovelleiro D et al. The HUPO proteomics standards initiative- mass spectrometry controlled vocabulary, Database (Oxford) 2013;2013:bat009.
- 63. Vizcaíno JA, Mayer G, Perkins S et al. The mzldentML Data Standard Version 1.2, Supporting Advances in Proteome Informatics, Mol Cell Proteomics 2017;16:1275-1285.
- 64. Ghali F, Krishna R, Lukasse P et al. Tools (Viewer, Library and Validator) that facilitate use of the peptide and protein identification standard format, termed mzldentML, Mol Cell Proteomics 2013;12:3026-3035.
- 65. Hoopmann MR, Mendoza L, Deutsch EW et al. An Open Data Format for Visualization and Analysis of Cross-Linked Mass Spectrometry Results, J Am Soc Mass Spectrom 2016.
- 66. Keller A, Eng J, Zhang N et al. A uniform proteomics MS/MS analysis platform utilizing open XML file formats, Mol Syst Biol 2005;1:2005 0017.
- 67. Fischer L, Rappsilber J. Quirks of Error Estimation in Cross-Linking/Mass Spectrometry, Anal Chem 2017;89:3829-3833.
- 68. Lima DB, de Lima TB, Balbuena TS et al. SIM-XL: A powerful and user-friendly tool for peptide cross-linking analysis, J Proteomics 2015;129:51-55.
- 69. Rost HL, Sachsenberg T, Aiche S et al. OpenMS: a flexible open-source software platform for mass spectrometry data analysis, Nat Methods 2016;13:741-748.
- 70. ME H. RNA Crosslinking Methods. In: ML (ed) Methods Enzymol.: Springer, 2009, 127-146.

- 71. Park AY, Robinson CV. Protein-nucleic acid complexes and the role of mass spectrometry in their structure determination, Critical Reviews in Biochemistry and Molecular Biology 2011;46:152-164.
- 72. Schmidt C, Kramer K, Urlaub H. Investigation of protein-RNA interactions by mass spectrometry--Techniques and applications, J Proteomics 2012;75:3478-3494.
- 73. Aoki T, Wolle D, Preger-Ben Noon E et al. Bi-functional cross-linking reagents efficiently capture protein-DNA complexes in Drosophila embryos, Fly (Austin) 2014;8:43-51.
- 74. Saliba AE, Vonkova I, Gavin AC. The systematic analysis of protein-lipid interactions comes of age, Nat Rev Mol Cell Biol 2015;16:753-761.
- 75. Gubbens J, Ruijter E, de Fays LE et al. Photocrosslinking and click chemistry enable the specific detection of proteins interacting with phospholipids at the membrane interface, Chem Biol 2009;16:3-14.
- 76. Yu SH, Bond MR, Whitman CM et al. Metabolic labeling of glycoconjugates with photocrosslinking sugars, Methods Enzymol 2010;478:541-562.
- 77. McCombs JE, Zou C, Parker RB et al. Enhanced Cross-Linking of Diazirine-Modified Sialylated Glycoproteins Enabled through Profiling of Sialidase Specificities, ACS Chem Biol 2016;11:185-192.
- 78. Bond MR, Zhang H, Vu PD et al. Photocrosslinking of glycoconjugates using metabolically incorporated diazirine-containing sugars, Nat Protoc 2009;4:1044-1063.
- 79. Toews J, Rogalski JC, Clark TJ et al. Mass spectrometric identification of formaldehyde-induced peptide modifications under in vivo protein cross-linking conditions, Anal Chim Acta 2008;618:168-183.
- 80. Arlt C, Gotze M, Ihling CH et al. Integrated Workflow for Structural Proteomics Studies Based on Cross-Linking/Mass Spectrometry with an MS/MS Cleavable Cross-Linker, Anal Chem 2016;88:7930-7937.
- 81. Mendes ML, Fischer L, Chen ZA et al. An integrated workflow for cross-linking/mass spectrometry, bioRxiv 2018:24.

Supplementary Table S1: Used vendor web page links.

Vendor	Web site (accessed 02/2020)	
Santa Cruz	http://www.scbt.com/chemicals-table-crosslinkers.html	
Biotechnology (SCB)		
SCB	http://www.scbio.de/chemicals-table-crosslinkers.html	
ChemCruz	http://www.scbio.de/chemicals-table-homobifunctional_crosslinkers.html	
Invitrogen	https://www.thermofisher.com/search/global/searchAction.action?query=crosslinkers	
Sigma Aldrich	http://www.sigmaaldrich.com/catalog/search?interface=All&term=crosslinker&N=0&page=2&mode=mode+matchpartialmax&focus=product	
Sigma Aldrich	http://www.sigmaaldrich.com/life-science/molecular-biology/molecular-biology-products.html?TablePage=9621740	
Clear Synth	http://www.clearsynth.ca/search.asp?category=Cross%20Linking%20Reagents	
ThermoFisher	https://www.thermofisher.com/search/results?query=crosslinkers&resultPage=1&resultsPerPage=60	
ThermoFisher	https://www.thermofisher.com/de/de/home/life-science/protein-biology/protein-labeling-crosslinking/protein-crosslinking/crosslinker-selection-	
	guide.html	
ThermoFisher	https://www.thermofisher.com/de/de/home/life-science/protein-biology/protein-labeling-crosslinking/protein-crosslinking.html	
ThermoFisher	https://www.thermofisher.com/de/de/home/life-science/protein-biology/protein-biology-learning-center/protein-biology-resource-library/pierce-	
	protein-methods/chemistry-crosslinking.html	
ThermoFisher	https://www.thermofisher.com/de/de/home/life-science/protein-biology/protein-labeling-crosslinking/protein-crosslinking/crosslinker-application-	
	<u>reference-guide.html</u>	
ThermoFisher	https://www.thermofisher.com/de/de/home/life-science/protein-biology/protein-labeling-crosslinking/protein-	
	labeling/biotinylation/photoreactive-biotinylation-reagents.html	
ThermoFisher	https://www.thermofisher.com/de/de/home/life-science/protein-biology/protein-labeling-crosslinking/protein-crosslinking/crosslinker-selection-	
	<u>tool.html</u>	
ThermoFisher	https://www.fishersci.com/us/en/products/I9C8L0FC/crosslinking-labeling-protein-modification.html	
CovaChem	http://www.covachem.com/protein_crosslinkers.html	
CreativeMolecules	http://www.creativemolecules.com/CM Products.htm	
ProteoChem	http://www.proteochem.com/proteincrosslinkers-c-1.html	
G Biosciences	http://www.gbiosciences.com/ResearchProducts/Protein-Cross-linking.aspx	
G Biosciences	http://www.genotech.com/protein-research/protein-crosslinking-modification.html	
G Biosciences	http://www.gbiosciences.com/Protein-Research/Cross-Linking-Modification/Protein-Cross-Linkers	
Chemometec	http://shop.chemometec.com/?s=cross-linker&post_type=product	
Interchim	http://www.interchim.eu	
Biocompare	http://www.biocompare.com/Protein-Biochemistry/11119-Chemicals/?search=crosslinker	
BioVision	http://www.biovision.com/search/results.html?keywords=crosslinker&x=0&y=0	
Abcam	http://www.abcam.com/products?keywords=crosslinker	
SoltecVentures	http://www.soltecventures.com/PhotoreactivePhenylAzides-c-8	
MolBio	http://www.molbio.com/Heterobi.htm	
MolBio	http://www.molbio.com/HomoNHS.htm	

VWR	https://de.vwr.com/store/product/18498564/cross-linkers-heterobifunctional
VWR	http://www.molbio.com/Homomaleim.htm
VWR	http://www.molbio.com/photoact.htm
GooglePatents	https://www.google.de/?tbm=pts&gws_rd=cr&ei=LCjIV6a7CYnWsAG2vIYw#tbm=pts&q=cross-linker
GooglePatents	https://www.google.de/?tbm=pts&gws_rd=cr&ei=LCjIV6a7CYnWsAG2vIYw#tbm=pts&q=crosslinker
ChemicalBook	http://www.chemicalbook.com/Search_EN.aspx?keyword=crosslinker
ChemicalBook	http://www.chemicalbook.com/ChemicalProductProperty_EN_CB9142806.htm
ChemicalBook	http://www.chemicalbook.com/Search_EN.aspx
ChemCD	http://www.chemcd.com/search?q=cross-linker
BioWorld	https://www.bio-world.com/index.php?main_page=advanced_search_result&search_in_description=7&keyword=cross-linker&x=29&y=16
VWR	https://de.vwr.com/store/product/18498564/cross-linkers-heterobifunctional
CaymanChem	https://www.caymanchem.com/Search?q=crosslinker
SoluLink	http://www.solulink.com/search/google?cx=001277416618431018690%3Ahgffxulp8iy&cof=FORID%3A11%3BNB%3A1&query=crosslinker&op=Go&
	form build id=form-
	7d49ffbea69c5b09c17b3c3c176a8bd0&form_token=fd2e01ee5d8c6941c678b5571f6972df&form_id=google_cse_searchbox_form
TRC (Toronto	http://www.trc-canada.com/products-listing/
Research	
Chemicals)	

Supplementary Table S2: Used vendor publications.

Vendor	Publication, available on the web as PDF files	web page (accessed 02/2020)
Creative Molecules Inc.	Crosslinking reagents test reaction procedure	http://www.creativemolecules.com
Creative Molecules Inc.	Isotopically-coded crosslinking reagents from Creative Molecules Inc., 2009	http://www.creativemolecules.com
Creative Molecules Inc.	Crosslinking reagents masses	http://www.creativemolecules.com
GBiosciences	Protein Cross Linker Selection Guide	http://www.gbiosciences.com
GBiosciences	Protein Cross-Linkers Handbook & Selection Guide	http://www.gbiosciences.com
GBiosciences	Protein Cross-Linking Techer's Guidebook (Cat. # BE-605)	http://www.gbiosciences.com
HiChrom	Gas Chromatography GC Derivatization Reagents	http://www.hichrom.com
Invitrogen	Molecular Probes Handbook, Chapter 5: Crosslinking and Photoactivatable Reagents, 11 th . ed., 2010	https://www.thermofisher.com
Interchim	Derivatization reagents	http://www.interchim.com
Mobitec	Chapter 5: Crosslinking and Photoreactive Reagents	http://www.mobitec.com/cms/index.html
Pierce	Cross-Linking Reagents, Technical Handbook, 2009	https://www.thermofisher.com
Regis	GC Derivatization Reagents	http://www.registech.com
Sigma Aldrich	Derivatization reagents: For selective response and Detection in Complex Matrices	https://www.sigmaaldrich.com
Sigma Aldrich	Guide to Derivatization Reagents for GC	https://www.sigmaaldrich.com
Thermo Scientific	Avidin-Biotin Technical Handbook	https://www.thermofisher.com
Thermo Scientific	Crosslinking Technical Handbook, 2012	https://www.thermofisher.com
Thermo Scientific	Thermo Scientific Pierce Reagents, Solvents and Accessories	https://www.cromlab.es

Supplementary Table S3: Used vendor-specific cross-linker selection tools resp. guides.

Vendor	Link to cross-linker selection tool / guide (accessed 02/2020)
Thermo Scientific	https://www.thermofisher.com/de/de/home/life-science/protein-biology/protein-labeling-crosslinking/protein-crosslinking/crosslinker-selection-
	<u>tool.html</u>
Thermo Scientific	https://www.thermofisher.com/de/de/home/life-science/protein-biology/protein-labeling-crosslinking/protein-crosslinking/crosslinker-selection-
	guide.html
GBiosciences	http://info.gbiosciences.com/blog/bid/160707/What-is-Protein-Cross-Linking-and-Which-Reagents-are-Used-in-it

Supplementary Table S4: Amino acid functional groups amenable to reactive groups of cross-linkers

Amino acid	Functional group	is reacting primarily with the following reactive groups of
	of amino acid	cross-linkers
S - serine		
T - threonine	hydroxyl (-OH)	isocyanate, vinyl sulfone
Y - tyrosine		
K - lysine		(form)aldehyde, aryl halide, carbodiimide, carboxyl,
N - asparagine	amine (–NH2),	hydroxybenzotriazole, phosphine, imidoester, N-
Q - glutamine	amides	hydroxyphthalimide, NHS ester, Sulfo-NHS ester, pentafluorophenyl, thioimidate, vinyl sulfone
R - arginine	guanidino	glyoxal, phenylglyoxal, formaldehyde
C - cysteine	sulfhydryl -SH	acryl, haloacetyl (bromo-, iodo-), platinum(II), maleimide, methanethiosulfonate, pyridinyldisulfide, vinyl sulfone
D – aspartic acid	carboxyl	amine, carbodiimide, dimethylether
E – glutamic acid	(-COOH)	
M - methionine	thioether (–S-)	platinum(II)
H - histidine	imidazole	platinum(II), formaldehyde
W - tryptophan	indole	formaldehyde
F - phenylalanine	phenyl	formaldehyde

Supplementary Table S5: Examples of typical spacer arm components, which are often used in cross-linking reagents

Spacer arms	
6-amino-caprionic acid	
DADPA (diaminodipropylamine)	
ethylenediamine	
hexanediamine	
(PEG) _n - polyethylenglycol	

Supplementary Table S6: Handles for detection, affinity enrichment or purification used in cross-linking reagents

Handle	Detection resp. enrichment / purification with		
Biotin	avidin, streptavidin or neutravidin		
alkyne	alkyne-azido click chemistry or immobilized azide linkers		
azide (-N3)	immobilized alkyne linkers or click chemistry using biotin-alkyne respective		
	biotin-phosphine		
CHCA (cyano-4-	UV-absorbing label for signal enhancement		
hydroxycinnamic acid)			
DNB (dinitrobenzole)	anti-DNB antibodies		

Supplementary Table S7: Term statistics of the XLMOD.obo ontology

Term type	Number of terms
Typedef's	16
Cross-linking attributes	16
Handles	10
Reactive group terms	39
Reactive group attributes (without reactivities)	8
Reactivity terms	37
Cross-linker related chemical modification terms	103
Zero-length cross-linker	6
Bifunctional cross-linker terms	248
Trifunctional cross-linker terms	13
Label transfer reagents	34
Amino acid derivatives	5
Amino acids	4
GC/MS derivatization reagents	174
LC/MS derivatization reagents	324
other terms	13
Nucleobases and their derivatives	
Total number of terms	1083

Supplementary Table S8: ID ranges for the terms in the XLMOD.obo ontology

Term type	ID range starting with
General terms, reactivities, handles, other properties	xlmod:00000
Concrete reactive groups and their attributes	xlmod:00100
Mono-functional cross-linker (cross-linker related chemical modifications)	xlmod:01000
Label transfer reagents	xlmod:01800
Analogues of amino acids	xlmod:01900
Bi-functional cross-linkers	xlmod:02000
Tri-functional cross-linkers	xlmod:03000
Derivatization reagents	xlmod:06000
Further reactivities (used for GC/MS and LC/MS reagents)	xlmod:06500
GC/MS derivatization reagents	xlmod:07000
LC/MS derivatization reagents	xlmod:09000