Ultra-Cold Cryogenic TEM with Liquid Helium and High Stability

Emily Rennich^{1,2,*}, Suk Hyun Sung^{1,3,*}, Nishkarsh Agarwal³, Maya Gates², Robert Kerns⁴, Robert Hovden^{3,5,†}, and Ismail El Baggari^{1,†}

¹The Rowland Institute at Harvard, Harvard University, Cambridge, MA 02142, United States. ²Department of Mechanical Engineering, University of Michigan, Ann Arbor, MI 48109, United States. ³Department of Materials Science and Engineering, University of Michigan, Ann Arbor, MI 48109, United States. ⁴Michigan Center for Materials Characterization, University of Michigan, Ann Arbor, MI 48109, United States. ⁵Applied Physics Program, University of Michigan, Ann Arbor, MI 48109, United States. ^{*}Contributed equally. [†]e-mail: hovden@umich.edu, ielbaggari@rowland.harvard.edu

Cryogenic transmission electron microscopy has revolutionized structural biology and materials science, but achieving temperatures below the boiling point of liquid nitrogen remains a long-standing aspiration. We introduce an ultra-cold liquid helium transmission electron microscope specimen holder, featuring continuous cryogen flow and vibration decoupling. This instrument is compatible with modern aberration-corrected microscopes and achieves sub-25 K base temperature, ±2 mK thermal stability over many hours, and atomic resolution—setting the stage for a new era of cryogenic electron microscopy.

Rapid advancements of cryogenic transmission electron microscopy (TEM) techniques have revolutionized the biological [1,2] and materials sciences [3–6] over the past decades. However, serious challenges persist in the realm of cryogenic TEM techniques, particularly when working at temperatures below 100 K. Ultra-cold temperatures below the boiling point of liquid nitrogen (77 K) in modern TEM are poised to improve the radiation tolerance and resolution in single-particle cryo-EM and cryo-electron tomography [7– 10], and provide access to quantum phases in materials science research [11–15].

Dedicated helium stage modules were constructed and implemented early in the history of TEM [16-20]; however, this approach is currently incompatible with highresolution objective lens pole pieces due to the strict design constraints of modern electron optics. In life science applications, cartridge-based liquid nitrogen microscopes have become popular even though they lack ultra-cold performance and in situ capabilities. As an alternative, dewarbased side-entry cryogenic holders have filled the need to cool specimens to low temperatures. Unfortunately, the rapid evaporation of the cryogen within the dewar introduces large mechanical vibrations and thermal instabilities that result in poor imaging resolution. Moreover, the relatively small volume of these dewars leads to short operation times, or hold times [13,14]. While some of these challenges can be mitigated for liquid nitrogen, liquid helium has a heat of vaporization \sim 60 times lower than that of liquid nitrogen, which dramatically increases instabilities and has made the use of liquid helium a persistent challenge for cryogenic TEM [14].

Here we introduce a side-entry specimen holder capable of liquid helium cooling while achieving atomic resolution. This ultra-cold cryogenic TEM specimen holder demonstrates millikelvin temperature stability maintained for hours near liquid helium temperatures (< 25 K at the specimen) in a modern aberration-corrected TEM. For the first time, a base temperature of 23 K is maintained over 4

hours. Moreover, the thermal stability is better than $\pm 2 \text{ mK}$, a significant (estimated $10 \times$) improvement over existing instruments. Thermal stability is directly associated with the reduction of specimen drift. By maintaining stable cryogenic temperatures with liquid helium for extended durations, scientists can conduct intricate experiments without grappling with unpredictable temperature fluctuations or recurrent dewar refills.

Figure 1a shows the atomic structure of gold at \sim 31 K using this ultra-cold TEM specimen holder inserted into a double-corrected JEOL 3100R05 operating at 300 kV. This holder (Fig. 1g) employs liquid cryogen flow through a copper heat exchanger as the cooling source. Liquid helium continuously flows from a large external dewar along a vacuum-insulated and radiation-shielded transfer line to the copper heat exchanger. The sample is coupled to the cooled heat exchanger via highly thermally conductive components running axially inside the long hollow metal rod of the holder, keeping it cold as long as liquid helium is flowing. The extended surface area of these cold components act as internal cold traps in the vacuum space, preventing ice accumulation on the specimen even after multiple hours of experiments.

Vibrations are mitigated by decoupling the specimen from the helium transfer line through two pairs of flexible ultra-high vacuum edge-welded bellows with a spring constant of 4.7 N/mm. Each bellow compresses 50A durometer rubber blocks that dampen vibrations emanating from the transfer line and room environment. The flexible bellows also accommodate movements of the microscope's goniometer, enabling the usual positional adjustments in the x, y and z axes, as well as axial rotation.

In addition to low base temperatures, this liquid helium cryogenic TEM holder allows precise control over a wide range of temperatures without compromising thermal stability. By adjusting the flow of helium and applying carefully controlled heating, the holder can be stabilized anywhere between room and base temperature with a stability




Fig. 1 | Ultra-cold cryogenic TEM imaging at atomic resolution with mK temperature stability. a), b Atomic resolution image of gold at 31 K, and 115 K, respectively. The top insets represent a zoom-in into the TEM images. Bottom insets show the Fourier transform of the atomically resolved TEM images. The scale bar represents 2 nm. c,d) Histograms of measured sensor temperatures at 11 K and 110 K, respectively. A temperature stability of ± 1.4 mK was measured at a heat exchanger base temperature of 11 K over a data collection period of 4 hours. A temperature stability of ± 2.0 mK was measured at a nitermediate heat exchanger temperature of 110 K over a data collection period of 4 hours. e),f) Plot of the measured heat exchanger temperature of 110 K over a data collection period of 4 hours. e),f) Plot of the measured heat exchanger temperature of 110 K over a data collection period of 4 hours. e),f) Plot of the measured heat exchanger temperature fluctuations as a function of time at 11 K and 100 K, respectively. Multi-hour long hold times with ≤ 2 mK stability were demonstrated at both temperatures. g) Rendered model of the instrument components including the cryogen-flow heat exchanger, the vibration damping module, and TEM specimen tip. Liquid helium flow through the heat exchanger is denoted using an arrow with relative holder temperature denoted by a blue-white gradient overlay. Temperature measurements were taken using a calibrated silicon diode sensor at the specimen tip and liquid flow heat exchanger. Temperature difference between heat exchanger and specimen tip is ~15 K.

of ± 2 mK or better. This capability has been demonstrated, reaching ± 1.4 mK and ± 2.0 mK stability at temperatures of 11 K and 110 K, respectively (see Fig. 1c-f). This instrument is also cryogen-agnostic: liquid nitrogen can be used instead of liquid helium if only temperatures above 110 K are needed.

Millikelvin stability over 10 hours is shown in the temperature plot (Fig. 1f) at an intermediate temperature of 110 K, though similarly long hold times are possible at any operating temperature. In fact, hold times are only limited by the available helium supply—this holder consumes helium at a rate of less than two litres per hour. Commonly available helium tanks are between 30–1000 liters, making long hold times easily achievable. These long hold times present a substantial opportunity for measurements that require long experimental durations including single particle reconstruction, chemical spectroscopy, tomography, and ptychography. The ability to easily reach intermediate temperatures also affords a great deal of flexibility: materials can now be reliably characterized at any desired temperature, not just near the boiling point of liquid nitrogen.

Beyond applications for the preservation of beamsensitive specimens, ultra-cold TEM opens up avenues for investigating quantum materials, crystalline solids in which novel electronic states emerge at low temperatures. An illustrative example is the observation of a structural phase transition below 33 K in the layered quantum material 2H-NbSe₂ (Fig. 2), which further validates ultra-cold TEM performance. 2H-NbSe₂ is a layered transition-metal dichalcogenide that exhibits charge density wave (CDW) ordering below 33 K and superconductivity below 7 K [21]. CDW formation is characterized by the appearance of superlattice peaks in reciprocal space. Figure 2a shows an electron diffraction pattern at 300 K from a thin exfoliated flake.




Fig. 2 | **The emergence of a charge density wave superlattice in 2H-NbSe₂ below 33 K.** a) TEM electron diffraction pattern captured at 300 K shows the room temperature phase containing only Bragg peaks corresponding to a c) normal crystal lattice. b) At temperatures below 33 K superlattice peaks (blue circles) appear, indicating emergence of charge density waves (denoted by blue intensities) with a threefold supercell in d). The 1/3 reciprocal lattice unit wavevector is represented by a blue arrow. c) Simplified model of the hexagonal atomic arrangement of Nb atoms (dark circles) in 2H-NbSe₂. Dashed lines represent the unit cell at high temperature. The lattice constant representing the Nb-Nb distance is 0.34 nm. d) Model of the atomic arrangement overlaid of the charge density wave superlattice (dashed lines). The light blue circles represent the extra charge density that forms a supercell at low temperature. e) Temperature-dependent resistance curve shows the emergence of electronic phase transitions including a CDW transition appears as a slight anomaly in the temperature-resistance curve whereas the superconducting transition appears as a sudden drop to zero resistance. The scale bars represent 2 nm⁻¹.

Bragg peaks with the hexagonal symmetry of the crystal are prominent. Accessing temperatures below 33 K with the ultra-cold holder, we observe the emergence of additional superlattice peaks (Fig. 2b, circles) indicating the formation of the CDW. These peaks are located at ¹/₃ reciprocal lattice units, the predicted wavevector for CDW in 2H-NbSe₂.

The results herein represent a significant advancement that addresses a longstanding desire to access ultra-cold temperatures within a modern aberration-corrected electron microscope, offering new opportunities for probing challenging biological and quantum materials. Cryogenic electron microscopy at ultra-cold temperatures is critical to improving radiation tolerance and imaging resolution of key materials underlying biological and renewable energy research. This development also enables a new field for characterization of quantum materials, including mapping of phase transitions such as superconductivity at temperatures unattainable with liquid nitrogen cooling.

References

- ¹W. Kühlbrandt, "The resolution revolution", Science **343**, 1443–1444 (2014).
- ²K. M. Yip, N. Fischer, E. Paknia, A. Chari, and H. Stark, "Atomicresolution protein structure determination by cryo-EM", Nature 587, 157–161 (2020).
- ³Y. Li, Y. Li, A. Pei, K. Yan, Y. Sun, C.-L. Wu, L.-M. Joubert, R. Chin, A. L. Koh, Y. Yu, et al., "Atomic structure of sensitive battery materials and

interfaces revealed by cryo-electron microscopy", Science 358, 506–510 (2017).

- ⁴M. J. Zachman, Z. Tu, S. Choudhury, L. A. Archer, and L. F. Kourkoutis, "Cryo-stem mapping of solid–liquid interfaces and dendrites in lithiummetal batteries", Nature 560, 345–349 (2018).
- ⁵I. El Baggari, B. H. Savitzky, A. S. Admasu, J. Kim, S.-W. Cheong, R. Hovden, and L. F. Kourkoutis, "Nature and evolution of incommensurate charge order in manganites visualized with cryogenic scanning transmission electron microscopy", Proc. Natl. Acad. Sci. 115, 1445–1450 (2018).
- ⁶T. Zhu, H. Sternlicht, Y. Ha, C. Fang, D. Liu, B. H. Savitzky, X. Zhao, Y. Lu, Y. Fu, C. Ophus, et al., "Formation of hierarchically ordered structures in conductive polymers to enhance the performances of lithium-ion batteries", Nat. Energy **8**, 129–137 (2023).
- ⁷R. Henderson, "Realizing the potential of electron cryo-microscopy", Q. Rev. Biophys. **37**, 3–13 (2004).
- ⁸J. Dubochet, F. P. Booy, R. Freeman, A. V. Jones, and C. A. Walter, "Low temperature electron microscopy", Annu. Rev. Biophys. **10**, 133– 149 (1981).
- ⁹K. Naydenova, A. Kamegawa, M. J. Peet, R. Henderson, Y. Fujiyoshi, and C. J. Russo, "On the reduction in the effects of radiation damage to two-dimensional crystals of organic and biological molecules at liquidhelium temperature", Ultramicroscopy 237, 113512 (2022).
- ¹⁰C. J. Russo, J. L. Dickerson, and K. Naydenova, "Cryomicroscopy in situ: what is the smallest molecule that can be directly identified without labels in a cell?", Faraday Discuss. 240, 277–302 (2022).
- ¹¹K. Harada, T. Matsuda, J. Bonevich, M. Igarashi, S. Kondo, G. Pozzi, U. Kawabe, and A. Tonomura, "Real-time observation of vortex lattices in a superconductor by electron microscopy", Nature 360, 51–53 (1992).
- ¹²W. Zhao, M. Li, C.-Z. Chang, J. Jiang, L. Wu, C. Liu, J. S. Moodera, Y. Zhu, and M. H. W. Chan, "Direct imaging of electron transfer and its influence on superconducting pairing at FeSe/SrTiO₃ interface", Sci. Adv. 4, eaao2682 (2018).
- ¹³A. M. Minor, P. Denes, and D. A. Muller, "Cryogenic electron microscopy for quantum science", MRS Bull. 44, 961–966 (2019).

- ¹⁴Y. Zhu, "Cryogenic electron microscopy on strongly correlated quantum materials", Acc. Chem. Res. 54, 3518–3528 (2021).
- ¹⁵K. Naydenova, P. Jia, and C. J. Russo, "Cryo-EM with sub–1 Å specimen movement", Science 370, 223–226 (2020).
- ¹⁶J. A. Venables, "Liquid helium cooled tilting stage for an electron microscope", Rev. Sci. Instrum. 34, 582–583 (1963).
- ¹⁷V. R. Matricardi, W. G. Lehmann, N. Kitamura, and J. Silcox, "Electron microscope observations of ferromagnetic domains in chromium tribromide", J. Appl. Phys. 38, 1297–1298 (1967).
- ¹⁸Y. Fujiyoshi, T. Mizusaki, K. Morikawa, H. Yamagishi, Y. Aoki, H. Kihara, and Y. Harada, "Development of a superfluid helium stage for highresolution electron microscopy", Ultramicroscopy 38, 241–251 (1991).
- ¹⁹W. Jiang, M. L. Baker, J. Jakana, P. R. Weigele, J. King, and W. Chiu, "Backbone structure of the infectious ε15 virus capsid revealed by electron cryomicroscopy", Nature 451, 1130–1134 (2008).
- ²⁰F. Börrnert, F. Kern, F. Harder, T. Riedel, H. Müller, B. Büchner, and A. Lubk, "The Dresden *in-situ* (S)TEM special with a continuous-flow liquid-helium cryostat", Ultramicroscopy **203**, 12–20 (2019).
- ²¹M. M. Ugeda, A. J. Bradley, Y. Zhang, S. Onishi, Y. Chen, W. Ruan, C. Ojeda-Aristizabal, H. Ryu, M. T. Edmonds, H.-Z. Tsai, A. Riss, S.-K. Mo, D. Lee, A. Zettl, Z. Hussain, Z.-X. Shen, and M. F. Crommie, "Characterization of collective ground states in single-layer NbSe₂", Nat. Phys. 12, 92–97 (2016).

Methods

Cryogenic Experiment

The holder incorporates liquid helium flow cooling and flexible thermal and vacuum connections to reduce mechanical vibration coupling. The internal cryogenic components operate under high vacuum conditions provided by the microscope column vacuum ($\sim 1.0 \times 10^{-7}$ torr). Due to the increase in internal surface area, the holder is pumped under vacuum for at least 30 minutes before imaging. Internally, a copper heat exchanger is connected to a copper braid and rod assembly that traverses the bore of the TEM holder.

A liquid helium transfer line connects the specimen holder to a large external helium source—here a 60 L stainless steel dewar was used. The dewar is depressurized when the helium transfer line is inserted and secured to the dewar. Pressure in the dewar is then raised and maintained at \sim 5 psi to ensure flow into the ultra-cold specimen holder.

For precise temperature control, we use a flexible heating element attached to the heat exchanger to maintain the desired temperature. The temperature is measured using a silicon diode with an excitation DC current of 10 μ A and controlled using a temperature controller (Cryo-Con Model 22C) with a proportional-integral-derivative (PID) scheme. A second temperature sensor is mounted at the tip, near specimen to monitor specimen temperature.

Specimen Preparation

Sputtered gold on an ultrathin carbon film covered a standard 3 mm copper mesh grid. 2H-NbSe₂ was exfoliated onto polydimethylsiloxane (PDMS) gel stamp and mechanically transferred using a home-built transfer system onto silicon-based TEM grids. The silicon based TEM grids contain a SiN_x membrane window TEM grid with 2 μ m holes and were produced by Norcada.

Temperature-Dependent Resistance Measurements

Temperature-dependent resistance measurements on 2H-NbSe₂ were performed in a 1.8 K cryostat with a superconducting magnet (Quantum Design PPMS DynaCool). Samples were prepared by exposing a fresh surface of 2H-NbSe₂ by scotch-tape exfoliation and thin gold wires were attached using silver paint. Resistance measurements were done in a 4-probe configuration with a 1 mA current being sourced from a Keithley 6220 Current Source and voltage measured via a Keithley 2812A Nanovoltmeter. The sample was cooled in helium exchange gas. The measurements shown here were performed under a 2 T out-of-plane magnetic field to match the magnetic field experienced by the sample in a TEM. Due to this magnetic field, the superconducting transition is suppressed to ~ 5 K instead of the \sim 7 K transition temperature under ambient conditions. The charge density wave transition is unaffected by the magnetic field.

Transmission Electron Microscopy

The in situ cooldown was performed on a double Cs-Corrected JEOL3100R05 STEM/TEM equipped with tungsten cold-FEG operated at 300 kV held under a $\sim 5.2 \times 10^{-11}$ torr vacuum. The microscope column pressure of $\sim 7.5 \times 10^{-8}$ torr was maintained around the specimen throughout the experiment. TEM images (Fig. 1a,b) were acquired on Gatan Ultrascan1000 CCD camera with 25 ms acquisition time. Electron diffraction patterns were collected using an 850 nm selected area aperture.

Acknowledgements

The first liquid helium cool-down occurred on July 1, 2022. Sub-30 K CDW phase transition 2H-NbSe₂ with hold times in excess of ten hours was observed on August 3, 2022. Atomic resolution images at 23 K where acquired on July 18, 2023. This research was supported by the Rowland Institute at Harvard. We acknowledge key contributions from Winfield Hill, Erik Madsen, Kal Banger, Alan Stern and Chris Stokes from the Rowland Institute at Harvard for their help with machining, electronics design, and cryogenics. Electron microscopy experiments were conducted using the Michigan Center for Materials Characterization $((MC)^2)$ at the University of Michigan with assistance from Tao Ma, Haiping Sun, and Allen Hunter. I.E. and E.R. acknowledge support from the Rowland Institute at Harvard. R.H. acknowledges support from the U.S. Department of Energy, Basic Energy Sciences, under award DE-SC0024147. E.R. acknowledges support from the National Science Foundation award 2228909.

Author contributions

I.E., R.H., E.R. conceived the design and experiments. E.R., S.H.S., N.A., M.G., R.H, and I.E. performed cooldown experiments and in situ microscopy measurements. I.E., S.H.S., N.A., M.G., E.R., R.K., and R.H. built the experimental microscopy setup and performed hardware assembly. S.H.S. and N.A. optimized temperature control. E.R., S.H.S., R.H., and I.E. prepared the manuscript. All authors reviewed and edited the manuscript.

Competing interests

E.R., R.H., I.E. are inventors on a filed patent application related to this work and licensed for commercialization. The other authors declare no competing interests.