MCL: Multi-view Enhanced Contrastive Learning for Chest X-ray Report Generation

Kang Liu¹, Zhuoqi Ma^{1,2,*}, Kun Xie¹, Zhicheng Jiao², Qiguang Miao^{1,*}

¹Xidian University

²Brown University

kangliu422@gmail.com, {zhuoqima,xiekun,qgmiao}@xidian.edu.cn,zhicheng_jiao@brown.edu

Abstract

Radiology reports are crucial for planning treatment strategies and enhancing doctor-patient communication, yet manually writing these reports is burdensome for radiologists. While automatic report generation offers a solution, existing methods often rely on single-view radiographs, limiting diagnostic accuracy. To address this problem, we propose MCL, a Multi-view enhanced Contrastive Learning method for chest X-ray report generation. Specifically, we first introduce multi-view enhanced contrastive learning for visual representation by maximizing agreements between multi-view radiographs and their corresponding report. Subsequently, to fully exploit patient-specific indications (e.g., patient's symptoms) for report generation, we add a transitional "bridge" for missing indications to reduce embedding space discrepancies caused by their presence or absence. Additionally, we construct Multi-view CXR and Two-view CXR datasets from public sources to support research on multi-view report generation. Our proposed MCL surpasses recent state-of-the-art methods across multiple datasets, achieving a 5.0% F₁ Rad-Graph improvement on MIMIC-CXR, a 7.3% BLEU-1 improvement on MIMIC-ABN, a 3.1% BLEU-4 improvement on Multi-view CXR, and an 8.2% F1,mic-14 CheXbert improvement on Two-view CXR.

Introduction

Radiology reports, crafted by experienced radiologists, meticulously document imaging findings from examinations such as X-rays, PET scans, and CTs, detailing abnormalities and initial diagnostic conclusions. These reports deliver vital imaging insights that enable physicians to develop efficient patient treatment strategies (Messina et al. 2022). However, the manual writing process is time-consuming and requires significant expertise, making it increasingly difficult to meet the demands of modern healthcare (Bannur et al. 2024), particularly in regions with limited medical resources.

Automatic chest X-ray report generation aims to produce detailed and accurate free-text reports from multi-view radiographs, helping radiologists improve diagnostic efficiency and consistency by providing high-quality draft reports. In clinical practice, limitations like X-ray equipment constraints and the complexity of human anatomical structures can prevent a single-view radiograph from achieving



Figure 1: Comparison of existing methods and our proposed method for chest X-ray report generation. Existing methods generate reports from single-view radiographs, while our approach utilizes multi-view radiographs.

optimal imaging quality and adequately displaying the overall anatomical structure. As a result, multi-view imaging examinations, such as postero-anterior (PA), antero-posterior (AP), and lateral views, are crucial for accurate diagnostics and personalized treatment. Consequently, the number of radiographs varies across studies. Typically, each study comprises a collection of radiographs, a corresponding report, and a patient-specific indication (which may sometimes be absent). This variability is also evident in public chest X-ray report generation datasets such as MIMIC-CXR (Johnson et al. 2019) and IU X-ray (Demner-Fushman et al. 2016). Such variability makes it challenging to effectively utilize multi-view radiographs from the same study to enhance the clinical accuracy of generated reports. One intuitive approach (Li et al. 2019; Chen et al. 2020; Liu et al. 2024c,b) treats each radiograph as an individual study, generating reports from single-view radiographs (see Figure 1). However, this scheme fails to fully exploit the rich anatomical information available in multi-view radiographs, potentially leading to inaccurate and inconsistent reports. In addition, since the IU X-ray dataset predominantly includes studies with two-view radiographs, several studies (Chen et al. 2020, 2021; Yang et al. 2023) have developed two-view re-

^{*}Corresponding author

port generation approaches, showing promise in producing informative reports. Despite this, these methods often struggle to integrate into clinical workflows due to the varying number of radiographs per study.

To address this challenge, we introduce a novel two-stage method, called Multi-view enhanced Contrastive Learning (MCL) for generating chest X-ray reports. In the representation learning stage, we propose multi-view enhanced contrastive learning for visual representation by leveraging semantic correspondences between multi-view radiographs within the same study and between these radiographs and their corresponding report. In the report generation stage, we utilize the cross-attention mechanism to exploit available indications fully, providing the model with patient background information. Additionally, we incorporate a transitional "bridge" for missing indications to reduce embedding space differences caused by the presence or absence of these indications. For multi-view enhanced contrastive learning, we first employ multi-positive contrastive learning to bring multi-view radiographs within the same study closer, improving the consistency of visual features. Then, we develop a multi-view fusion module to integrate varying numbers of radiographs per study, producing fused visual features for subsequent cross-modal alignment. Finally, we apply contrastive learning for instance-wise and token-wise alignment, maximizing agreement between multi-view radiographs and their corresponding report. We evaluate our proposed MCL on the MIMIC-CXR (Johnson et al. 2019), MIMIC-ABN (Ni et al. 2020; Hou et al. 2023a), and our curated Multi-view CXR and Two-view CXR datasets. Experiment results demonstrate the effectiveness of MCL, showing a 5.0% F₁ RadGraph improvement on MIMIC-CXR, a 7.3% BLEU-1 improvement on MIMIC-ABN, a 3.1% BLEU-4 improvement on Multi-view CXR, and an 8.2% F1,mic-14 CheXbert improvement on Two-view CXR. Our key contributions are outlined as follows:

- We propose a novel multi-view enhanced contrastive learning that facilitates cross-modal alignment between multi-view radiographs and their corresponding report, addressing the issue that existing algorithms cannot handle varying numbers of views.
- We incorporate a transitional "bridge" for missing indications to reduce embedding space differences caused by their presence or absence, thereby enhancing the capture of patient background information.
- We curate Multi-view CXR and Two-view CXR datasets from two public sources, ensuring that each study includes multiple radiographs. This supports research on multi-view report generation, particularly for scenarios involving varying numbers of views or two-view setups.

Related Work

Image Captioning

Image captioning (Li et al. 2022; Chen, ZHANG, and Hinton 2023) and report generation (Hou et al. 2023b; Bu et al. 2024) both convert visual information into textual descriptions. Advances in image captioning, particularly with the encoder-decoder framework, have influenced report generation. However, report generation differs by requiring the capture of subtle differences between radiographs and producing detailed and lengthy reports.

Medical Image Analysis Meets Multi-view Learning

Multi-view learning (Xu, Tao, and Xu 2013; Zhao et al. 2017) aims to enable models to capture shared and complementary information from different perspectives of the same scene, leading to a more comprehensive understanding. In representation learning, Paul et al. (Paul et al. 2021) develop a multi-view semantic embedding extracted from X-ray reports, CT reports, and visual traits to enhance the clinical accuracy of chest X-ray diagnosis. REFERS (Zhou et al. 2022) introduces a view fusion module that integrates fixedperspective visual features, enhancing cross-modal alignment. In medical report generation, FMVP (Liu et al. 2024f) incorporates single radiographs and auxiliary inputs (i.e., disease tags and medical concepts) as multi-view information, using two shared and continuous cross-attention mechanisms (Vaswani et al. 2017) to assist in generating informative reports. Given that the IU X-ray dataset (Demner-Fushman et al. 2016) primarily comprises studies with twoview radiographs, several approaches (Chen et al. 2020, 2021; Yang et al. 2023) have been developed to handle twoview report generation. These methods enhance the clinical efficacy of generated reports and offer valuable insights for exploring multi-view report generation. However, these approaches often struggle with the varying number of views encountered in real-world scenarios. To tackle this issue, we introduce a multi-view enhanced contrastive learning method that accommodates an unfixed number of views.

Medical Image Analysis Meets Contrastive Learning

Contrastive learning is pivotal in medical image analysis, particularly in tasks related to medial visual representation learning (Wang et al. 2022a) and report generation (Shen et al. 2024). It aligns radiographs with their corresponding reports by minimizing the distance between positive pairs (radiographs and associated reports) while maximizing the distance between negative pairs (radiographs and unrelated reports). Significant advancements have been made in this area. For instance, MedCLIP (Wang et al. 2022b) adopts a semantic matching loss for global alignment between decoupled radiographs and reports. ARL (Chen, Li, and Wan 2022) regards knowledge as an intermediary to facilitate semantic alignment between radiographs and their reports. MGCA (Wang et al. 2022a) introduces multi-grained crossmodal alignment at the instance, pathological region, and disease levels for generalized medical visual representation. SEI (Liu et al. 2024c) utilizes global and local cross-modal alignment between radiographs and factual serialization in reports to generate radiology reports. While these methods effectively align radiographs with their reports by extracting visual features from single-view radiographs, they often overlook the benefits of incorporating multi-view information from the same study. In response, our work introduces



Figure 2: Illustration of our proposed MCL, which comprises a visual encoder, a text encoder, and a text generator. MCL is a two-stage method: 1) Representation learning using multi-view enhanced contrastive learning. 2) Report generation based on patient-specific indications. Model inference solely involves stage 2, excluding the language modeling (LM) loss.

multi-view enhanced contrastive learning, integrating multiview visual features to improve cross-modal alignment.

Method

The overall architecture of our proposed MCL, shown in Figure 2, consists of two stages. In the representation learning stage, we propose multi-view enhanced contrastive learning to align multi-view radiographs with their corresponding report. In the report generation stage, we fully exploit patient-specific indications to provide the text generator with patient background information.

Problem Formulation

Let $D_{tr} = \{ (x_{i,k}, x_{i,\backslash k}, z_i, r_i) | 1 \le i \le n, 1 \le k \le m_i \}$ be the training set, where *n* is the number of studies. Each study *i* includes m_i radiographs (views), a patient-specific indication z_i (which may be absent), and a corresponding report r_i . Here, $x_{i,k}$ refers to the anchor scan, while $x_{i,\backslash k} = \{x_{i,j} | j \ne k, 1 \le j \le m_i\}$, which may be empty, represents the set of auxiliary references excluding the k^{th} view. The indication z_i , derived from the radiology report, provides the patient's reasons for the examination. In this paper, we aim to learn the mapping $F_{\theta}(\cdot)$ that maps the input $(x_{i,k}, x_{i,\backslash k}, z_i)$ to output r_i on the training set D_{tr} , namely $F_{\theta}: (x_{i,k}, x_{i,\backslash k}, z_i) \rightarrow r_i$.

Representation Learning with Multi-view Enhanced Contrastive Learning

We introduce multi-view enhanced contrastive learning for visual representation by utilizing semantic correspondences both among multi-view radiographs within the same study and between these radiographs and their corresponding report. First, we employ multi-positive contrastive learning to bring multi-view radiographs closer, thereby enhancing visual feature consistency across them. Afterward, we develop a multi-view fusion module to integrate varying numbers of radiographs per study, generating fused visual features for subsequent cross-modal alignment. Finally, we apply contrastive learning for instance-wise and token-wise cross-modal alignment, maximizing agreement between the radiographs and their report.

Visual Features Extraction. We use ResNet101 (He et al. 2016), pre-trained on ImageNet, as the visual encoder. The feature maps from the last convolutional layer of ResNet101 are regarded as the visual features of radiographs, formulated as $\boldsymbol{V} \in \mathbb{R}^{M \times p \times d_1}$. Here, $M = \sum_{i=1}^{B} m_i$ denotes the number of radiographs in the batch. B, p, and d_1 represent the batch size, the dimensions of the feature map, and the number of channels, respectively.

Textual Features Extraction. Drawing inspiration from (Yan et al. 2023; Liu et al. 2024d), we first extract factual serialization from reports using the structural entities approach (Liu et al. 2024d) to reduce the noise in the following cross-modal alignment. Factual serialization in a report, shown in Figure 2, is a brief sentence composed of clinically significant keyword groups extracted from reports. Following this, factual serialization is input into the text encoder, a six-layer pre-trained SciBERT (Beltagy, Lo, and Cohan 2019), to generate textual features $T \in \mathbb{R}^{B \times k \times d_2}$. Here, k and d_2 denote the number of textual tokens and dimensionality per token, respectively. Note that textual features T are exclusively utilized in the representation learning stage.

Multi-positive Contrastive Learning between Multiview Radiographs. To capture semantic correspondences between multi-view radiographs within the same study, we employ multi-positive contrastive learning (Tian et al. 2023). This method aligns radiographs from the same study while distinguishing them from those in different studies, thereby enhancing the consistency of visual features. We begin by excluding studies with only a single-view radiograph from the batch (Notably, these studies are still used for subsequent cross-modal alignment). For each anchor scan $x_{i,a}$, we compute the contrastive categorical distribution q to measure the similarity between $x_{i,a}$ and all other scans $x_{\backslash a}$. Here, $x_{\backslash a}$ comprises both auxiliary references of $x_{i,a}$ and radiographs from different studies. We denote the global visual features, after ℓ_2 normalization, as $oldsymbol{v}_{i,a}$ for anchor scans and $oldsymbol{v}_{ackslash a}$ for the other scans. The $q \in \mathbb{R}^{K \times (K-1)}$ is formulated as:

$$\boldsymbol{q}_{i} = \frac{\exp\left(\boldsymbol{v}_{i,a} \cdot \left(\boldsymbol{v}_{\backslash a}\right)^{T} \middle/ \tau_{1}\right)}{\sum_{j=1}^{K-1} \exp\left(\boldsymbol{v}_{i,a} \cdot \left(\boldsymbol{v}_{\backslash a}^{j}\right)^{T} \middle/ \tau_{1}\right)}, \qquad (1)$$

where $\tau_1 \in \mathcal{R}^+$ is the temperature parameter. The number of radiographs in a batch is $K = \sum_{j=1}^{B} m_j$, where $m_j > 1$ for all j. Next, we construct the ground-truth categorical distribution $\boldsymbol{p} \in \mathbb{R}^{K \times (K-1)}$, denoted as:

$$\boldsymbol{p}_{i} = \frac{\mathbb{I}_{match}\left(\boldsymbol{x}_{i,a}, \boldsymbol{x}_{\backslash a}\right)}{\sum_{j=1}^{K-1} \mathbb{I}_{match}\left(\boldsymbol{x}_{i,a}, \boldsymbol{x}_{\backslash a}^{j}\right)},$$
(2)

here, $\mathbb{I}_{match}(\cdot, \cdot)$ is an indicator function that determines whether two radiographs belong to the same study. Finally, the multi-positive contrastive loss is calculated using the cross-entropy between q and p, formulated as:

$$\mathcal{L}_{MPC} = -\frac{1}{K} \sum_{i=1}^{K} \boldsymbol{p}_i \log \boldsymbol{q}_i.$$
(3)

Despite the varying number of radiographs per study, the different number of non-zero elements in p_i account for this variability. Applying cross-entropy between q_i and p_i pushes the multi-view radiographs in the i^{th} study closer together in the embedding space.

Multi-view Fusion Module. Effective integration of multi-view visual features is crucial for aligning multi-view radiographs with their corresponding report. However, the varying number of radiographs per study complicates simple concatenation, leading to inconsistent channel dimensions in the fused visual features and hindering cross-modal alignment. Although averaging and maximization methods ensure channel consistency, they result in a loss of detail. To tackle this challenge, we propose a multi-view fusion module that integrates multi-view visual features using the scaled dot-product cross-attention mechanism (Vaswani et al. 2017), denoted as ATTN(Q, K, V). In this module, the anchor scan $x_{i,a}$ functions as queries, while auxiliary references $x_{i,\backslash a}$ serve as both keys and values. A skip connection (He et al. 2016) and layer normalization (LN) are then applied to produce the fused visual features for subsequent cross-modal alignment:

$$\tilde{\boldsymbol{V}}_{i,a} = LN\left(\boldsymbol{V}_{i,a} + ATTN\left(\boldsymbol{V}_{i,a}, \boldsymbol{V}_{i,\backslash a}, \boldsymbol{V}_{i,\backslash a}\right)\right), \quad (4)$$

where $\tilde{V} = \left\{ \tilde{V}_{i,a} | 1 \leq i \leq B \right\} \in \mathbb{R}^{B \times p \times d_1}$ ensures consistent channel dimensions across studies while preserving information from both anchor scans and auxiliary references. Notably, we instruct the module to focus on one study at a time, adeptly handling the variability in the number of radiographs per study.

Instance-wise Alignment Loss. Inspired by (Cheng et al. 2023; Wang et al. 2022a; Liu et al. 2024d), we propose the instance-wise alignment loss to maximize the agreement between multi-view radiographs and their corresponding report in the embedding space. We begin by projecting the fused visual features $ilde{V}$ and textual features T into a unified embedding space d using the two-layer convolutional projection head, respectively. We denote the results as $\bar{V} \in \mathbb{R}^{B \times p \times d}$ and $\bar{T} \in \mathbb{R}^{B \times k \times d}$, respectively. Due to the similarity in radiological features, radiology reports from different studies may be identical. Therefore, besides treating the visual-textual feature pair from the same study as a positive pair, we also consider reports with identical contents from different studies as additional positive pairs, along with the current visual features. Next, the image-totext contrastive categorical distribution $q^{V
ightarrow T} \in \mathbb{R}^{B imes B}$ is computed to estimate the similarity between multi-view radiographs and their corresponding report:

$$\boldsymbol{q}^{V \to T} = \frac{\exp\left(\boldsymbol{\bar{v}} \cdot \boldsymbol{\bar{t}}^{T} / \tau_{2}\right)}{\sum_{j=1}^{B} \exp\left(\boldsymbol{\bar{v}} \cdot \left(\boldsymbol{\bar{t}}^{j}\right)^{T} / \tau_{2}\right)},$$
(5)

here, \bar{v} and \bar{t} are global visual and textual features processed by ℓ_2 normalization. Likewise, we obtain the symmetric text-to-image contrastive categorical distribution $q^{T \to V}$. The global ground-truth categorical distribution p^g is:

$$\boldsymbol{p}_{i,j}^{g} = \frac{\mathbb{I}_{equal}\left(r_{i}, r_{j}\right)}{\sum_{k=1}^{B} \mathbb{I}_{equal}\left(r_{i}, r_{k}\right)},\tag{6}$$

where $p^g \in \mathbb{R}^{B \times B}$. $\mathbb{I}_{equal}(\cdot, \cdot)$ denotes an indicator function that determines whether two reports are identical. The instance-wise alignment loss is formulated as:

$$\mathcal{L}_{global} = -\frac{1}{B} \sum_{i=1}^{B} \left(\boldsymbol{p}_{i}^{g} \log \boldsymbol{q}_{i}^{V \to T} + \boldsymbol{p}_{i}^{g} \log \boldsymbol{q}_{i}^{T \to V} \right).$$
(7)

Notably, our method differs from previous works (Cheng et al. 2023; Wang et al. 2022a) in two key ways: we utilize factual serialization in reports for alignment instead of complete reports, and we employ multi-positive contrastive loss (Tian et al. 2023) to learn semantic correspondences between radiographs and reports, rather than single-positive contrastive loss (van den Oord, Li, and Vinyals 2019).

Token-wise Alignment Loss. Inspired by (Wang et al. 2022a; Liu et al. 2024d), we employ the token-wise alignment loss \mathcal{L}_{local} (Liu et al. 2024d) to learn fine-grained visual features. This loss is achieved through single-positive contrastive learning between token embeddings of visual and textual data. In conclusion, the training objective during the representation learning stage is formulated as:

$$\mathcal{L}_{pretrain} = \mathcal{L}_{MPC} + \mathcal{L}_{global} + \mathcal{L}_{local}.$$
 (8)

Report Generation Based on Patient-specific Indications

Indication Features Extraction. Patient-specific indications, as shown in Figure 2, are collected before examinations and provide critical background information, such as the patient's condition. However, de-identification in public datasets often introduces noise into indications, such as ____, @, and -year-old. We preprocess the data as described in (Liu et al. 2024c), removing the noise and standardizing gender expressions. The cleaned indications are then fed into a text encoder, initialized with the pre-trained model from stage 1, to extract indication features.

Report Generation. Indications may sometimes be absent in certain studies, as detailed in Table 1 in the Appendix. As a result, some existing methods either directly ignore indications (Chen et al. 2020; Wang et al. 2023; Shen et al. 2024) or fail to exploit them fully (Tian et al. 2019; Nguyen et al. 2023). To address this issue, we utilize the Transformer Decoder (Chen, Li, and Wan 2022) to exploit available indications fully, providing the text generator with patient background information. For cases where indications are missing, we incorporate a transitional "bridge" (i.e., Transformer Encoder (Chen, Li, and Wan 2022)) to mitigate embedding space differences caused by the presence or absence of indications. The text generator then produces reports in an autoregressive manner, optimized by minimizing the negative log-likelihood as follows:

$$\mathcal{L}_{LM_i} = -\sum_{t=1}^{L} \log P\left(\left.\tilde{w}_i^t\right| x_{i,a}, x_{i,\backslash a}, z_i, \tilde{w}_i^{j$$

where \tilde{w}_i^t denotes the t^{th} predicted word. $\tilde{w}_i^{j < t}$ is the preceding word sequence generated by the text generator. $P(\cdot)$, z_i , and L represent the likelihood function, the indication, and the maximum number of generated tokens, respectively.

Experiments

Experimental Settings

Datasets. 1) MIMIC-CXR (Johnson et al. 2019) is a largescale chest X-ray dataset with free-text radiology reports and varying numbers of radiographs per study. 2) MIMIC-ABN (Ni et al. 2020), a subset of MIMIC-CXR, focuses on abnormal sentences within radiology reports. 3) Our curated Multi-view CXR aggregates studies with multiple views from both MIMIC-CXR and IU X-ray (Demner-Fushman et al. 2016). 4) Our curated Two-view CXR is a variant of Multi-view CXR that includes only two views per study, with other settings unchanged. We treat the "FINDINGS" section of the radiology report as the reference report, in line with established practices (Chen et al. 2020; Wang et al. 2023; Shen et al. 2024). MIMIC-ABN contains only abnormal sentences in this section. MIMIC-CXR and MIMIC-ABN datasets adhere to their official splits. For Multi-view CXR, we follow the MIMIC-CXR split for its MIMIC-CXR component and baseline split (Chen et al. 2020; Yang et al. 2023; Shen et al. 2024) for its IU X-ray. Statistics for these datasets are summarized in Table 1 in the Appendix.

Evaluation Metrics. We evaluate model performance using natural language generation (NLG) and clinical efficacy (CE) metrics. NLG metrics assess the lexical similarities between generated and reference reports, including BLEU*n* (B-*n*, where $n \in \{1, 2, 3, 4\}$), METEOR (MTR), and ROUGE-L (R-L). CE metrics evaluate the clinical correctness of generated reports, including F₁ RadGraph (RG) (Jain et al. 2021), F_{1,mic-5} CheXbert (CX5), and F_{1,mic-14} CheXbert (CX14) (Smit et al. 2020). RG measures the overlap in clinical entities and their relationships, which better aligns with radiologists' evaluations than B-3 and CX14 metrics (Yu et al. 2023). CX5 and CX14 evaluate the ability to describe 5 and 14 observations (e.g., *Edema*) using the multi-label classification micro-averaged F₁ score.

Implementation Details. The batch size, τ_1 , and τ_2 are set to 32, 0.5, and 0.5, respectively. The maximum number of generated tokens, L, is fixed at 100 for all datasets. The optimal model is selected based on the combined scores of B-4, RG, and CX14 on the validation set. Experimental results from the chosen optimal model are reported on the test set (please refer to the Appendix for more details).

Results and Analyses

Comparison with State-of-the-art Methods. We compare our proposed MCL with 11 state-of-the-art (SOTA) methods: R2Gen (Chen et al. 2020), CMN (Chen et al. 2021), CGPT2 (Nicolson, Dowling, and Koopman 2023), MET (Wang et al. 2023), KiUT (Huang, Zhang, and Zhang 2023), SA (Yan et al. 2023), FMVP (Liu et al. 2024f), MAN (Shen et al. 2024), PMRG (Jin et al. 2024), Med-LLM (Liu et al. 2024e), and HERGen (Wang, Du, and Yu 2024). Results are presented in Table 1, where "Input Size" denotes the image resolution processed by the visual encoder (e.g., 224 refers to 224×224 pixels). MCL outperforms all baselines across every metric on four datasets, as evidenced by a 5.0% F_1 RadGraph improvement on MIMIC-CXR, a 7.3% BLEU-1 improvement on MIMIC-ABN, a 3.1% BLEU-4 improvement on Multi-view CXR, and an 8.2% F_{1,mic-14} CheXbert improvement on Two-view CXR. These results highlight the effectiveness of MCL in generating clinically accurate reports. Additionally, at a resolution of 384×384 , MCL consistently surpasses its performance at 224×224 resolution on most metrics across these datasets, demonstrating that higher image resolution can significantly enhance the model's performance.

Effect of Multi-view Enhanced Contrastive Learning. We evaluate the impact of our multi-view enhanced contrastive learning in stage 1 on model performance, as shown in Table 2. Results show positive effects from multi-positive contrastive loss (MCL vs. (f)), instance-wise alignment loss (MCL vs. (e)), and token-wise alignment loss (MCL vs. (d)). Please see the Appendix for more details.

Effect of Report Generation Based on Indications. We assess the impact of indications in stage 2 on model performance, as detailed in Table 2. Both single-view ((b) vs. (a)) and multi-view methods (MCL vs. (g)) exhibit significant improvements across all metrics, emphasizing the importance of indications in generating accurate reports. Additionally, MCL shows a 3.1% increase in the sum of all

Deteret	Mathad	Veen	Input			NLG M	fetrics ↑			C	E Metric	s ↑
Dataset	Method	rear	Size	B-1	B-2	B-3	B-4	MTR	R-L	RG	CX5	CX14
Comparison with	n single-view r	nethod	S									
	R2Gen [♭]	2020	224	0.353	0.218	0.145	0.103	0.142	0.277	0.207	0.340	0.340
	CMN [♭]	2021	224	0.353	0.218	0.148	0.106	0.142	0.278	0.220	0.461	0.391
	CGPT2 [♭]	2023	384	0.393	0.248	0.171	0.127	0.155	0.286	-	-	0.442
	MET [♭]	2023	-	0.386	0.250	0.169	0.124	0.152	0.291	-	-	0.311
	KiUT [♭]	2023	224	0.393	0.243	0.159	0.113	0.160	0.285	-	-	0.321
	SA♭	2023	256	-	0.184	-	-	-	-	0.228	-	0.394
MIMIC-CXR	FMVP [♭]	2023	224	0.389	0.236	0.156	0.108	0.150	0.284	-	-	0.336
	MAN [♭]	2024	224	0.396	0.244	0.162	0.115	0.151	0.274	-	-	0.389
	PMRG [♭]	2024	224	0.398	-	-	0.112	0.157	0.268	-	-	0.476
	Med-LLM [♭]	2024	224	-	-	-	0.128	0.161	0.289	-	-	0.395
	HERGen [♭]	2024	384	0.395	0.248	0.169	0.122	0.156	0.285	-	-	-
	MCL(Ours)	-	224	0.395	0.262	<u>0.190</u>	0.147	0.167	0.311	0.276	0.557	0.499
	MCL(Ours)	-	384	0.408	0.271	0.197	0.151	0.171	0.313	0.278	0.578	0.517
	$\Delta \uparrow$	-	-	+1.0%	+2.1%	+2.6%	+2.3%	+1.0%	+2.2%	+5.0%	+11.7%	+4.1%
	R2Gen [♯]	2020	224	0.253	0.144	0.092	0.063	0.106	0.229	0.179	0.501	0.442
	CMN♯	2021	224	0.256	0.147	0.095	0.066	0.110	0.230	0.183	0.528	0.460
MIMIC-ABN	MCL(Ours)	-	224	<u>0.310</u>	<u>0.185</u>	<u>0.125</u>	<u>0.090</u>	0.127	0.246	<u>0.214</u>	<u>0.535</u>	<u>0.482</u>
	MCL(Ours)	-	384	0.329	0.196	0.131	0.093	0.134	0.255	0.220	0.545	0.503
	$\Delta \uparrow$	-	-	+7.3%	+4.9%	+3.6%	+2.7%	+2.4%	+2.5%	+3.7%	+1.7%	+4.3%
	R2Gen [♯]	2020	224	0.359	0.225	0.155	0.114	0.143	0.297	0.255	0.431	0.384
	CMN♯	2021	224	0.404	0.252	0.170	0.122	0.160	0.311	0.279	0.475	0.416
Multi-view CXR	MCL(Ours)	-	224	<u>0.413</u>	0.274	<u>0.199</u>	<u>0.152</u>	<u>0.174</u>	<u>0.335</u>	<u>0.328</u>	<u>0.515</u>	<u>0.487</u>
	MCL(Ours)	-	384	0.415	0.276	0.200	0.153	0.177	0.336	0.329	0.557	0.508
	$\Delta \uparrow$	-	-	+1.1%	+2.4%	+3.0%	+3.1%	+1.7%	+2.5%	+5.0%	+8.2%	+9.2%
Comparison with	n two-view me	thods										
	R2Gen [♯]	2020	224	0.346	0.219	0.153	0.113	0.141	0.302	0.267	0.413	0.400
	CMN♯	2021	224	0.387	0.241	0.166	0.122	0.151	0.310	0.268	0.437	0.425
Two-view CXR	MCL(Ours)	-	224	<u>0.393</u>	0.256	<u>0.184</u>	<u>0.140</u>	<u>0.165</u>	<u>0.322</u>	0.304	<u>0.531</u>	<u>0.501</u>
	MCL(Ours)	-	384	0.411	0.270	0.195	0.150	0.172	0.326	<u>0.302</u>	0.547	0.507
	$ \Delta \uparrow$	-	-	+2.4%	+2.9%	+2.9%	+2.8%	+2.1%	+1.6%	+3.6%	+11.0%	+8.2%

Table 1: Comparisons with SOTA methods on MIMIC-CXR, MIMIC-ABN, Multi-view CXR, and Two-view CXR datasets. Δ denotes the performance difference between MCL and the best peer methods. ^b means results cited from the original work, while [#] indicates results reproduced using official codes. The best and second-best values are in **bold** and <u>underlined</u>, respectively.

Model	M/S	Representation Learning				Report Generation		NLG M	letrics ↑	CE Metrics †			
WIGUCI		F/R	MPC	G	L	indications("bridge")	B-2	B-4	MTR	R-L	RG	CX5	CX14
R2Gen	S						0.218	0.103	0.142	0.277	0.207	0.340	0.340
(a)	S	F		\checkmark	\checkmark		0.213	0.102	0.142	0.278	0.233	0.512	0.430
(b)	S	F		\checkmark	\checkmark	$\checkmark (\checkmark)$	0.242	0.130	0.158	0.297	0.256	0.548	0.484
(c)	Μ		\checkmark			$\checkmark (\checkmark)$	0.261	0.146	0.165	0.309	0.270	0.551	0.482
(d)	Μ	F	\checkmark	\checkmark		$\checkmark (\checkmark)$	0.263	0.145	0.166	0.305	0.268	0.540	0.480
(e)	Μ	F	\checkmark		\checkmark	$\checkmark (\checkmark)$	0.246	0.136	0.157	0.305	0.255	0.513	0.450
(f)	Μ	F		\checkmark	\checkmark	$\checkmark (\checkmark)$	0.252	0.142	0.162	0.307	0.266	0.536	0.477
(g)	M	F	\checkmark	\checkmark	\checkmark		0.226	0.108	0.148	0.279	0.230	0.539	0.477
(h)	Μ	R	\checkmark	\checkmark	\checkmark	$\checkmark (\checkmark)$	0.258	0.142	0.164	0.307	0.266	0.522	0.470
(i)	Μ	F	\checkmark	\checkmark	\checkmark	√ (×)	0.259	0.145	0.164	0.312	0.272	0.548	0.488
MCL	М	F	\checkmark	\checkmark	\checkmark	$\checkmark (\checkmark)$	0.262	0.147	0.167	0.311	0.276	0.557	0.499

Table 2: Ablation study on MIMIC-CXR with a resolution of 224×224 . "M/S" indicates a multi-view or single-view method. "F/R" denotes cross-modal alignment using either fact serialization or report. "MPC", "G", and "L" indicate multi-positive contrastive loss, instance-wise alignment loss (G), and token-wise alignment loss (L), respectively. 'Indications("bridge")' denotes whether indications are used in stage 2 and if a transitional "bridge" is added for missing indications.

B-4↑ RG↑ CX14↑ 0.696/0.160 0.857/0.182 0.000/0.000	Reference report: frontal and lateral views of the chest are obtained . the lungs remain hyperinflated , suggesting chronic obstructive pulmonary disease . no focal consolidation , pleural effusion , or evidence of pneumothorax is seen . the cardiac and mediastinal silhouettes are stable and unremarkable . hilar contours are also stable .	Generated report by R2Gen: pa and lateral views of the chest provided . there is no focal consolidation effusion or pneumothorax . the cardiomediastinal silhouette is normal . imaged osseous structures are intact . no free air below the right hemidiaphragm is seen .	Generated report by MCL (Ours): frontal and lateral views of the chest were obtained . the lungs are hyperinflated with flattening of the diaphragms, suggesting chronic obstructive pulmonary disease . no focal consolidation , pleural effusion , or evidence of pneumothorax is seen . the cardiac and mediastinal silhouettes are unremarkable .
B-4 [†] RG [†] CX14 [†] 0.658/0.000 0.762/0.000 0.667/0.500	Reference report: in comparison with the study of, the monitoring and support devices remain in place . diffuse bilateral pulmonary opacifications are essentially unchanged .	Generated report by R2Gen: a single portable ap chest radiograph was obtained . an endotracheal tube terminates in the mid thoracic trachea . a left port-a- cath terminates in the right atrium . a right chest tube terminates in the right apex . a right pleural effusion is small . lung volumes are low . there is no focal consolidation or pneumothorax . cardiac and mediastinal contours are normal .	Generated report by MCL (Ours): in comparison with the study of, the monitoring and support devices remain in place . diffuse bilateral pulmonary opacifications are again seen , consistent with widespread pneumonia and pulmonary edema .

Figure 3: Examples of generated reports on the MIMIC-CXR test set with a resolution of 224×224 . The cell "A/B" represents "MCL/R2Gen". Sentences in the reference report are shown in different colors. Each sentence in generated reports is highlighted in matching colors corresponding to those in the reference report. Failure sentences in MCL are <u>underlined</u>.

metrics compared to (i), confirming the effectiveness of employing a transitional "bridge" for missing indications.

Model Benefits from Multi-view Radiographs and Indications. Table 4 in the Appendix shows performance variations on subsets of the test set with and without multi-view radiographs (indications). Significant improvements in NLG and RG metrics are observed for studies with multi-view radiographs (indications). This suggests that multi-view radiographs offer comprehensive visual features, while indications provide essential background information, enhancing the accuracy and completeness of clinical reports.

Clinical Accuracy of 14 Observations. Tables 2 and 3 in the Appendix illustrate the clinical accuracy of 14 observations related to thoracic diseases and support devices. Results reveal that: 1) MCL significantly improves most F_1 compared to R2Gen (Chen et al. 2020). 2) Although MCL is not specifically designed for imbalanced observations, it slightly outperforms the baseline, improving F_1 for *Pleural Other* and *Fracture* by 12% and 2.3% in MIMIC-CXR, and by 1.7% and 7% in Two-view CXR, respectively.

Case Study. Figure 3 presents two examples from the MIMIC-CXR test set to illustrate our generated reports, with additional examples available in Figure 1 of the Appendix. A color distribution that closely matches the reference reports indicates comprehensive coverage of clinical findings, while longer color bars in generated reports suggest more detailed descriptions. Results show that MCL produces higher-quality draft reports for radiologists compared to R2Gen (Chen et al. 2020). For instance, in case 1, radiologists only needed to add *"hilar contours are also stable"*. However, MCL and R2Gen struggle to effectively capture disease progression due to the lack of consideration for patient temporal data, as seen with *"unchanged"* in case 2.

Evaluation via Large Language Models. We evaluate the performance of our generated reports using the GREEN (Ostmeier et al. 2024) score, which employs a fine-tuned LLaMA-2 (Touvron et al. 2023) model to detect clinically significant errors in radiology reports (details in the Ap-



Figure 4: Comparisons with baselines on the MIMIC-CXR test set. "# Matched Findings" refers to the number of matched findings between generated and reference reports. "MCL-s" denotes our MCL model at a resolution of s.

pendix). Figure 4 shows that MCL surpasses R2Gen (Chen et al. 2020), CMN (Chen et al. 2021), and CGPT2 (Nicolson, Dowling, and Koopman 2023) in both "# Matched Findings" and GREEN score, demonstrating MCL's superior ability to generate clinically accurate reports.

Conclusion

In this paper, we have proposed MCL, a two-stage method for generating chest X-ray reports. We first introduced multiview enhanced contrastive learning to align multi-view radiographs with their corresponding report, addressing the challenge of handling varying numbers of views in existing methods. Subsequently, we fully exploited available indications for report generation. Experiments on MIMIC-CXR, MIMIC-ABN, Multi-view CXR, and Two-view CXR datasets suggest that incorporating multi-view radiographs, indications, and high-resolution radiographs improves the model's accuracy in describing clinical findings. In the future, we plan to utilize patient temporal data (Hou et al. 2023b) to capture disease progression and enhance model reliability through prediction uncertainty (Kang et al. 2024).

References

Bannur, S.; Bouzid, K.; Castro, D. C.; Schwaighofer, A.; Bond-Taylor, S.; Ilse, M.; Pérez-García, F.; Salvatelli, V.; Sharma, H.; Meissen, F.; Ranjit, M.; Srivastav, S.; Gong, J.; Falck, F.; Oktay, O.; Thieme, A.; Lungren, M. P.; Wetscherek, M. T.; Alvarez-Valle, J.; and Hyland, S. L. 2024. MAIRA-2: Grounded Radiology Report Generation. arXiv:2406.04449.

Beltagy, I.; Lo, K.; and Cohan, A. 2019. SciBERT: A Pretrained Language Model for Scientific Text. In *Findings of the Association for Computational Linguistics: EMNLP* 2019, 3615–3620.

Bu, S.; Li, T.; Yang, Y.; and Dai, Z. 2024. Instancelevel Expert Knowledge and Aggregate Discriminative Attention for Radiology Report Generation. In *Proceedings of the IEEE/CVF Conference on Computer Vision and Pattern Recognition*, 14194–14204.

Chen, T.; ZHANG, R.; and Hinton, G. 2023. Analog Bits: Generating Discrete Data using Diffusion Models with Self-Conditioning. In *The Eleventh International Conference on Learning Representations*.

Chen, X.; Fang, H.; Lin, T.-Y.; Vedantam, R.; Gupta, S.; Dollar, P.; and Zitnick, C. L. 2015. Microsoft COCO Captions: Data Collection and Evaluation Server. arXiv:1504.00325.

Chen, Z.; Li, G.; and Wan, X. 2022. Align, reason and learn: Enhancing medical vision-and-language pre-training with knowledge. In *Proceedings of the 30th ACM International Conference on Multimedia*, 5152–5161.

Chen, Z.; Shen, Y.; Song, Y.; and Wan, X. 2021. Crossmodal Memory Networks for Radiology Report Generation. In *Proceedings of the 59th Annual Meeting of the Association for Computational Linguistics and the 11th International Joint Conference on Natural Language Processing* (Volume 1: Long Papers), volume 1, 5904–5914.

Chen, Z.; Song, Y.; Chang, T.-H.; and Wan, X. 2020. Generating Radiology Reports via Memory-driven Transformer. In *Findings of the Association for Computational Linguistics: EMNLP 2020*, 1439–1449.

Cheng, P.; Lin, L.; Lyu, J.; Huang, Y.; Luo, W.; and Tang, X. 2023. PRIOR: Prototype Representation Joint Learning from Medical Images and Reports. In *Proceedings of the IEEE/CVF International Conference on Computer Vision*, 21361–21371.

Demner-Fushman, D.; Kohli, M. D.; Rosenman, M. B.; Shooshan, S. E.; Rodriguez, L.; Antani, S.; Thoma, G. R.; and McDonald, C. J. 2016. Preparing a collection of radiology examinations for distribution and retrieval. *Journal of the American Medical Informatics Association*, 23(2): 304–310.

He, K.; Zhang, X.; Ren, S.; and Sun, J. 2016. Deep residual learning for image recognition. In *Proceedings of the IEEE/CVF Conference on Computer Vision and Pattern Recognition*, 770–778.

Hou, W.; Cheng, Y.; Xu, K.; Li, W.; and Liu, J. 2023a. RE-CAP: Towards Precise Radiology Report Generation via Dynamic Disease Progression Reasoning. In *Findings of the* Association for Computational Linguistics: EMNLP 2023, 2134–2147.

Hou, W.; Xu, K.; Cheng, Y.; Li, W.; and Liu, J. 2023b. ORGAN: Observation-Guided Radiology Report Generation via Tree Reasoning. In *Proceedings of the 61st Annual Meeting of the Association for Computational Linguistics (Volume 1: Long Papers)*, 8108–8122.

Huang, Z.; Zhang, X.; and Zhang, S. 2023. KiUT: Knowledge-injected U-Transformer for Radiology Report Generation. In *Proceedings of the IEEE/CVF Conference on Computer Vision and Pattern Recognition*, 19809–19818.

Irvin, J.; Rajpurkar, P.; Ko, M.; Yu, Y.; Ciurea-Ilcus, S.; Chute, C.; Marklund, H.; Haghgoo, B.; Ball, R.; Shpanskaya, K.; Seekins, J.; Mong, D. A.; Halabi, S. S.; Sandberg, J. K.; Jones, R.; Larson, D. B.; Langlotz, C. P.; Patel, B. N.; Lungren, M. P.; and Ng, A. Y. 2019. Chexpert: A large chest radiograph dataset with uncertainty labels and expert comparison. In *Proceedings of the AAAI conference on artificial intelligence*, volume 33, 590–597.

Jain, S.; Agrawal, A.; Saporta, A.; Truong, S.; Duong, D. N. D. N.; Bui, T.; Chambon, P.; Zhang, Y.; Lungren, M.; Ng, A.; Langlotz, C.; Rajpurkar, P.; and Rajpurkar, P. 2021. Radgraph: Extracting clinical entities and relations from radiology reports. In *Thirty-fifth Conference on Neural Information Processing Systems Datasets and Benchmarks Track (Round 1)*, volume 1.

Jin, H.; Che, H.; Lin, Y.; and Chen, H. 2024. PromptMRG: Diagnosis-Driven Prompts for Medical Report Generation. In *Proceedings of the AAAI Conference on Artificial Intelligence*, volume 38, 2607–2615.

Johnson, A. E. W.; Pollard, T. J.; Greenbaum, N. R.; Lungren, M. P.; ying Deng, C.; Peng, Y.; Lu, Z.; Mark, R. G.; Berkowitz, S. J.; and Horng, S. 2019. MIMIC-CXR-JPG, a large publicly available database of labeled chest radiographs. arXiv:1901.07042.

Kang, X.; Ma, Z.; Liu, K.; Li, Y.; and Miao, Q. 2024. Multi-scale Information Sharing and Selection Network with Boundary Attention for Polyp Segmentation. arXiv:2405.11151.

Li, C. Y.; Liang, X.; Hu, Z.; and Xing, E. P. 2019. Knowledge-Driven Encode, Retrieve, Paraphrase for Medical Image Report Generation. In *Proceedings of the AAAI Conference on Artificial Intelligence*, volume 33, 6666– 6673.

Li, J.; Li, D.; Xiong, C.; and Hoi, S. 2022. BLIP: Bootstrapping Language-Image Pre-training for Unified Vision-Language Understanding and Generation. In *International Conference on Machine Learning*, volume 162, 12888– 12900.

Liu, A.; Guo, Y.; Yong, J.-h.; and Xu, F. 2024a. Multigrained Radiology Report Generation with Sentence-level Image-language Contrastive Learning. *IEEE Transactions on Medical Imaging*, 1–1.

Liu, C.; Tian, Y.; Chen, W.; Song, Y.; and Zhang, Y. 2024b. Bootstrapping Large Language Models for Radiology Report Generation. In *Proceedings of the AAAI Conference on Artificial Intelligence*, volume 38, 18635–18643. Liu, K.; Ma, Z.; Kang, X.; Zhong, Z.; Jiao, Z.; Baird, G.; Bai, H.; and Miao, Q. 2024c. Structural Entities Extraction and Patient Indications Incorporation for Chest X-Ray Report Generation. In *Medical Image Computing and Computer Assisted Intervention – MICCAI 2024*, 433–443. Cham: Springer Nature Switzerland. ISBN 978-3-031-72384-1.

Liu, K.; Ma, Z.; Liu, M.; Jiao, Z.; Kang, X.; Miao, Q.; and Xie, K. 2024d. Factual Serialization Enhancement: A Key Innovation for Chest X-ray Report Generation. arXiv:2405.09586.

Liu, R.; Li, M.; Zhao, S.; Chen, L.; Chang, X.; and Yao, L. 2024e. In-Context Learning for Zero-shot Medical Report Generation. In *ACM Multimedia 2024*.

Liu, Z.; Zhu, Z.; Zheng, S.; Zhao, Y.; He, K.; and Zhao, Y. 2024f. From Observation to Concept: A Flexible Multi-View Paradigm for Medical Report Generation. *IEEE Transactions on Multimedia*, 26: 5987–5995.

Messina, P.; Pino, P.; Parra, D.; Soto, A.; Besa, C.; Uribe, S.; Andía, M.; Tejos, C.; Prieto, C.; and Capurro, D. 2022. A survey on deep learning and explainability for automatic report generation from medical images. *ACM Computing Surveys (CSUR)*, 54(10s): 1–40.

Nguyen, D.; Chen, C.; He, H.; and Tan, C. 2023. Pragmatic Radiology Report Generation. In *Machine Learning for Health (ML4H)*, volume 225, 385–402.

Ni, J.; Hsu, C.; Gentili, A.; and McAuley, J. J. 2020. Learning Visual-Semantic Embeddings for Reporting Abnormal Findings on Chest X-rays. In *Findings of the Association for Computational Linguistics: EMNLP 2020*, 1954–1960.

Nicolson, A.; Dowling, J.; and Koopman, B. 2023. Improving chest X-ray report generation by leveraging warm starting. *Artificial Intelligence in Medicine*, 144: 102633.

Ostmeier, S.; Xu, J.; Chen, Z.; Varma, M.; Blankemeier, L.; Bluethgen, C.; Michalson, A. E.; Moseley, M.; Langlotz, C.; Chaudhari, A. S.; and Delbrouck, J.-B. 2024. GREEN: Generative Radiology Report Evaluation and Error Notation. arXiv:2405.03595.

Paul, A.; Shen, T. C.; Lee, S.; Balachandar, N.; Peng, Y.; Lu, Z.; and Summers, R. M. 2021. Generalized Zero-Shot Chest X-Ray Diagnosis Through Trait-Guided Multi-View Semantic Embedding With Self-Training. *IEEE Transactions on Medical Imaging*, 40(10): 2642–2655.

Shen, H.; Pei, M.; Liu, J.; and Tian, Z. 2024. Automatic Radiology Reports Generation via Memory Alignment Network. In *Proceedings of the AAAI Conference on Artificial Intelligence*, volume 38, 4776–4783.

Smit, A.; Jain, S.; Rajpurkar, P.; Pareek, A.; Ng, A.; and Lungren, M. 2020. Combining Automatic Labelers and Expert Annotations for Accurate Radiology Report Labeling Using BERT. In *Findings of the Association for Computational Linguistics: EMNLP 2020.*

Tian, J.; Zhong, C.; Shi, Z.; and Xu, F. 2019. Towards Automatic Diagnosis from Multi-modal Medical Data. In *Interpretability of Machine Intelligence in Medical Image Computing and Multimodal Learning for Clinical Decision Support*, volume 11797, 67–74. Tian, Y.; Fan, L.; Isola, P.; Chang, H.; and Krishnan, D. 2023. StableRep: Synthetic Images from Text-to-Image Models Make Strong Visual Representation Learners. In *Advances in Neural Information Processing Systems*, volume 36, 48382–48402.

Touvron, H.; Martin, L.; Stone, K.; et al. 2023. Llama 2: Open Foundation and Fine-Tuned Chat Models. arXiv:2307.09288.

van den Oord, A.; Li, Y.; and Vinyals, O. 2019. Representation Learning with Contrastive Predictive Coding. arXiv:1807.03748.

Vaswani, A.; Shazeer, N.; Parmar, N.; Uszkoreit, J.; Jones, L.; Gomez, A. N.; Kaiser, L. u.; and Polosukhin, I. 2017. Attention is All you Need. In *Advances in Neural Information Processing Systems*, volume 30.

Wang, F.; Du, S.; and Yu, L. 2024. HERGen: Elevating Radiology Report Generation with Longitudinal Data. arXiv:2407.15158.

Wang, F.; Zhou, Y.; Wang, S.; Vardhanabhuti, V.; and Yu, L. 2022a. Multi-granularity cross-modal alignment for generalized medical visual representation learning. In *Advances in Neural Information Processing Systems*, volume 35, 33536–33549.

Wang, Z.; Liu, L.; Wang, L.; and Zhou, L. 2023. ME-Transformer: Radiology Report Generation by Transformer with Multiple Learnable Expert Tokens. In *Proceedings of the IEEE/CVF Conference on Computer Vision and Pattern Recognition*, 11558–11567.

Wang, Z.; Wu, Z.; Agarwal, D.; and Sun, J. 2022b. Med-CLIP: Contrastive Learning from Unpaired Medical Images and Text. In *Findings of the Association for Computational Linguistics: EMNLP 2022*, 3876–3887.

Xu, C.; Tao, D.; and Xu, C. 2013. A Survey on Multi-view Learning. arXiv:1304.5634.

Yan, B.; Liu, R.; Kuo, D. E.; Adithan, S.; Reis, E. P.; Kwak, S.; Venugopal, V. K.; O'Connell, C.; Saenz, A.; Rajpurkar, P.; and Moor, M. 2023. Style-aware radiology report generation with radgraph and few-shot prompting. In *Findings of the Association for Computational Linguistics: EMNLP 2023*, 14676–14688.

Yang, S.; Wu, X.; Ge, S.; Zheng, Z.; Zhou, S. K.; and Xiao, L. 2023. Radiology report generation with a learned knowledge base and multi-modal alignment. *Medical Image Analysis*, 86: 102798.

Yu, F.; Endo, M.; Krishnan, R.; Pan, I.; Tsai, A.; Reis, E. P.; Fonseca, E. K. U. N.; Lee, H. M. H.; Abad, Z. S. H.; Ng, A. Y.; Langlotz, C. P.; Venugopal, V. K.; and Rajpurkar, P. 2023. Evaluating progress in automatic chest x-ray radiology report generation. *Patterns*, 4(9): 100802.

Zhao, J.; Xie, X.; Xu, X.; and Sun, S. 2017. Multi-view learning overview: Recent progress and new challenges. *Information Fusion*, 38: 43–54.

Zhou, H.-Y.; Chen, X.; Zhang, Y.; Luo, R.; Wang, L.; and Yu, Y. 2022. Generalized radiograph representation learning via cross-supervision between images and free-text radiology reports. *Nature Machine Intelligence*, 4(1): 32–40.

Supplementary Material

Experimental Settings

Implementation Details. 1) MIMIC-CXR: In stage 1, we train the model for 50 epochs using the AdamW optimizer with a learning rate of 5e-5. In stage 2, we use the RAdam optimizer for 50 epochs, with learning rates of 5e-6 for the pre-trained model and 5e-5 for other parameters. 2) MIMIC-ABN and Multi-view CXR: Since most radiographs in these datasets originate from MIMIC-CXR (Johnson et al. 2019), we skip stage 1 and directly fine-tune the pre-trained model from the MIMIC-CXR dataset using the RAdam optimizer for 50 epochs, with learning rates of 5e-6 for MIMIC-ABN (Ni et al. 2020) and 5e-7 for Multiview CXR. To ensure fairness, we fine-tuned the pre-trained models of comparable methods on the MIMIC-CXR dataset for performance evaluation. 3) Two-view CXR: We utilize the same optimizer, learning rate, and epoch settings as in MIMIC-CXR. Peer methods are implemented according to their official code on the MIMIC-CXR, utilizing concatenation for visual feature fusion. 4) Other settings: For MCL with a resolution of 224×224 , radiograph pre-processing follows the procedure used in R2Gen (Chen et al. 2020). For MCL at a resolution of 384×384 , radiographs are resized to 448×448 , randomly cropped to 384×384 , rotated by 5 degrees, and then normalized. Additionally, we exclude studies with empty reference reports or those lacking clinical significance (e.g., reference report contains only Portable AP upright chest film ____ at 09:31 is submitted.). Metrics for natural language generation (NLG) and clinical efficacy (CE) are calculated using the pycocoevalcap (Chen et al. 2015), radgraph (Jain et al. 2021), and f1chexbert (Smit et al. 2020) libraries.

Results and Analyses

Effect of Multi-view Enhanced Contrastive Learning. Key findings from Table 2 in the main text are summarized as follows: 1) MCL enhances performance across all metrics, especially in CE metrics compared to (h), demonstrating the effectiveness of incorporating factual serialization in reports for cross-modal alignment. 2) A comparison between MCL and (f) reveals that multi-positive contrastive loss (MPC) effectively improves the consistency of visual features across multi-view radiographs by maximizing the agreement between these views. Notably, visual representation learning using MPC loss alone, as demonstrated in (c), also yields satisfactory results, further validating its effectiveness. 3) MCL exhibits improvements of 1.6%, 2.1%, and 4.9% in the B-2, RG, and CX14 metrics, respectively, compared to (e), indicating the critical role of instance-wise alignment. 4) Although MCL exhibits a slight decrease in the B-2 metric compared to (d), it shows improvement in other metrics, especially CE metrics, highlighting the significance of token-wise alignment loss. As a result, these results demonstrate the positive contributions of our multiview enhanced contrastive learning to model performance.

Model Benefits from Indications and Multi-view Radiographs. Table 6 in the Appendix reveals performance variations on subsets of the MIMIC-CXR test set with and

without multi-view radiographs (indications). Due to the lack of single-view radiographs in the MIMIC-ABN, Multiview CXR, and Two-view CXR datasets, the benefits of multi-view radiographs cannot be evaluated. Therefore, the experiment is conducted only on the MIMIC-CXR dataset. Significant improvements in NLG and RG metrics are observed for studies with multi-view radiographs (indications). However, the CX5 and CX14 metrics show a notable decline. This may be because these metrics cover only 14 specific observations, as shown in Table 4, which may not capture all relevant details (Figure 3 in the main text and Figure 5 in the Appendix). To explore this, we calculate failure rates for these metrics, denoted as CX5-F and CX14-F. We define the the CX14 (or CX5) metric as invalid when B-4 > 0.1, RG > 0.2, and CX14 = 0 (or CX5 = 0). These thresholds ensure certain similarities between the generated and reference reports, as shown in Figure 5 in the Appendix. The high failure rates of CX5 and CX14 contribute to lower performance in studies with multi-view radiographs (indications). Despite these limitations, CX5 and CX14 remain valuable for interpreting diagnostic performance based on these 14 specific observations.

Clinical Accuracy of 14 Observations. Tables 4 and 5 in the Appendix present the clinical accuracy of 14 observations related to thoracic diseases and support devices. To ensure a fair comparison, MCL uses the same 224×224 resolution as the baseline (Chen et al. 2020). Each report is automatically labeled with 14 types of observations using the CheXpert labeler (Irvin et al. 2019). The No Findings observation is classified as either present or absent, while other observations include the additional categories of uncertain and not-mentioned. Following the default settings of the f1chexbert library (Smit et al. 2020), not-mentioned and absent are categorized as negative, with other categories classified as positive. Our findings show that: 1) MCL significantly improves most F1 compared to R2Gen (Chen et al. 2020). 2) However, performance for some rare observations, such as Lung Lesion, Pleural Other, and Fracture, remains poor, consistent with previous works (Chen et al. 2020; Liu et al. 2024a). This is primarily due to the limited number of training samples for these patterns.

Case Study. Figure 3 in the main text and Figure 5 in the Appendix present five examples from the MIMIC-CXR test set, depicting the reports generated by our model. Key findings are as follows: 1) Our proposed MCL consistently outperforms the baseline in these cases, delivering effective draft reports for radiologists to review. 2) Although the CX5 and CX14 metrics may not capture all relevant details, they provide valuable interpretability by assessing diagnostic performance based on 14 specific and common observations (see Tables 4 and 5 in the Appendix).

Evaluation via Large Language Models. The GREEN score (Ostmeier et al. 2024) identifies six types of clinically significant errors: (a) False report of a finding in the candidate; (b) Missing a finding present in the reference; (c) Misidentification of a finding's anatomic location/position; (d) Misassessment of the severity of a finding; (e) Mentioning a comparison that isn't in the reference; (f) Omitting a comparison detailing a change from a prior study. Addi-

Snlit	MIMIC-CXR			MIMIC-ABN			Multi-view CXR			Two-view CXR			
Spiit	# Img	# Rpt	% Ind	# Img	# Rpt	% Ind	# Img	# Rpt	% Ind	# Img	# Rpt	% Ind	
Train	269,239	150,957	66.4%	69,526	34,763	64.6%	220,978	100,505	67.1%	165,056	82,528	67.9%	
Val	2,113	1,182	65.4%	526	263	62.7%	2,299	1,057	71.2%	1,778	889	72.1%	
Test	3,852	2,343	57.3%	756	378	56.3%	3,947	1,805	67.6%	3,000	1,500	68.9%	

Table 3: Statistics for the training, validation, and test sets across the four datasets. "# Img" and "# Rpt" denote the number of radiographs and reports, respectively. "% Ind" is the ratio of indications.

		Two-view CXR												
Observation	$\mathcal{O}_{\alpha}(n)$	R2Gen			M	CL (Ou	ırs)	$\mathcal{O}_{\alpha}(n)$		R2Gen	l	M	CL (Ou	rs)
	-70 (n)	Р	R	\mathbf{F}_1	P	R	\mathbf{F}_1	<i>// (ii)</i>	Р	R	\mathbf{F}_1	P	R	\mathbf{F}_1
ECM	10.2% (1,096)	0.377	0.353	0.365	0.396	0.345	0.369	9.8% (293)	0.371	0.280	0.319	0.351	0.300	0.324
Cardiomegaly	14.7% (1,573)	0.567	0.347	0.431	0.628	0.566	0.595	14.1% (423)	0.590	0.409	0.483	0.578	0.567	0.573
Lung Opacity	13.3% (1,422)	0.457	0.231	0.307	0.578	0.352	0.437	13.1% (391)	0.614	0.069	0.124	0.567	0.335	0.421
Lung Lesion	2.6% (281)	0.481	0.046	0.084	0.217	0.036	0.061	2.8% (84)	0.500	0.012	0.023	0.357	0.060	0.102
Edema	8.3% (882)	0.445	0.254	0.323	0.514	0.447	0.478	5.8% (174)	0.523	0.259	0.346	0.475	0.333	0.392
Consolidation	3.3% (356)	0.188	0.096	0.127	0.240	0.154	0.188	3.5% (104)	0.300	0.115	0.167	0.279	0.163	0.206
Pneumonia	4.6% (490)	0.238	0.204	0.220	0.264	0.186	0.218	3.7% (112)	0.276	0.071	0.113	0.349	0.268	0.303
Atelectasis	11.1% (1,184)	0.421	0.249	0.313	0.496	0.519	0.507	10.2% (306)	0.463	0.248	0.323	0.483	0.464	0.473
Pneumothorax	1% (102)	0.333	0.029	0.054	0.291	0.157	0.204	0.8% (23)	0.000	0.000	0.000	0.222	0.087	0.125
Pleural Effusion	12.5% (1,340)	0.804	0.196	0.316	0.716	0.669	0.691	10.5% (315)	0.787	0.375	0.508	0.709	0.667	0.687
Pleural Other	1.7% (181)	0.000	0.000	0.000	0.217	0.083	0.120	1.9% (56)	0.500	0.018	0.034	0.091	0.036	0.051
Fracture	2% (215)	0.000	0.000	0.000	0.071	0.014	0.023	2.3% (70)	0.000	0.000	0.000	0.188	0.043	0.070
Support Devices	12.5% (1,333)	0.723	0.489	0.583	0.755	0.713	0.734	9.5% (285)	0.706	0.488	0.577	0.722	0.646	0.681
No Finding	2.2% (238)	0.130	0.626	0.215	0.223	0.660	0.333	12.0% (358)	0.387	0.950	0.550	0.547	0.813	0.654
micro avg	-	0.432	0.280	0.340	0.538	0.465	0.499	-	0.484	0.341	0.400	0.539	0.469	0.501
macro avg	-	0.369	0.223	0.238	0.400	0.350	0.354	-	0.430	0.235	0.255	0.423	0.342	0.362

Table 4: Clinical accuracy of our proposed MCL on the MIMIC-CXR and Two-view CXR datasets. ECM stands for *Enlarged Cardiomediastinum*. P, R, and F₁ denote the precision, recall, and F₁-score, respectively.

	MIMIC-ABN							Multi-view CXR						
Observation	0/2 (n)	R2Gen				CL (Ou	rs)	Ø (n)	R2Gen				CL (Ou	rs)
	<i>70</i> (<i>n</i>)	Р	R	\mathbf{F}_1	Р	R	F1	<i>70</i> (<i>n</i>)	Р	R	F1	Р	R	\mathbf{F}_1
ECM	5.4 (94)	0.197	0.149	0.170	0.205	0.181	0.192	10.7 (884)	0.342	0.189	0.243	0.403	0.369	0.385
Cardiomegaly	12.5 (220)	0.523	0.723	0.607	0.571	0.636	0.602	14.3 (1,179)	0.629	0.415	0.500	0.639	0.465	0.538
Lung Opacity	20.2 (354)	0.596	0.463	0.521	0.679	0.514	0.585	13.1 (1,079)	0.564	0.073	0.130	0.608	0.320	0.419
Lung Lesion	4.9 (86)	0.000	0.000	0.000	0.500	0.093	0.157	3.1 (258)	0.400	0.008	0.015	0.154	0.023	0.040
Edema	7.1 (124)	0.417	0.242	0.306	0.470	0.444	0.456	6.4 (530)	0.540	0.292	0.379	0.459	0.385	0.419
Consolidation	3.4 (60)	0.255	0.233	0.243	0.339	0.350	0.344	3.2 (267)	0.198	0.097	0.131	0.292	0.169	0.214
Pneumonia	6.2 (108)	0.254	0.139	0.180	0.300	0.278	0.288	4.1 (339)	0.286	0.083	0.128	0.288	0.221	0.250
Atelectasis	10.7 (188)	0.341	0.633	0.443	0.406	0.527	0.458	10.1 (833)	0.444	0.202	0.277	0.480	0.449	0.464
Pneumothorax	0.0 (0)	0.000	0.000	0.000	0.000	0.000	0.000	0.7 (55)	0.263	0.091	0.135	0.182	0.109	0.136
Pleural Effusion	9.2 (162)	0.660	0.623	0.641	0.598	0.735	0.659	10.5 (861)	0.770	0.469	0.583	0.719	0.631	0.672
Pleural Other	3.1 (54)	0.250	0.019	0.034	0.350	0.130	0.189	2.1 (170)	0.333	0.006	0.012	0.179	0.071	0.101
Fracture	3.9 (68)	0.000	0.000	0.000	0.000	0.000	0.000	2.6 (211)	0.000	0.000	0.000	0.182	0.019	0.034
Support Devices	10.6 (186)	0.748	0.608	0.671	0.691	0.699	0.695	9.1 (747)	0.657	0.522	0.582	0.680	0.632	0.655
Ño Finding	2.9 (50)	0.150	0.460	0.227	0.135	0.300	0.186	10.0 (823)	0.329	0.968	0.491	0.473	0.943	0.630
micro avg	-	0.456	0.429	0.442	0.496	0.469	0.482	-	0.462	0.329	0.384	0.525	0.454	0.487
macro avg	-	0.314	0.307	0.289	0.375	0.349	0.344	-	0.411	0.244	0.258	0.410	0.343	0.354

Table 5: Clinical accuracy of our proposed MCL on the MIMIC-ABN and Multi-view CXR datasets. ECM stands for *Enlarged Cardiomediastinum*. P, R, and F₁ denote the precision, recall, and F₁-score, respectively.

B-4† RG† CX14† 0.491/0.125 0.732/0.310 0.857/0.333	Reference report: lung volumes are low	Generated report by R2Gen: frontal and	Generated report by MCL (Ours): lung
	compared to the previous study . left-sided aicd	lateral views of the chest were obtained . dual-	volumes are low . left - sided aicd device is
	device is noted with single lead terminating in	lead left-sided pacemaker is again seen with	noted with single lead terminating in the
	unchanged position in the right ventricle . heart	leads extending to the expected positions of	right ventricle . heart size is mildly enlarged .
	size appears at least mildly enlarged . the	the right atrium and right ventricle . there are	mediastinal contours are unremarkable .
	mediastinal and hilar contours are unremarkable .	low lung volumes . there are relatively low	there is crowding of the bronchovascular
	there is crowding of the bronchovascular structures	lung volumes . there is elevation of the right	structures with likely mild pulmonary
	without overt pulmonary edema . streaky opacities	hemidiaphragm and overlying atelectasis . no	vascular congestion . patchy bibasilar
	in the lung bases likely reflect areas of atelectasis .	definite focal consolidation is seen . there is no	airspace opacities likely reflect atelectasis .
	no pleural effusion or pneumothorax is present .	pleural effusion or pneumothorax . the cardiac	no pleural effusion or pneumothorax is seen .
	there are no acute osseous abnormalities .	and mediastinal silhouettes are stable .	there are no acute osseous abnormalities .
B-41 RG1 CX141 0.591/0.245 0.919/0.421 0.000/0.000	Reference report: frontal and lateral views of	Generated report by R2Gen: pa and lateral	Generated report by MCL (Ours): frontal
	the chest were obtained . there are relatively low	views of the chest provided . lung volumes are	and lateral views of the chest were obtained .
	lung volumes . mild elevation of the right	low limiting assessment . allowing for this	there are relatively low lung volumes . there
	hemidiaphragm persists . there is persistent right	there is no focal consolidation effusion or	is mild elevation of the right hemidiaphragm .
	base atelectasis . no new focal consolidation is	pneumothorax . the cardiomediastinal	no focal consolidation , pleural effusion , or
	seen . there is no pleural effusion or	silhouette is normal . imaged osseous	evidence of pneumothorax is seen . the
	pneumothorax . the cardiac and mediastinal	structures are intact . no free air below the	cardiac and mediastinal silhouettes are
	silhouettes are unremarkable .	right hemidiaphragm is seen .	unremarkable .
B-41 RG1 CX141 0.100/0.001 0.235/0.074 0.000/0.000	Reference report: the right loculated pleural effusion has decreased in size from there is a little to no remaining effusion on the left and there has been interval resolution of the mild interstitial pulmonary edema . there are no focal consolidations to suggest active infectious process . there is no pneumothorax . the hilar and cardiomediastinal contours are normal . the hemidiaphragms are flattened and the ap diameter is increased .	Generated report by R2Gen: pa and lateral chest views were obtained with patient in upright position . analysis is performed in direct comparison with the next preceding similar study of . the heart size remains unchanged and is within normal limits . unchanged appearance of thoracic aorta . the pulmonary vasculature is not congested . there is now a right-sided port-a-cath system seen to terminate in the mid portion of the svc . no pneumothorax is identified in the apical area . similar as on the preceding examination there is a diffuse interstitial pattern in both lungs with a peripheral distribution in the periphery	Generated report by MCL (Ours): frontal and lateral chest radiographs demonstrate a right chest port with the tip terminating in the low svc . there is a small right apical pneumothorax without evidence of tension . post - surgical changes are noted in the right lower lung . there is no focal consolidation or left pleural effusion . the cardiomediastinal silhouette is normal .

Figure 5: Examples of generated reports on the MIMIC-CXR test set with a resolution of 224×224 . The cell "A/B" represents "MCL/R2Gen". Sentences in the reference report are shown in different colors. Each sentence in generated reports is highlighted in matching colors corresponding to those in the reference report.

Reference Report: the cardiac, mediastinal, and hilar contours are normal. pulmonary vascularity is normal, and the lungs are clear. no pleural effusion or pneumothorax is present. there are no acute osseous abnormalities. clips are noted within the upper abdomen compatible with prior cholecystectomy. Generated Report by MCL (Ours): the cardiac, mediastinal and hilar contours are normal. the lungs are clear and the pulmonary vascularity is normal. no pleural effusion or pneumothorax is present. there are no acute osseous abnormalities. several clips are noted within the left upper quadrant of the abdomen.

[Explanation]:

GREEN Response The candidate report is largely accurate, but there is a misidentification of the anatomic location of the clips. The reference report mentions that the clips are in the upper abdomen, while the candidate report specifies them as being in the left upper quadrant of the abdomen. This could potentially lead to confusion or misinterpretation in a clinical setting.

[Clinically Significant Errors]:

(a) False report of a finding in the candidate: 0.

(b) Missing a finding present in the reference: 0.

(c) Misidentification of a finding's anatomic location/position: 1. The candidate report incorrectly identifies the location of the clips as being in the left upper quadrant of the abdomen, while the reference report mentions them as being in the upper abdomen.

(d) Misassessment of the severity of a finding: 0.

(e) Mentioning a comparison that isn't in the reference: 0.

(f) Omitting a comparison detailing a change from a prior study: 0.

[Clinically Insignificant Errors]:

(a) False report of a finding in the candidate: 0.

(b) Missing a finding present in the reference: 0.

(c) Misidentification of a finding's anatomic location/position: 0.

(d) Misassessment of the severity of a finding: 0.

(e) Mentioning a comparison that isn't in the reference: 0.

(f) Omitting a comparison detailing a change from a prior study: 0.

[Matched Findings]:

5. The cardiac, mediastinal and hilar contours are normal; The lungs are clear; The pulmonary vascularity is normal; No pleural effusion or pneumothorax is present; There are no acute osseous abnormalities.

Figure 6: Clinically significant errors and matched findings in a generated report.

Metric	w/ MV	w/o MV	w/ Ind	w/o Ind
%	70.7%	29.3%	57.8%	42.2%
п	2,724	1,128	2,225	1,627
B-1	0.406	0.369	0.413	0.372
B-2	0.270	0.242	0.280	0.238
B-3	0.196	0.176	0.206	0.169
B-4	0.151	0.136	0.161	0.128
MTR	0.170	0.157	0.174	0.157
R-L	0.316	0.298	0.321	0.298
RG	0.287	0.249	0.299	0.244
CX5	0.547	0.573	0.548	0.566
CX14	0.474	0.545	0.490	0.509
CX5-F	24.0%	12.5%	24.0%	16.0%
CX14-F	24.0%	12.5%	24.0%	16.0%

Table 6: Breakdown of MCL's metrics on the MIMIC-CXR test set using a resolution of 224×224 , categorized by: (i) the presence of multi-view radiographs (MV), and (ii) the inclusion of an indication (Ind), in all studies. CX5-F and CX14-F denote the failure rates of the CX5 and CX14 metrics, respectively.

tionally, GREEN calculates the number of matched findings between generated and reference reports, denoted as "# Matched Findings", with higher values indicating better alignment. The GREEN (Ostmeier et al. 2024) score is formulated as:

 $GREEN = \frac{\# \text{ Matched Findings}}{\# \text{ Matched Findings} + \sum_{i=(a)}^{(f)} \# \text{ Error}_{\text{sig.}, i}},$ (10)

where $\# \operatorname{Error}_{\operatorname{sig.}, i}$ refers to the count of the i^{th} clinically significant error. A higher GREEN score indicates better report quality. Clinically significant errors and matched findings in a generated report are presented in Figure 6 of the Appendix.