EyecareGPT: Boosting Comprehensive Ophthalmology Understanding with Tailored Dataset, Benchmark and Model

Sijing Li¹, Tianwei Lin¹, Lingshuai Lin², Wenqiao Zhang¹, Jiang Liu¹, Xiaoda Yang¹, Juncheng Li¹, Yucheng He³, Xiaohui Song¹, Jun Xiao¹, Yueting Zhuang¹, Beng Chin Ooi⁴

¹Zhejiang University, ²Harbin Institute of Technology, ³The First People's Hospital of Chenzhou, ⁴National University of Singapore



Code

Figure 1: Overview of the Eyecare-100K. Eyecare-100K aggregates real-world ophthalmic data across 8 modalities, 15+ anatomical structures and 100+ eye diseases, supporting multi-modal report generation and fine-grained visual QA tasks.

Abstract

Medical Large Vision-Language Models (Med-LVLMs) demonstrate significant potential in healthcare, but their reliance on general medical data and coarse-grained global visual understanding limits them in intelligent ophthalmic diagnosis. Currently, intelligent ophthalmic diagnosis faces three major challenges: (i) Data. The lack of deeply annotated, high-quality, multi-modal ophthalmic visual instruction data; (ii) Benchmark. The absence of a comprehensive and systematic benchmark for evaluating diagnostic performance; (iii) Model. The difficulty of

adapting holistic visual architectures to fine-grained, region-specific ophthalmic lesion identification. In this paper, we propose the **Eyecare Kit**, which systematically tackles the aforementioned three key challenges with the tailored dataset, benchmark and model: First, we construct a multi-agent data engine with real-life ophthalmology data to produce **Eyecare-100K**, a high-quality ophthalmic visual instruction dataset. Subsequently, we design **Eyecare-Bench**, a benchmark that comprehensively evaluates the overall performance of LVLMs on intelligent ophthalmic diagnosis tasks across multiple dimensions. Finally, we develop the **EyecareGPT**, optimized for fine-grained ophthalmic visual understanding thoroughly, which incorporates an adaptive resolution mechanism and a layer-wise dense connector. Extensive experimental results indicate that the EyecareGPT achieves state-of-the-art performance in a range of ophthalmic tasks, underscoring its significant potential for the advancement of open research in intelligent ophthalmic diagnosis. Our project is available at https://github.com/DCDmllm/EyecareGPT.

1 Introduction

Large Vision-Language Models (LVLMs) (Liu et al. 2024b,a; Lin et al. 2024; Team et al. 2024; Hurst et al. 2024) achieve remarkable progress in open-world visual understanding tasks (Ren et al. 2024; Zhang et al. 2019b), demonstrating potential in medical scenarios. In recent years, Medical Large Vision-Language Models (Med-LVLMs), such as LLaVA-Med (Li et al. 2023), HuatuoGPT-Vision(Chen et al. 2024b), and HealthGPT (Lin et al. 2025) trained on extensive medical visual instruction data, advanced medical tasks including pathological diagnosis and knowledge reasoning. However, due to the lack of deep modeling of modality-specific features and domain-specific expertise in vertical medical fields, existing Med-LVLMs still exhibit significant limitations in fine-grained visual understanding and region-specific intelligent diagnosis within specific disciplines. Taking ophthalmology as an example, the field involves complex medical imaging modalities and highly specialized clinical requirements (Balas et al. 2024; Xu et al. 2024b), where current ophthalmic foundation models (Shi et al. 2024) and Med-LVLMs (Saab et al. 2024; AlSaad et al. 2024; Lin et al. 2025) fail to provide effective support. Therefore, developing advanced Med-LVLM specifically tailored for ophthalmology, equipped with fine-grained visual understanding and reasoning capabilities, becomes an urgent need to drive intelligent ophthalmology research and applications.

Effectively transferring the Med-LVLM paradigm to ophthalmology requires a systematic analysis of the domain's unique requirements and challenges. First, existing ophthalmic datasets primarily focus on single imaging modality and pathological classification task, lacking deeply annotated, high-quality, heterogeneous multi-modal visual instruction data that cover temporal sequence data (e.g., FFA, ICGA), spatial sequence data (e.g., UBM, CT), and complex modalities (e.g., Fundus, OCT, Slit-Lamp). Second, the absence of a comprehensive benchmark for intelligent ophthalmic diagnosis hinders accurate evaluation of Med-LVLMs on fine-grained visual understanding and report generation tasks, and limits guidance for model optimization. Finally, current Med-LVLM architectures, relying on coarse-grained global features, often overlook critical fine-grained priors and local details in ophthalmic imaging, failing to meet the precision requirements of intelligent diagnosis. To address these challenges, we propose the **Eyecare Kit**, which systematically advances the adaptability and performance of Med-LVLMs in ophthalmology through innovations in three key aspects: **Dataset**, **Benchmark**, and **Model**.

(i) Dataset. To address the scale, modality, and task diversity of ophthalmic data, we collect real-world data from 13 public datasets, 3 hospitals, and 3 public medical case libraries, covering 8 imaging modalities, over 15 anatomical structures, and more than 100 types of eye diseases (Figure 1). A multi-agent data engine is developed for information extraction, cleaning, formatting, and expert review, resulting in a comprehensive dataset named Eyecare-100K, containing approximately 102,000 visual question answering (VQA) pairs. Eyecare-100K is the first comprehensive ophthalmic dataset to simultaneously encompass multiple modalities, tasks, and diseases, and is expected to serve as a key resource for advancing multi-modal intelligent understanding in ophthalmology.

(ii) Benchmark. To deeply evaluate the comprehensive ophthalmology understanding capability of a Med-LVLM, we develop a benchmark named Eyecare-Bench, which includes three key clinical metrics: closed QA, open QA, and report generation. The data instances in Eyecare-Bench are drawn from the designated test set of Eyecare-100K, comprising about 15,000 carefully sampled examples across tasks, modalities, and disease categories to ensure balanced and representative evaluation. Notably, we design multidimensional evaluation metrics for different tasks and introduce GPT-4 to provide a more comprehensive evaluation of the report generation capabilities of models. Eyecare-Bench provides significant reference value for the open community to research more accurate and reliable intelligent systems for eye diseases.

(iii) Model. To address the fine-grained and regionspecific demands of intelligent ophthalmic diagnosis and to validate the effectiveness of the Eyecare Kit in supporting high-quality data and systematic evaluation, we propose the EyecareGPT model. We employ SigLIP (Zhai et al. 2023) as a high-resolution visual feature extractor to enhance local lesion perception. To accommodate variable resolutions in clinical ophthalmic imaging, we design an adaptive resolution mechanism (You et al. 2024; Guo and Huang 2025; Zhang et al. 2022) for dynamic adjustment, improving consistency across multi-modal images. Additionally, we introduce a Layer-wise Dense Connector (LDC) to densely integrate multi-scale visual features and preserve fine-grained structural information. Finally, we provide two scalable EyecareGPT variants to enable flexible deployment in diverse real-world settings.

Experimental results show that the Eyecare Kit provides high-quality, deeply annotated data for intelligent ophthalmic diagnosis and establishes a comprehensive evaluation benchmark, effectively supporting the optimization and development of Med-LVLMs. The main contributions of this work are as follows:

- **High-Quality Dataset.** We propose the first comprehensive ophthalmic visual instruction dataset.
- **Comprehensive Benchmark.** We build a systematic benchmark to evaluate the clinical performance of LVLMs on 3 core tasks: closed QA, open QA and report generation.
- Adapted LVLM Architecture. We introduce an LVLM architecture adapted to complex, heterogeneous oph-thalmic clinical imaging, achieving SoTA performance.
- Facilitating Open Research. We will fully open-source the dataset, benchmark, and model to facilitate research on intelligent ophthalmic diagnosis.

2 Related Work

2.1 Medical Large-Vision Language Models

Med-LVLMs achieve groundbreaking progress in processing and understanding medical imaging, offering new possibilities for clinical diagnosis and treatment (Chen et al. 2024a; Xu et al. 2024a). Med-Flamingo (Moor et al. 2023) leverages multimodal knowledge across medical disciplines for pre-training, extending the Flamingo (Alayrac et al. 2022) framework into the medical domain. Models such as LLaVA-Med (Li et al. 2023) and UMIT (Yu et al. 2025) adopt a two-stage training strategy combining pre-training and fine-tuning, enhancing vision-text alignment and multitask adaptability. To address language adaptability and dataset specificity challenges, HuatuoGPT-Vision (Chen et al. 2024b) introduces the PubMedVision dataset, comprising 1.3 million high-quality medical samples and markedly improving model adaptability. Specialized LVLMs like LLava-Rad (Zambrano Chaves et al. 2025) focus on radiology image understanding, actively exploring report generation tasks aligned with clinical practice. However, existing specialized ophthalmic models exhibit limited generalization; for instance, DeepDR-LLM (Li et al. 2024) is restricted to auxiliary diagnosis of diabetic retinopathy, and Ophtha-LLaMA2 (Zhao et al. 2023) uses only three ophthalmic modalities for fine-tuning. Currently, Med-LVLMs are gradually evolving from general medical tasks toward clinical practicality, yet the scarcity of high-quality datasets in specific subfields continues to hinder their development.

2.2 Ophthalmic Multi-Modal Datasets

High-quality ophthalmic datasets hold significant clinical and societal value in advancing intelligent diagnostic models. Currently, publicly available datasets primarily consist of ocular images labeled with classification tags or diagnostic keywords. For ocular disease classification, typical fundus image datasets include IDRID (Porwal et al. 2018), ACRIMA (Ovreiu, Paraschiv, and Ovreiu 2021), JSIEC (Cen et al. 2021), ODIR (on Ocular Disease Intelligent Recognition 2019), Harvard-GDP (Luo et al. 2023), MuRed (Rodríguez, AlMarzouqi, and Liatsis 2022), and DeepDRiD (Liu et al. 2022). Similarly, OCT2017 (Kermany, Zhang, and Goldbaum 2018a), Kermany (Kermany, Zhang, and Goldbaum 2018b), OCTID (Gholami et al. 2020), and OCTDL (Kulyabin et al. 2024) provide ocular OCT images with corresponding classification labels. Although these datasets contribute significantly to their respective tasks, they exhibit notable limitations in imaging modality diversity and fine-grained annotation, restricting their use in more complex intelligent diagnostic applications. The recently proposed multimodal ophthalmic dataset LMOD (Qin et al. 2025) covers five imaging modalities, partially addressing the limitations of single-modality datasets. However, the lack of medical visual instruction datasets for LVLM training highlights the need to develop larger, more diverse, and finely annotated heterogeneous multimodal ophthalmic datasets to support model training and evaluation.

3 Eyecare Kit: Eyecare-100K

3.1 Data Collation and Organization

Existing ophthalmic datasets are typically limited to a single task or modality and suffer from inconsistent standards and uneven distributions. These limitations make them unsuitable for constructing high-quality visual instruction datasets, thereby restricting the development of Med-LVLMs in intelligent ophthalmic applications. To address this, we propose **Eyecare-100K**, a comprehensive ophthalmic visual instruction dataset that covers diverse heterogeneous multi-modal imaging, aiming to provide standardized data support for intelligent understanding in ophthalmology.

To address the lack of available data for key clinical modalities, we collaborate with three large public hospitals to collect and annotate real-world ophthalmic cases with anonymization procedures. We also systematically collect and standardize multiple public ophthalmic datasets across different modalities and labels (see Appendix for details). To further expand diversity, Eyecare-100K incorporates examination cases from public medical repositories such as Radiopaedia (Knipe 2005) and MedPix (Siragusa et al. 2024), as well as professional slit-lamp textbooks (Liang and Zhang 2022). In total, Eyecare-100K integrates 58,485 ophthalmic images from 13 public datasets, 3 hospitals, and 3 medical case repositories. The dataset covers 8 imaging modalities (see Fig. 2): (1) Fluorescein Angiography (FA), (2) Indocyanine Green Angiography (ICGA), (3) Optical Coherence Tomography (OCT), (4) Fundus Photography, (5) Ultrasound Biomicroscopy (UBM), (6) Slit-Lamp, (7) Fluorescein Staining Imaging, and (8) Computed Tomography (CT), spanning 15 anatomical structures and over 100 ophthalmic diseases and rare conditions, significantly enhancing dataset diversity and comprehensiveness.

Considering that traditional medical data primarily consist of classification labels, segmentation annotations,



Figure 2: Data statistics of Eyecare-100K.

or brief textual descriptions and lack the visual instruction structures needed for fine-tuning Med-LVLMs, we develop a multi-agent data engine to extract, clean, standardize, and perform expert review on large-scale raw data. Ultimately, the data are organized into three types of VQA tasks: closed QA (multiple-choice questions), open QA (shortform questions), and report generation (long-text answers), to equip the models with fine-grained ophthalmic visual understanding and reasoning capabilities.

3.2 Multi-Agent Data Engine

We develop an automated **multi-agent data engine** to create Eyecare-100K, converting categorized labels and raw reports into structured VQA pairs. As shown in Figure 3, the agent engine comprises 6 components as follows.

Analyzer for Description Extraction. Given that a large volume of raw clinical data is stored in PDF format, containing imaging modalities, diagnostic results, and other details. Qwen2-VL-2B-OCR (Wang et al. 2024) is adapted as our analyzer to automatically extract key information and anonymize sensitive information.

Intelligent Collector for Medical Cases. We construct an intelligent collector to extract metadata from authorized public medical case repositories. The extracted metadata includes imaging modalities, anatomical structures, and diagnostic descriptions, enriching the diversity of data sources and expanding the coverage of clinical cases.

Translator for Data Sanitization. To address the common issues of complex abbreviations and mixed-language expressions in medical descriptions, we integrate Gemini-2.0-Flash (Google AI) as an automated translator. This module accurately interprets the contextual meaning of medical terms and converts raw descriptions into clear, standardized professional language, thereby enhancing the consistency of the dataset.

Template Library of QA and Prompt. To support finegrained tasks (closed QA, open QA, report generation), we designed a diverse VQA and prompt template library (see Appendix). For single-label data, we apply VQA templates to create open tasks involving modality identification, disease diagnosis, and lesion grading. For data containing detailed pathological information, we use prompt templates



Figure 3: Framework of Multi-Agent Data Engine.

to generate all three types of tasks. This design facilitates model training across multiple dimensions—understanding, reasoning, and generation—enhancing overall performance in real-world clinical scenarios.

Rewriter for Generating VQAs. Claude-3.7 serves as the rewriter, using prompts to extract key information from processed descriptions and construct reliable VQA pairs with the extracted information as answers. It supports generating VQA types such as modality identification, disease diagnosis and lesion localization, etc. For report generation prompts, it automatically organizes three components—examination types, imaging findings, and diagnostic recommendations—from the processed descriptions and generates a standardized Markdown format report.

Human Preference-Based Reviewer. To improve the accuracy and standardization of automatically generated data, we randomly select 10% constructed VQA instructions and report samples, and introduce five human reviewers to inspect the data quality. Medical professionals concentrate on semantic validity, medical consistency, and format standardization. Each data entry undergoes two rounds of independent review, effectively ensuring the fine-grained data quality of Eyecare-100K.

4 Eyecare Kit: Eyecare-Bench

We propose **Eyecare-Bench** to systematically evaluate the performance of Med-LVLMs in intelligent ophthalmic diagnosis. The data instances in Eyecare-Bench are drawn from the designated test set of Eyecare-100K, comprising about 15,000 examples. Sampling is carefully designed across all task types, imaging modalities, and ophthalmic disease categories to ensure representative proportions within each class, maintaining the balance and comprehensiveness of the test set. To the best of our knowledge, Eyecare-Bench is the most comprehensive benchmark to date for evaluating LVLMs in ophthalmology.

Table 1: Ten-Criteria evaluation framework.

Definition of Indicators	Weight
Scoring Indicators	
A) The number of abnormal features in candidate report that are not mentioned in the reference report.	1
B) The number of times the candidate report describes the disease severity incorrectly.	4
C) The number of times the candidate report describes the disease location incorrectly.	4
D) The number of missing key findings compared to the reference report.	6
E) Whether the diagnosis or suspected diagnosis is included.	2
F) Whether the description of the examination type exists and is correct.	2
G) Whether there is a treatment recommendation.	2
H) Whether the report structure is clear.	5
I) Whether the candidate outcome contains particularly serious clinical errors.	15
Correct Rate Calculation Indicator	
J) Whether the diagnosis is similar or approximately correct.	-

4.1 Multi-Dimensional Evaluation Suite

To systematically evaluate model performance on multitask and multi-format ophthalmic VQA tasks, we design a multi-dimensional evaluation suite, **EyeEval**. EyeEval defines fine-grained evaluation metrics for three task types, comprehensively covering aspects such as generation quality, factual consistency, and linguistic faithfulness.

VQA Evaluation Metrics In closed QA tasks, questions are presented in a multiple-choice format, aiming to assess the response accuracy of models. Therefore, we adopt accuracy as the primary evaluation metric. For open QA tasks, we focus on the factual consistency and linguistic coherence of generated answers. Specifically, we use F1-RadGraph (Yu et al. 2023), BERTScore-F1 (Zhang et al. 2019a), and F1-BioBert (Lee et al. 2020) to evaluate factual accuracy, BLEU (Papineni et al. 2002) and ROUGE-L (Lin and Hovy 2003) to assess surface similarity and language faithfulness.

Report Generation Evaluation Metrics. We recognize that traditional evaluation methods based on n-grams (Culy and Riehemann 2003) or semantic embeddings (Bakarov 2018) often introduce bias due to the inherent diversity and uncertainty of long-form text generation and correlate poorly with expert assessments. Therefore, in addition to the commonly used NLP metrics, we find five authoritative experts to develop a Ten-criteria evaluation framework (see Table 1) covering four key dimensions: **accuracy**, **completeness**, **structural rationality**, and **clinical practicability**. Each indicator is assigned a weight according to its clinical importance, and the total score of report is capped at 100. Evaluations are conducted using GPT-4 based on this refined rubric.

Specifically, indicators (A)–(D) are quantitative metrics assessing the accuracy and completeness of abnormal find-

ings in the report. The next five are Boolean indicators: (E)-(H) evaluate the structural coherence of the report, and (I) assesses the presence of critical errors that could affect clinical decision-making. Indicator (J) independently evaluates the diagnostic accuracy of the report. For quantitative indicators, deductions are applied proportionally based on the number of errors and corresponding weights. For Boolean indicators, points are deducted if the condition is not met.

According to the scoring criteria, we grade the reports as follows: (i) Excellent Report (90-100). (ii) Usable Report (80-90). (iii) Report Under Review (60-80). (iv) Unusable Report (Below 60). In the subsequent experiments, we validate the consistency between this evaluation framework and human expert assessments, thereby demonstrating the reliability of the framework and clinical relevance in quantitatively measuring report quality.

5 Eyecare Kit: EyecareGPT

5.1 Large Vision-Language Models

The input of LVLMs typically consists of an image $\boldsymbol{x}^{\text{img}} \in \mathbb{R}^{c \times h \times w}$ and a discrete text sequence $\boldsymbol{x}^{\text{txt}}$. Specifically, a vision encoder E^{img} and a text encoder E^{txt} are employed to transfer each individual image and the text sequence into a sequence of visual tokens $\mathbf{V} = (v_1, v_2, \ldots, v_{N_{\text{img}}})$ and textual tokens $\mathbf{T} = (t_1, t_2, \ldots, t_{N_{\text{txt}}})$, respectively. Subsequently, the visual tokens and the text tokens are fused to form a multi-modal embedding representation $\mathbf{U} = (\mathbf{V}, \mathbf{T})$, which is then fed into a pre-trained large language model $M_{\text{LLM}}(\cdot | \theta_{\text{LLM}})$ for conditional generation. The joint probability of the output response sequence $\mathbf{R} = (r_1, r_2, \ldots, r_{N_{\text{res}}})$ is modeled as:

$$P_{\theta_{\text{LLM}}}(\mathbf{R} \mid \mathbf{U}) = \prod_{i=1}^{N_{\text{res}}} P_{\theta_{\text{LLM}}}(r_i \mid \mathbf{U}, r_{< i}), \qquad (1)$$

where r_i is conditioned on the multi-modal input embedding U and the previously generated tokens $r_{\langle i}$.

The optimization objective during training is to minimize the cross-entropy loss between the generated sequence and the reference sequence $R^* = (r_1^*, r_2^*, \ldots, r_{N_{\text{res}}}^*)$, which is formulated as:

$$\mathcal{L}(\mathbf{R}^* \mid \mathbf{U}) = -\sum_{i=1}^{N_{\text{res}}} \log P_{\theta_{\text{LLM}}}(r_i^* \mid \mathbf{U}, r_{< i}).$$
(2)

5.2 Layer-wise Dense Connector

To better preserve fine-grained structural information critical for ophthalmic understanding, we propose a **Layerwise Dense Connector** (LDC) that densely aggregates multi-scale visual features from different layers of the vision encoder.

Specifically, given an input image, the vision encoder $M_{\text{SLP}}(\cdot | \theta_{\text{SLP}})$ produces a sequence of hierarchical features



Figure 4: Model architecture of EyecareGPT.

 $H = (h_1, h_2, \ldots, h_l)$. Conventional approaches typically utilize only the final feature h_l , risking the loss of low-level details important for fine-grained tasks. To address this, LDC performs average pooling over shallow-to-intermediate layers and deep layers separately to obtain fine-grained features F^{fine} and coarse-grained features F^{coarse} :

$$F^{\text{fine}} = \text{AvgPooling}(h_1, \dots, h_k), \qquad (3)$$

$$F^{\text{coarse}} = \text{AvgPooling}(h_{k+1}, \dots, h_{l-1}).$$
(4)

The final visual representation is obtained by concatenating F^{fine} , F^{coarse} , and h_l along the channel dimension:

$$F = \text{Concatenate}(F^{\text{fine}}, F^{\text{coarse}}, h_l, \text{dim=channel}).$$
 (5)

The fused feature F is then projected to align with the textual embeddings and fed into the LLMs for multi-modal reasoning.

This simple yet effective design enhances fine-grained visual understanding while maintaining overall model efficiency.

5.3 Adaptive Anyres Mechanism

Ophthalmic clinical imaging exhibits diverse resolution distributions, complicating consistent visual modeling. To address this, we propose an **Adaptive Anyres mechanism** that dynamically adapts to varying input scales while preserving fine-grained details.

Specifically, for an input image x^{img} , if its resolution exceeds a predefined threshold, it is partitioned into four subimages centered around the center-symmetric point and concatenated with the original image:

$$x_{\text{anyres}}^{\text{img}} = (x^{\text{img}}, x_1^{\text{sub}}, x_2^{\text{sub}}, x_3^{\text{sub}}, x_4^{\text{sub}}) \,. \tag{6}$$

The original and sub-images are encoded by the LDC-based vision encoder $E^{\rm LDC},$ and adaptive pooling is applied to the

features of the sub-images before concatenating them with the original feature:

$$V = \text{Concat}\left(E^{\text{LDC}}(x^{\text{img}}), \text{ AdaptPooling}\left(E^{\text{LDC}}(x^{\text{sub}}_{i})_{i=1}^{4}\right)\right)$$
(7)

This simple mechanism enables the model to dynamically adapt to heterogeneous resolution distributions, improving the consistency and robustness of multi-modal ophthalmic modeling with minimal computational overhead.

6 Experiments

6.1 Data and Experimental Setup

Data Details. We follow a three-stage training paradigm, first performing alignment and supervised fine-tuning on a mixture of LLaVA (Liu et al. 2023) and PubMedVision (Chen et al. 2024b) data, and then continuing training on Eyecare-100K to enhance domain-specific performance. Detailed settings are provided in the Appendix. We systematically evaluate our model on both the proposed Eyecare-Bench and the existing ophthalmology-related benchmark, OmniMedVQA(Hu et al. 2024), ensuring a comprehensive assessment of its generalization ability and diagnostic performance.

Model Details. We conduct a zero-shot evaluation on 12 representative LVLMs, including eight open-world LVLMs (e.g., LLaVA-v1.5 (Liu et al. 2023), Qwen2.5-VL (Bai et al. 2025), InternVL2.5 (Chen et al. 2025), mPLUG-Owl3 (Ye et al. 2024), Yi-VL (AI et al. 2025), MiniCPM-V2.6 (Yao et al. 2024), gemma-3 (Team et al. 2025), Claude3.5) and four Med-LVLMs (e.g., Med-Flamingo (Moor et al. 2023), LLaVA-Med (Li et al. 2023), MedVLM-R1 (Pan et al. 2025), HealthGPT (Lin et al. 2025)). Models (LLaVA-v1.5, mPLUG-Owl3, gemma-3, Med-Flamingo) are not included in the evaluation when they fail to respond to the report generation instructions correctly.

6.2 Main Experiments

Visual Question Answering. The experimental results for closed and open VQA tasks are presented in Table 2 and Table 3, respectively. The main observations are as follows: **(i) Superior Performance:** EyecareGPT-3.8B achieves SoTA performance across all evaluation metrics with a relatively small parameter size. For closed VQA, it attains an average accuracy of 84.56%, significantly outperforming the second-best model at 62.95% (Qwen2.5-VL-7B). In open VQA, EyecareGPT-3.8B achieves the best F1-Bio score of 90.55, demonstrating outstanding fine-grained ophthalmic image understanding. EyecareGPT-7B, further scaled from a different base model, shows continued performance gains, validating the scalability of the approach. **(ii) Limitations of Existing Medical LVLMs:** Due to the lack of high-quality, domain-specific data, existing Med-LVLMs

Table 2: Performance comparison between EyecareGPT-3.8B and other baseline methods on the *closed QA* task from the opensourced OmniMedVQA and our proposed Eyecare-Bench. Where FS. denotes Fluorescein Staining Images modality. We use **bold** text to indicate the best results and underline to indicate the second-best results.

Model			Ey	ecare -Ben	ch			Omni	OmniMedVQA		
	FS.	Slit-Lamp	OCT	Fundus	FA-ICGA	UBM	СТ	ОСТ	Fundus	8	
			Generalis	st Models							
LLaVA-1.5-7B (Liu et al. 2023)	20.43	65.22	30.52	12.58	6.84	20.26	19.01	51.70	26.40	28.22	
Qwen2.5-VL-7B (Bai et al. 2025)	31.74	75.71	57.86	44.90	75.79	68.66	74.65	68.74	68.46	62.95	
InternVL-2.5-8B (Chen et al. 2025)	32.61	58.57	52.31	37.88	73.62	62.26	61.97	78.67	77.36	59.47	
mPLUG-Owl3-7B (Ye et al. 2024)	16.09	41.43	55.75	30.07	60.10	52.45	71.83	63.56	36.66	47.55	
Yi-VL-6B (AI et al. 2025)	36.52	50.00	50.07	20.40	55.26	58.64	59.86	63.84	36.12	47.86	
MiniCPM-V2.6-8B (Yao et al. 2024)	25.22	58.33	59.05	16.95	79.11	73.47	66.90	86.81	77.31	60.35	
Gemma-3-4B (Team et al. 2025)	22.17	71.67	46.10	25.71	60.10	63.33	33.10	53.48	57.95	48.18	
Claude-3.5	35.22	70.97	64.07	32.28	68.28	63.11	52.82	78.96	63.07	58.75	
			Medical	Models							
Med-Flamingo-8.3B (Moor et al. 2023)	34.78	34.48	33.16	19.39	40.57	40.94	31.69	26.96	29.11	32.34	
LLaVA-Med-7B (Li et al. 2023)	12.61	26.67	37.25	12.98	39.73	31.98	20.42	26.81	29.38	26.54	
MedVLM-R1-2B (Pan et al. 2025)	31.14	64.41	59.13	42.52	55.02	56.72	63.83	71.17	76.76	57.86	
HealthGPT-M3-3.8B (Lin et al. 2025)	41.30	63.33	63.28	20.66	77.80	61.19	69.72	75.11	63.86	59.58	
EyecareGPT-3.8B	<u>60.87</u>	77.03	<u>89.76</u>	<u>75.10</u>	<u>91.43</u>	<u>81.66</u>	85.21	100.00	100.00	<u>84.56</u>	
EyecareGPT-7B	61.43	77.64	90.09	82.25	92.96	86.78	84.33	<u>99.26</u>	<u>99.56</u>	86.03	

Table 3: Performance comparison between EyecareGPT-3.8B and other baseline methods on the *open QA* task from our proposed Eyecare-Bench. We use **bold** text to indicate the best results and underline to indicate the second-best results.

Model	OCT		Fundus		FA-	ICGA	U	ВМ	(СТ		
	F1-Bio	Rouge-L	F1-Bio	Rouge-L	F1-Bio	Rouge-L	F1-Bio	Rouge-L	F1-Bio	Rouge-L		
Generalist Models												
LLaVA-1.5-7B (Liu et al. 2023)	52.60	15.35	67.57	18.38	22.57	7.51	62.29	15.90	27.88	10.18		
Qwen2.5-VL-7B (Bai et al. 2025)	66.55	23.21	81.63	28.82	62.29	21.28	81.74	16.74	69.20	18.22		
InternVL-2.5-8B (Chen et al. 2025)	63.51	20.94	71.44	22.92	46.84	15.86	64.24	21.03	48.78	16.68		
mPLUG-Owl3-7B (Ye et al. 2024)	42.19	19.86	79.27	25.99	31.08	9.69	56.90	22.34	60.36	18.03		
Yi-VL-6B (AI et al. 2025)	56.71	20.60	71.15	22.24	17.77	7.90	59.86	20.46	31.52	14.44		
MiniCPM-V2.6-8B (Yao et al. 2024)	63.60	26.88	78.13	26.92	42.71	11.95	69.20	24.10	62.83	18.88		
Gemma-3-4B (Team et al. 2025)	60.29	20.45	74.48	24.93	38.67	9.20	80.96	23.83	64.04	23.16		
Claude-3.5	62.96	21.20	85.93	28.87	42.06	12.57	78.75	26.48	62.22	16.75		
			Med	ical Models								
Med-Flamingo-8.3B (Moor et al. 2023)	29.13	11.47	45.96	14.45	32.31	10.76	34.90	10.65	38.30	10.59		
LLaVA-Med-7B (Li et al. 2023)	51.79	23.25	82.38	26.03	32.35	9.31	68.92	20.87	66.95	17.03		
MedVLM-R1-2B (Pan et al. 2025)	60.29	19.26	76.46	25.10	50.40	18.13	63.53	22.72	62.83	18.42		
HealthGPT-M3-3.8B (Lin et al. 2025)	51.45	13.15	61.55	16.97	56.24	17.07	71.05	21.07	57.80	15.87		
EyecareGPT-3.8B	95.53	49.68	<u>90.79</u>	37.03	86.75	49.21	<u>95.78</u>	47.61	83.90	36.11		
EyecareGPT-7B	96.26	50.10	90.88	38.13	87.86	51.24	96.60	47.26	87.27	36.70		

show no significant advantage over general LVLMs in ophthalmic diagnostic tasks. The consistent performance of the EyecareGPT models further highlights the critical role of Eyecare-100K in enhancing fine-grained domain-specific visual understanding. (iii) Multidimensional Evaluation Metrics: Compared to the existing benchmark OmniMed-VQA, Eyecare-Bench covers a broader range of imaging modalities and task types, posing greater challenges and practical evaluation. The results demonstrate that Eyecare-Bench effectively reveals performance bottlenecks and deficiencies of LVLMs in ophthalmology, offering valuable insights for model optimization.

Report Generation. Table 4 presents the experimental results for the report generation task: (i) The EyecareGPT achieves the best performance across all evaluation metrics. Under both GPT-based and traditional evaluations, the diagnostic reports generated by EyecareGPT exceed 50% accuracy across three imaging modalities and show excellent performance on structural and linguistic consistency metrics such as F1-RadGraph and ROUGE-L. These results demonstrate that EyecareGPT accurately understands complex ophthalmic images and generates high-quality professional reports. (ii) Although some general models (e.g., Qwen2.5, Claude3.5) and medical models (e.g., MedVLM, MedVLM-R1) perform reasonably well in open VQA tasks, they show significant deficiencies in structured medical report generation, failing to maintain scores within a usable range. This highlights the current limitations of LVLMs in handling multi-modal heterogeneous ophthalmic data and the urgent need for targeted optimization.

To further systematically reveal performance differences among models, we introduce a GPT-4-based multi-

Table 4: Performance comparison between EyecareGPT-3.8B and other baseline methods on the *report generation* task from our proposed Eyecare-100K benchmark. We use **bold** text to indicate the best results and <u>underline</u> to indicate the second-best results.

Model		FA-ICGA	1		UBM			СТ			
mouel	Acc ^{GPT}	F1-Rad	Rouge-L	Acc ^{GPT}	F1-Rad	Rouge-L	Acc ^{GPT}	F1-Rad	Rouge-L		
			Generalist I	Models							
Qwen2.5-VL-7B (Bai et al. 2025)	17.00	6.91	15.54	19.54	4.26	7.21	42.86	11.37	18.23		
InternVL-2.5-8B (Chen et al. 2025)	5.92	5.19	8.51	0.00	3.53	8.39	0.00	8.00	12.68		
Yi-VL-6B (AI et al. 2025)	2.26	5.12	9.13	0.00	2.14	8.03	3.85	6.73	14.68		
MiniCPM-V2.6-8B (Yao et al. 2024)	3.34	6.12	8.56	0.00	3.78	7.33	0.00	5.01	11.88		
Claude-3.5	14.53	6.37	12.96	2.98	9.56	14.38	25.00	11.05	16.23		
			Medical M	[odels							
LLaVA-Med-7B (Li et al. 2023)	0.14	3.53	12.64	0.00	4.63	8.32	0.00	1.02	12.15		
MedVLM-R1-2B (Pan et al. 2025)	11.28	5.76	4.38	2.22	4.05	6.41	32.14	11.87	10.64		
HealthGPT-M3-3.8B (Lin et al. 2025)	14.41	7.30	12.86	3.02	5.53	10.30	17.35	14.47	18.21		
EyecareGPT-3.8B	52.62	25.04	47.91	58.05	42.83	57.04	50.00	35.39	42.73		
EyecareGPT-7B	53.91	26.04	48.32	60.06	42.98	58.43	52.43	36.19	43.54		



Figure 5: GPT-4-based evaluationon results for report generation task in (a) FA, (b) UBM, and (c) CT modalities.

dimensional evaluation approach, with results shown in Figure 5. EyecareGPT consistently demonstrates stable and superior performance across all imaging modalities, with more than 50% of its generated reports rated as clinically usable, showing particularly strong potential in FA and UBM modalities. In contrast, other models perform relatively better on CT tasks but show clear declines in FA and UBM, reflecting the imbalance of modality distribution in current public datasets. Overall, these results further validate the robustness of EyecareGPT in multi-modal tasks and demonstrate the scientific value of our proposed multidimensional evaluation framework in guiding future LVLM optimization.

6.3 Ablation and In-Depth Study

Effect of LDC and Anyres. We integrate the Layer-wise Dense Connector (LDC) and the Adaptive Anyres Mechanism to enhance the ability of the model to capture finegrained structural information and multi-scale visual fea-



Figure 6: Results after fine-tuning on Eyecare-100K.

tures, and validate the contribution of each module through ablation studies. As shown in Table 5, LDC, as a plug-andplay module, consistently improves model performance across all three task types by integrating multi-scale visual features. Similarly, the Adaptive Anyres Mechanism strengthens the model's capability for fine-grained region recognition, achieving notable gains over the baseline. Furthermore, the combined use of LDC and Anyres synergistically balances multi-scale visual modeling and resolution adaptability, significantly boosting diagnostic accuracy and report generation quality, thereby demonstrating the effectiveness and generalizability of the architecture in intelligent ophthalmic diagnosis.

Ablation on Eyecare-100K. We validate Eyecare-100K's effectiveness in enhancing visual understanding of complex ophthalmic scenarios by comparing model performance before and after fine-tuning, as shown in Figure 6. In closed VQA tasks, the average accuracy of the model improves from 65.30% to 84.56%, reaching 100% accuracy on the OmniMedVQA subtask. In open VQA and report generation tasks, the average F1-Bio score increases from 64.89 to 90.55, and the GPT-evaluated diagnostic accuracy rises from 11.48% to 53.56%. These significant improvements demonstrate the high quality and broad applicability of Eyecare-100K in constructing multi-modal, multi-task instruction datasets.

Task	AnvRes	LD	C]	Eyecare-10	0K				C)mniMeo	łVQA	Avg.
14011		22	-	FS	s. s	L.	ОСТ	Fundus	FA-ICGA	UB	м	СТ	0	ОСТ	Fundus	
βĄ	-	-		60.0	00 68	3.92	87.19	73.77	88.62	78.2	25	82.39	9	9.26	98.65	81.89
	•	-		60.7	<u>71 76</u>	6.47	88.61	76.42	90.78	80.5	5	83.95	9	9.41	100.00	84.10
Se	-	C		60.0	00 70	0.23	87.58	74.78	89.20	79.9	96	83.45	<u>9</u>	9.50	<u>99.56</u>	82.79
CIC	۲	C	>	60.8	87 77	.03	89.76	75.10	91.43	81.6	66	85.21	10	00.00	100.00	84.56
Task	AnvRes	LD	ОСТ			Fur	ndus	FA	-ICGA		τ	JBM		СТ		
	,			F1-F	Bio Rou	ıge-L	F1-Bio	Rouge-L	F1-Bio	Roug	e-L F	1-Bio	Ro	uge-L	F1-Bio	Rouge-L
A	-	-		95.	12 48	3.32	87.15	36.46	85.89	46.7	7	93.48	4	5.25	79.38	32.23
Ŏ	•	-		<u>95.3</u>	<u>38 49</u>	0.17	89.76	36.90	86.34	48.1	.5	95.66	4	6.94	82.23	33.43
per	-	C	>	95.2	20 48	8.98	88.90	36.77	85.95	47.7	'5	94.73	4	5.79	80.98	32.93
Ō	Ø	C		95.	53 49	.68	90.79	37.03	86.75	49.2	21	95.78	4	7.61	83.90	36.11
Task	AnvRes	LDC	1		FA-l	CGA	UBM								СТ	
Tubit	1111/1000	220	Ac	c ^{GPT}	Score ^{Avg}	F1-Rad	Rouge	L Acc ^{GPT}	ScoreAvg	F1-Rad	Rouge	L Ac	c ^{GPT}	ScoreAvg	F1-Rad	Rouge-L
Gen	-	-	45	5.37	83.02	23.56	46.71	32.48	72.36	39.14	53.78	42	2.71	76.71	29.42	40.47
L L	Ø	-	<u>51</u>	1.37	85.13	24.76	47.80	54.86	<u>79.07</u>	<u>41.36</u>	56.10	45	5.88	78.35	32.27	<u>41.52</u>
īodi	-	Ø	47	7.83	84.85	24.62	46.81	44.42	77.45	40.29	54.43	44	1.97	77.25	30.33	40.65
Re	V	Ø	52	2.62	85.49	25.04	47.91	58.05	79.90	42.83	57.04	50	0.00	78.64	35.39	42.73

Table 5: Ablation study of the effect of the individual module for three tasks.



Figure 7: Physician preference for generated reports (a) and EyeEval reliability (b, c).

Expert Physician Review. We further conduct an expert physician review of the report generation task in Eyecare-Bench. We recruited ten clinicians to rank the responses from EyecareGPT-3.8B, Qwen-VL-7B, Claude 3.5, LLaVA-Med, MedVLM-R1, and HealthGPT-M3, and to select the report with the highest satisfaction. We randomly sample 500 report generation VOA pairs along with the answers generated by the aforementioned six models and randomly order them for the doctors to choose from. The final results of the doctors' selections are shown in Figure 7 (a), indicating that the reports generated by EyecareGPT are more satisfactory to clinicians and better meet clinical needs. Simultaneously, we also asked the doctors to score the reports generated by EvecareGPT according to our proposed EveEval evaluation system. The report scores and deduction details across four dimensions evaluated by GPT-4 and doctors are shown in Figure 7 (b) and (c) respectively, demonstrating that EyeEval and doctor diagnoses have high consistency and reliability.

7 Conclusion

We propose **Eyecare Kit**, a comprehensive framework that addresses critical challenges in intelligent ophthalmic di-

agnosis through tailored dataset, benchmark, and model. Experiments validate the effectiveness of **Eyecare-100K**, **Eyecare-Bench**, and **EyecareGPT**, with EyecareGPT achieving SOTA results. We believe Eyecare Kit lays a solid foundation for future advances in domain-specific Med-LVLMs and ophthalmic AI applications.

References

AI, .; :; Young, A.; Chen, B.; Li, C.; Huang, C.; Zhang, G.; Zhang, G.; Wang, G.; Li, H.; Zhu, J.; Chen, J.; Chang, J.; Yu, K.; Liu, P.; Liu, Q.; Yue, S.; Yang, S.; Yang, S.; Xie, W.; Huang, W.; Hu, X.; Ren, X.; Niu, X.; Nie, P.; Li, Y.; Xu, Y.; Liu, Y.; Wang, Y.; Cai, Y.; Gu, Z.; Liu, Z.; and Dai, Z. 2025. Yi: Open Foundation Models by 01.AI. arXiv:2403.04652.

Alayrac, J.-B.; Donahue, J.; Luc, P.; Miech, A.; Barr, I.; Hasson, Y.; Lenc, K.; Mensch, A.; Millican, K.; Reynolds, M.; Ring, R.; Rutherford, E.; Cabi, S.; Han, T.; Gong, Z.; Samangooei, S.; Monteiro, M.; Menick, J.; Borgeaud, S.; Brock, A.; Nematzadeh, A.; Sharifzadeh, S.; Binkowski, M.; Barreira, R.; Vinyals, O.; Zisserman, A.; and Simonyan, K. 2022. Flamingo: a Visual Language Model for Few-Shot Learning. arXiv:2204.14198.

AlSaad, R.; Abd-Alrazaq, A.; Boughorbel, S.; Ahmed, A.; Renault, M.-A.; Damseh, R.; and Sheikh, J. 2024. Multimodal large language models in health care: applications, challenges, and future outlook. *Journal of medical Internet research*, 26: e59505.

Bai, S.; Chen, K.; Liu, X.; Wang, J.; Ge, W.; Song, S.; Dang, K.; Wang, P.; Wang, S.; Tang, J.; Zhong, H.; Zhu, Y.; Yang, M.; Li, Z.; Wan, J.; Wang, P.; Ding, W.; Fu, Z.; Xu, Y.; Ye, J.; Zhang, X.; Xie, T.; Cheng, Z.; Zhang, H.; Yang, Z.; Xu, H.; and Lin, J. 2025. Qwen2.5-VL Technical Report. *arXiv* preprint arXiv:2502.13923.

Bakarov, A. 2018. A survey of word embeddings evaluation methods. *arXiv preprint arXiv:1801.09536*.

Balas, M.; Ramalingam, V.; Pandya, B.; Abdelaal, A.; and Shi, R. B. 2024. Adaptive optics imaging in ophthalmology: redefining vision research and clinical practice. *JFO Open Ophthalmology*, 100116.

Cen, L.-P.; Ji, J.; Lin, J.-W.; Ju, S.-T.; Lin, H.-J.; Li, T.-P.; Wang, Y.; Yang, J.-F.; Liu, Y.-F.; Tan, S.; et al. 2021. Automatic detection of 39 fundus diseases and conditions in retinal photographs using deep neural networks. *Nature communications*, 12(1): 4828.

Chen, J.; Gui, C.; Ouyang, R.; Gao, A.; Chen, S.; Chen, G.; Wang, X.; Cai, Z.; Ji, K.; Wan, X.; et al. 2024a. Towards injecting medical visual knowledge into multimodal llms at scale. In *Proceedings of the 2024 Conference on Empirical Methods in Natural Language Processing*, 7346–7370.

Chen, J.; Gui, C.; Ouyang, R.; Gao, A.; Chen, S.; Chen, G. H.; Wang, X.; Zhang, R.; Cai, Z.; Ji, K.; et al. 2024b. Huatuogptvision, towards injecting medical visual knowledge into multimodal llms at scale. *arXiv preprint arXiv:2406.19280*.

Chen, Z.; Wang, W.; Cao, Y.; Liu, Y.; Gao, Z.; Cui, E.; Zhu, J.; Ye, S.; Tian, H.; Liu, Z.; Gu, L.; Wang, X.; Li, Q.; Ren, Y.; Chen, Z.; Luo, J.; Wang, J.; Jiang, T.; Wang, B.; He, C.; Shi, B.; Zhang, X.; Lv, H.; Wang, Y.; Shao, W.; Chu, P.; Tu, Z.; He, T.; Wu, Z.; Deng, H.; Ge, J.; Chen, K.; Zhang, K.; Wang, L.; Dou, M.; Lu, L.; Zhu, X.; Lu, T.; Lin, D.; Qiao, Y.; Dai, J.; and Wang, W. 2025. Expanding Performance Boundaries of Open-Source Multimodal Models with Model, Data, and Test-Time Scaling. arXiv:2412.05271.

Culy, C.; and Riehemann, S. Z. 2003. The limits of N-gram translation evaluation metrics. In *Proceedings of Machine Translation Summit IX: Papers*.

Gholami, P.; Roy, P.; Parthasarathy, M. K.; and Lakshminarayanan, V. 2020. OCTID: Optical coherence tomography image database. *Computers & Electrical Engineering*, 81: 106532.

Google AI. ???? Gemini 2.0 Flash. Technical Blog Post. Accessed April 11, 2025.

Guo, Y.; and Huang, W. 2025. LLaVA-NeXT-Med: Medical Multimodal Large Language Model. In 2025 Asia-Europe Conference on Cybersecurity, Internet of Things and Soft Computing (CITSC), 474–477. IEEE.

Hu, Y.; Li, T.; Lu, Q.; Shao, W.; He, J.; Qiao, Y.; and Luo, P. 2024. Omnimedvqa: A new large-scale comprehensive evaluation benchmark for medical lvlm. In *Proceedings of the IEEE/CVF Conference on Computer Vision and Pattern Recognition*, 22170–22183.

Hurst, A.; Lerer, A.; Goucher, A. P.; Perelman, A.; Ramesh, A.; Clark, A.; Ostrow, A.; Welihinda, A.; Hayes, A.; Radford, A.; et al. 2024. Gpt-4o system card. *arXiv preprint arXiv:2410.21276*.

Kermany, D.; Zhang, K.; and Goldbaum, M. 2018a. Labeled optical coherence tomography (OCT) and chest Xray images for classification (2018). *Mendeley Data*, *v2 https://doi. org/10.17632/rscbjbr9sj https://nihcc. app. box. com/v/ChestXray-NIHCC*.

Kermany, D.; Zhang, K.; and Goldbaum, M. 2018b. Large Dataset of Labeled Optical Coherence Tomography (OCT) and Chest X-Ray Images.

Knipe, H. 2005. Radiopaedia's mission is to create the best radiology reference the world has ever seen and to make it available for free, forever, for all. Website. https://radiopaedia.org.

Kulyabin, M.; Zhdanov, A.; Nikiforova, A.; Stepichev, A.; Kuznetsova, A.; Ronkin, M.; Borisov, V.; Bogachev, A.; Korotkich, S.; Constable, P. A.; et al. 2024. Octdl: Optical coherence tomography dataset for image-based deep learning methods. *Scientific data*, 11(1): 365.

Lee, J.; Yoon, W.; Kim, S.; Kim, D.; Kim, S.; So, C. H.; and Kang, J. 2020. BioBERT: a pre-trained biomedical language representation model for biomedical text mining. *Bioinformatics*, 36(4): 1234–1240.

Li, C.; Wong, C.; Zhang, S.; Usuyama, N.; Liu, H.; Yang, J.; Naumann, T.; Poon, H.; and Gao, J. 2023. Llava-med: Training a large language-and-vision assistant for biomedicine in one day. *Advances in Neural Information Processing Systems*, 36: 28541–28564.

Li, J.; Guan, Z.; Wang, J.; Cheung, C. Y.; Zheng, Y.; Lim, L.-L.; Lim, C. C.; Ruamviboonsuk, P.; Raman, R.; Corsino, L.; et al. 2024. Integrated image-based deep learning and language models for primary diabetes care. *Nature medicine*, 30(10): 2886–2896.

Liang, Q.; and Zhang, Y. 2022. *Interpretation of Slit Lamp Images of Anterior Segment Diseases*. People's Medical Publishing House Co., Ltd. ISBN 9787117330626.

Lin, C.-Y.; and Hovy, E. 2003. Automatic evaluation of summaries using n-gram co-occurrence statistics. In *Proceedings of the 2003 human language technology conference of the North American chapter of the association for computational linguistics*, 150–157.

Lin, J.; Yin, H.; Ping, W.; Molchanov, P.; Shoeybi, M.; and Han, S. 2024. Vila: On pre-training for visual language models. In *Proceedings of the IEEE/CVF conference on computer vision and pattern recognition*, 26689–26699.

Lin, T.; Zhang, W.; Li, S.; Yuan, Y.; Yu, B.; Li, H.; He, W.; Jiang, H.; Li, M.; Song, X.; et al. 2025. HealthGPT: A Medical Large Vision-Language Model for Unifying Comprehension and Generation via Heterogeneous Knowledge Adaptation. *arXiv preprint arXiv:2502.09838*.

Liu, H.; Li, C.; Li, Y.; and Lee, Y. J. 2024a. Improved baselines with visual instruction tuning. In *Proceedings of the IEEE/CVF Conference on Computer Vision and Pattern Recognition*, 26296–26306. Liu, H.; Li, C.; Li, Y.; Li, B.; Zhang, Y.; Shen, S.; and Lee, Y. J. 2024b. Llavanext: Improved reasoning, ocr, and world knowledge.

Liu, H.; Li, C.; Wu, Q.; and Lee, Y. J. 2023. Visual Instruction Tuning. In *NeurIPS*.

Liu, R.; Wang, X.; Wu, Q.; Dai, L.; Fang, X.; Yan, T.; Son, J.; Tang, S.; Li, J.; Gao, Z.; et al. 2022. DeepDRiD: diabetic retinopathy-grading and image quality estimation challenge. Patterns 3 (6), 100512 (2022).

Luo, Y.; Shi, M.; Tian, Y.; Elze, T.; and Wang, M. 2023. Harvard glaucoma detection and progression: A multimodal multitask dataset and generalization-reinforced semi-supervised learning. In *Proceedings of the IEEE/CVF International Conference on Computer Vision*, 20471–20482.

Moor, M.; Huang, Q.; Wu, S.; Yasunaga, M.; Dalmia, Y.; Leskovec, J.; Zakka, C.; Reis, E. P.; and Rajpurkar, P. 2023. Med-flamingo: a multimodal medical few-shot learner. In *Machine Learning for Health (ML4H)*, 353–367. PMLR.

on Ocular Disease Intelligent Recognition, P. U. I. C. 2019. Ocular Disease Intelligent Recognition (ODIR-2019). Accessed [Date you accessed the dataset].

Ovreiu, S.; Paraschiv, E.-A.; and Ovreiu, E. 2021. Deep learning & digital fundus images: Glaucoma detection using DenseNet. In 2021 13th international conference on electronics, computers and artificial intelligence (ECAI), 1–4. IEEE.

Pan, J.; Liu, C.; Wu, J.; Liu, F.; Zhu, J.; Li, H. B.; Chen, C.; Ouyang, C.; and Rueckert, D. 2025. Medvlm-r1: Incentivizing medical reasoning capability of vision-language models (vlms) via reinforcement learning. *arXiv preprint arXiv:2502.19634*.

Papineni, K.; Roukos, S.; Ward, T.; and Zhu, W.-J. 2002. Bleu: a method for automatic evaluation of machine translation. In *Proceedings of the 40th annual meeting of the Association for Computational Linguistics*, 311–318.

Porwal, P.; Pachade, S.; Kamble, R.; Kokare, M.; Deshmukh, G.; Sahasrabuddhe, V.; and Meriaudeau, F. 2018. Indian diabetic retinopathy image dataset (IDRiD): a database for diabetic retinopathy screening research. *Data*, 3(3): 25.

Qin, Z.; Yin, Y.; Campbell, D.; Wu, X.; Zou, K.; Tham, Y.-C.; Liu, N.; Zhang, X.; and Chen, Q. 2025. LMOD: A Large Multimodal Ophthalmology Dataset and Benchmark for Large Vision-Language Models. arXiv:2410.01620.

Ren, T.; Liu, S.; Zeng, A.; Lin, J.; Li, K.; Cao, H.; Chen, J.; Huang, X.; Chen, Y.; Yan, F.; et al. 2024. Grounded sam: Assembling open-world models for diverse visual tasks. *arXiv* preprint arXiv:2401.14159.

Rodríguez, M. A.; AlMarzouqi, H.; and Liatsis, P. 2022. Multi-label retinal disease classification using transformers. *IEEE Journal of Biomedical and Health Informatics*.

Saab, K.; Tu, T.; Weng, W.-H.; Tanno, R.; Stutz, D.; Wulczyn, E.; Zhang, F.; Strother, T.; Park, C.; Vedadi, E.; et al. 2024.

Capabilities of gemini models in medicine. *arXiv preprint arXiv:2404.18416*.

Shi, D.; Zhang, W.; Yang, J.; Huang, S.; Chen, X.; Yusufu, M.; Jin, K.; Lin, S.; Liu, S.; Zhang, Q.; et al. 2024. EyeCLIP: A visual-language foundation model for multi-modal oph-thalmic image analysis. *arXiv preprint arXiv:2409.06644*.

Siragusa, I.; Contino, S.; La Ciura, M.; Alicata, R.; and Pirrone, R. 2024. Medpix 2.0: a comprehensive multimodal biomedical dataset for advanced AI applications. *arXiv* preprint arXiv:2407.02994.

Team, G.; Georgiev, P.; Lei, V. I.; Burnell, R.; Bai, L.; Gulati, A.; Tanzer, G.; Vincent, D.; Pan, Z.; Wang, S.; et al. 2024. Gemini 1.5: Unlocking multimodal understanding across millions of tokens of context. *arXiv preprint arXiv:2403.05530.*

Team, G.; Kamath, A.; Ferret, J.; Pathak, S.; Vieillard, N.; Merhej, R.; Perrin, S.; Matejovicova, T.; Ramé, A.; Rivière, M.; et al. 2025. Gemma 3 Technical Report. *arXiv preprint arXiv:2503.19786.*

Wang, P.; Bai, S.; Tan, S.; Wang, S.; Fan, Z.; Bai, J.; Chen, K.; Liu, X.; Wang, J.; Ge, W.; Fan, Y.; Dang, K.; Du, M.; Ren, X.; Men, R.; Liu, D.; Zhou, C.; Zhou, J.; and Lin, J. 2024. Qwen2-VL: Enhancing Vision-Language Model's Perception of the World at Any Resolution. *arXiv preprint arXiv:2409.12191*.

Xu, D.; Chen, Y.; Wang, J.; Huang, Y.; Wang, H.; Jin, Z.; Wang, H.; Yue, W.; He, J.; Li, H.; et al. 2024a. Mlevlm: Improve multi-level progressive capabilities based on multi-modal large language model for medical visual question answering. In *Findings of the Association for Computational Linguistics ACL 2024*, 4977–4997.

Xu, P.; Chen, X.; Zhao, Z.; and Shi, D. 2024b. Unveiling the clinical incapabilities: a benchmarking study of GPT-4V (ision) for ophthalmic multimodal image analysis. *British Journal of Ophthalmology*, 108(10): 1384–1389.

Yao, Y.; Yu, T.; Zhang, A.; Wang, C.; Cui, J.; Zhu, H.; Cai, T.; Li, H.; Zhao, W.; He, Z.; et al. 2024. MiniCPM-V: A GPT-4V Level MLLM on Your Phone. *arXiv preprint arXiv:2408.01800*.

Ye, J.; Xu, H.; Liu, H.; Hu, A.; Yan, M.; Qian, Q.; Zhang, J.; Huang, F.; and Zhou, J. 2024. mPLUG-Owl3: Towards Long Image-Sequence Understanding in Multi-Modal Large Language Models. arXiv:2408.04840.

You, K.; Zhang, H.; Schoop, E.; Weers, F.; Swearngin, A.; Nichols, J.; Yang, Y.; and Gan, Z. 2024. Ferret-ui: Grounded mobile ui understanding with multimodal llms. In *European Conference on Computer Vision*, 240–255. Springer.

Yu, F.; Endo, M.; Krishnan, R.; Pan, I.; Tsai, A.; Reis, E. P.; Fonseca, E. K. U. N.; Lee, H. M. H.; Abad, Z. S. H.; Ng, A. Y.; et al. 2023. Evaluating progress in automatic chest x-ray radiology report generation. *Patterns*, 4(9). Yu, H.; Yi, S.; Niu, K.; Zhuo, M.; and Li, B. 2025. UMIT: Unifying Medical Imaging Tasks via Vision-Language Models. *arXiv preprint arXiv:2503.15892.*

Zambrano Chaves, J. M.; Huang, S.-C.; Xu, Y.; Xu, H.; Usuyama, N.; Zhang, S.; Wang, F.; Xie, Y.; Khademi, M.; Yang, Z.; et al. 2025. A clinically accessible small multimodal radiology model and evaluation metric for chest X-ray findings. *Nature Communications*, 16(1): 3108.

Zhai, X.; Mustafa, B.; Kolesnikov, A.; and Beyer, L. 2023. Sigmoid loss for language image pre-training. In *Proceedings of the IEEE/CVF international conference on computer vision*, 11975–11986.

Zhang, T.; Kishore, V.; Wu, F.; Weinberger, K. Q.; and Artzi, Y. 2019a. Bertscore: Evaluating text generation with bert. *arXiv preprint arXiv:1904.09675*.

Zhang, W.; Lin, T.; Liu, J.; Shu, F.; Li, H.; Zhang, L.; Wanggui, H.; Zhou, H.; Lv, Z.; Jiang, H.; et al. 2024. Hyperllava: Dynamic visual and language expert tuning for multimodal large language models. *arXiv preprint arXiv:2403.13447*.

Zhang, W.; Tang, S.; Cao, Y.; Pu, S.; Wu, F.; and Zhuang, Y. 2019b. Frame augmented alternating attention network for video question answering. *IEEE Transactions on Multimedia*, 22(4): 1032–1041.

Zhang, W.; Zhu, L.; Hallinan, J.; Zhang, S.; Makmur, A.; Cai, Q.; and Ooi, B. C. 2022. Boostmis: Boosting medical image semi-supervised learning with adaptive pseudo labeling and informative active annotation. In *Proceedings of the IEEE/CVF Conference on Computer Vision and Pattern Recognition*, 20666–20676.

Zhao, H.; Ling, Q.; Pan, Y.; Zhong, T.; Hu, J.-Y.; Yao, J.; Xiao, F.; Xiao, Z.; Zhang, Y.; Xu, S.-H.; Wu, S.-N.; Kang, M.; Wu, Z.; Liu, Z.; Jiang, X.; Liu, T.; and Shao, Y. 2023. Ophtha-LLaMA2: A Large Language Model for Ophthalmology. arXiv:2312.04906.

Appendix

This is the Appendix for "EyecareGPT: Boosting Comprehensive Ophthalmology Understanding with Tailored Dataset, Benchmark and Model". This Appendix is organized as follows:

- Section A presents the details of the experimental implementation, the training process of EyecareGPT, the construction details of Eyecare-100K, and the specific information of Eyecare-Bench.
- Section B presents our more detailed ablation experimental results and a brief experimental analysis.
- Section C shows typical data examples in Eyecare-100K.

A Implementation Details

A.1 Model Details

We use SigLIP-SO400M-Patch14-384 as the visual feature extractor to capture multi-scale visual features, ensuring the model's ability to recognize and understand local details. The visual features are aligned with text embeddings through an MLP and jointly fed into the large language model for conditional output.

EyecareGPT offers two versions: **EyecareGPT-3.8B** and **EyecareGPT-7B**, which are based on Phi-3.5-mini-Instruct and Qwen2.5-7B-Instruct as the pre-trained LLMs, respectively. Table 6 shows the details.

Table 6: Overview of the components of EyecareGPT.

Model	ViT	Adapter	LLM	Params	LoRA Rank
EyecareGPT-3.8B	SigLIP-SO400M	2-layer MLP	Phi-3.5-mini-Instruct	3.8B	64
EyecareGPT-7B	SigLIP-SO400M	2-layer MLP	Qwen2.5-7B-Instruct	7B	64

A.2 Training Details

This study adopts a three-stage training strategy to progressively build the model's visual understanding and intelligent ophthalmic diagnostic capabilities. In the first stage, we train the model on data for aligning from LLaVA-558K and PubMedVision to enhance image description and basic vision-language alignment capabilities. In the second stage, we use supervised finetuning data from LLaVA-665K and PubMedVision to further strengthen the model's visual instruction following and general medical understanding. In the third stage, we perform specialized fine-tuning on Eyecare-100K, focusing on three core task types in clinical ophthalmic diagnosis while avoiding noise introduced by other data sources, thereby optimizing the model's domain-specific adaptability and accuracy. Hyperparameter configurations for each training stage are detailed in Table 7 to ensure training efficiency and convergence.

Table 7: Overview of hyperparameter configurations.

II	Eye	ecareGPT-3.	8B	EyecareGPT-7B					
пуреграгатет	Stage-1	Stage-2	Stage-3	Stage-1	Stage-2	Stage-3			
Optimizer	AdamW	AdamW	AdamW	AdamW	AdamW	AdamW			
Adapter LR	1e-3	2e-5	2e-5	1e-3	2e-5	2e-5			
Learning Rate	/	2e-4	2e-4	/	2e-4	2e-4			
Global Batch Size	256	128	32	256	128	32			
Weight Decay	0	0	0	0	0	0			
Dropout Rate	0	0.05	0.05	0	0.05	0.05			
LR Scheduler	Warm Up	Warm Up	Constant	Warm Up	Warm Up	Constant			
Max Sequence Length	2048	2048	2048	2048	2048	2048			

A.3 Construction details of Eyecare-100K

Data Source Details: In the data collection phase, we gathered eye report data with four modalities – Fluorescein Angiography (FA), Indocyanine Green Angiography (ICGA), Ultrasound Biomicroscopy (UBM), and Computed Tomography (CT) – from our hospital. Specifically, this included 2081 CT images, 3425 UBM images, 15048 FA images, and 2394 ICGA images.

Furthermore, to fully leverage real-world data from existing public datasets, we collected 10 previously published ophthalmological datasets with Fundus and Optical Coherence Tomography (OCT) modalities, as detailed in Table 8. These are all single-modality datasets containing disease classification or grading labels and corresponding images. We also collected three publicly available datasets from Kaggle, including one fluorescence-stained image dataset, one slit-lamp dataset, and one OCT dataset.

Dataset Name	Modality Type	Source
IDRID	Fundus	Aravind Eye Hospital, Madurai, India
ACRIMA	Fundus	University of Jaén, Spain
JSIEC	Fundus	Joint Shantou University - Eye Center, China
ODIR	Fundus	Multiple hospitals in China
MuReD	Fundus	Multiple hospitals in China
DeepDRID	Fundus	Multiple hospitals in India
OCT2017	OCT	Shiley Eye Institute, University of California San Diego, USA
OCTID	OCT	Zhongshan Ophthalmic Center, Sun Yat-sen University, China
OCTDL	OCT	University of Tuebingen, Germany
Kermany	OCT	Multiple sources (publicly available)

Table 8: Overview of Existing Eye Publicly Available Datasets Collected

Rewrite Prompt: For the report data collected from the hospital, we processed it and used Claude 3.7 to construct three types of data. Through discussions with doctors, we identified three key components in the reports: Image Type, Imaging Findings, and Diagnostic Suggestions. Therefore, we designed prompts to guide Claude in generating reports with a unified structure. Figure 8 illustrates the prompt template we designed for the CT modality.

Rewrite Prompt
System Prompt: Transform the following original CT ocular imaging report into a standardized, professional report format. The rewritten version should employ precise and formal medical language, follow a consistent structure, and be suitable for clinical documentation and communication. Ensure clarity, accuracy, and the use of standardized imaging terminology throughout.
Instruction: Original Report Content: (RAW_REPORT)
Rewrite Requirements:
1. Examination Type: Ocular CT Examination
2. Image Findings: Provide detailed, structured descriptions under the following
headings:
 Eyeball: Describe the size, shape, symmetry, and contour of both eyes. Intraocular Structures: Report on the lens, vitreous humor, and any internal ocular structures. Describe any abnormalities such as calcifications, dislocations, or opacities.
- Orbital Tissues: Assess the extraocular muscles, optic nerve, retrobulbar fat, and bony orbit. Note any thickening, displacement, or bone destruction.
- Additional Findings: Include any pathology involving adjacent paranasal sinuses, soft tissue masses, foreign bodies, or other incidental findings. State size, location, density, and nature (if identifiable).
3. Diagnostic Impression & Recommendation:
 Preliminary Impression: Provide a likely diagnosis or differential based on the imaging appearance. If inconclusive, state: "Further correlation with clinical findings is recommended."
- Recommendation: Indicate if follow-up imaging, clinical evaluation, or referral to a specialist is advised.

Figure 8: Rewrite Prompt

Eva	aluation Prompt
Syst the g are t indic corre and Inst quar A)T the r B)T inco C)T inco D)T.	tem Prompt: You are a professional eye doctor. The diagnostic accuracy of generated image report was judged according to the reference report. There ten evaluation indicators, and the calculation method and examples of each cator are given below. For each quantity metric, the number of error types esponding to each metric in the candidate report needs to be counted. For true false indica-tors, 0 is negative and 1 is positive. rruction: There are ten evaluation aspects, of which the first four are ntitative indicators and the last six are true and false indicators: he number of abnormal features in candidate report that are not mentioned in reference report. he number of times the candidate report describes the disease severity brrectly. he number of times the candidate report describes the disease location brrectly. The number of missing key findings compared to the reference report.
 Desi Elev (G, 0	ired output format: yen indicators: [(A, n_A), (B, n_B), (C, n_C), (D, n_D), (E, 0 or 1), (F, 0 or 1), 0 or 1), (H, 0 or 1), (I, 0 or 1), (J, 0 or 1)]

Figure 9: Evaluation Prompt

QA Templates: For the aforementioned datasets that only contain classification or grading labels, we analyzed the data characteristics of their labels and designed different Question-Answering (QA) templates for each. This allowed us to transform the original data into open-ended Question-Answering pairs. Examples of the QA templates are shown in the Table 9.

A.4 GPT-4 Evaluation Prompt

We designed an evaluation system called EyeEval and introduced GPT-4 for the evaluation process. The template used for GPT-4's evaluation is shown in Figure 9. According to the scoring criteria, we grade the reports as follows:

- Excellent Report (90-100): The report contains virtually no errors, the information is relatively complete, the structure is clear, and it does not contain serious clinical errors.
- Usable Report (80-90): The report may contain some minor errors, but overall the information is complete, the structure is clear, and it does not contain serious clinical errors.
- **Report Under Review (60-80):** The report contains numerous errors or missing information, the diagnosis may be inaccurate, or the report structure is disorganized, requiring further review.
- Unusable Report (Below 60): The report contains a large number of errors, severely missing information, diagnostic errors, or contains serious clinical errors, making it unsuitable for clinical decision-making.

B Supplemental Experimental Results

In this section, we include additional experiments to demonstrate the superiority of Eyecare Kit.

B.1 Additional Evaluation Metrics

This section provides a detailed overview of the supplementary evaluation metrics employed in the main experiment. Table 10 outlines the specific metrics used to assess the performance of the open question answering task. Similarly, Table 11 presents the additional evaluation metrics utilized for the report generation task. By including these supplementary evaluations, we aim to provide a more holistic and nuanced understanding of the models' capabilities and limitations in generating free-form textual responses. The results demonstrate that Eyexpert achieved the best performance across both tasks on the supplementary semantic similarity metrics (F1-Rad and BertScore-F1) and the text similarity metrics (BLEU-1 and BLEU-4).

B.2 Eyecare-100K Fine-tuning Ablation Study Results

In the main text, we only present the experimental performance of EyecareGPT-3.8B before and after fine-tuning on Eyecare-100K. The specific evaluation results are shown in the table 12. The results demonstrate a significant improvement in the

Table 9: Sample Question Answering (QA) Templates for Data Conversion.

Question1:

- 1. Is the eye in this picture diseased?.
- 2. Does the eye shown in the image have any disease?
- 3. Is there any sign of illness in the eye in this photo?
- 4. Does this eye image show any signs of abnormalities?
- 5. Does the eye in the image show signs of disease?
- 6. Is there evidence of a disorder in the eye in this picture?
- 7. Are there any visible abnormalities in the eye image?

Positive ;condition;:

- 1. Yes, the eye in the picture has {condition}.
- 2. Yes, the image reveals the presence of {condition} in the eye.
- 3. Yes, the eye shown in this image is impacted by {condition}.
- 4. Yes, this image depicts an eye presenting {condition}.
- 5. Yes, the eye in this image shows evidence of {condition}.
- 6. Yes, the image illustrates an eye with {condition}.

Negative ¡condition¿:

- 1. No, very healthy.
- 2. No, the eye appears healthy in the image.
- 3. No. This image shows that the retina looks normal, with no hemorrhages, exudates or other signs of abnormality.
- 4. No, the eye image appears normal.
- 5. No, the findings from the retinal image suggest a normal and healthy eye.
- 6. No, there are no indications of disease in the image.
- 7. No, the retinal image indicates a healthy eye, with no signs of hemorrhages, exudates, or other pathological changes.
- 8. No significant abnormalities were detected in the eye image.

Question2:

- 1. What ocular disease is evident in this image?
- 2. What eye condition is visible in this picture?
- 3. What condition is affecting the eye shown in the image?
- 4. What issue is apparent in the eye shown here?
- 5. What is wrong with the eye in the image?
- 6. Which disease can be seen in the eye from this picture?
- 7. What health issue is present in the eye in this image?
- 8. What health concern is evident in the eye in this image?
- 9. What problem does the eye shown in the image have?

Positive ;condition;:

- 1. The eye in the image exhibits signs of {condition}.
- 2. {condition} is evident in the eye depicted in the image.
- 3. The image reveals the presence of {condition} in the eye.
- 4. In this picture, the eye appears to be affected by {condition}.
- 5. This image shows an eye with {condition}.
- 6. The eye in the photograph shows signs of {condition}.
- 7. {condition} is visible in the eye from this picture.

Negative ;condition;:

- 1. The eye in this image is very healthy.
- 2. This picture shows a perfectly healthy eye with no signs of disease.
- 3. The eye depicted in the image is completely healthy, showing no illness.
- 4. There is no indication of disease in the eye shown by this image. It's very healthy.
- 5. According to this image, the eye is very healthy and free from any disease.
- 6. The photo indicates a very healthy eye with no presence of disease.

Model	OCT				Fundus			FA-ICGA			СТ			UBM		
	F1-Rad	BLEU-1	BLEU-4	F1-Rad	BLEU-1	BLEU-4	F1-Rad	BLEU-1	BLEU-4	F1-Rad	BLEU-1	BLEU-4	F1-Rad	BLEU-1	BLEU-4	
						(Generalist N	1 odels								
LLaVA-1.5	8.50	11.20	2.18	6.76	28.57	2.44	3.33	3.01	0.26	7.48	6.58	0.86	15.69	12.69	1.93	
Qwen2.5-VL	13.39	22.23	5.06	20.46	36.45	10.21	6.12	11.12	2.44	11.37	16.28	2.95	15.91	6.13	1.15	
InternVL-2.5	12.90	20.06	4.43	16.75	30.09	7.30	4.38	10.49	1.38	9.39	17.02	3.32	17.75	25.34	4.50	
mPLUG-Owl3	10.57	16.63	3.05	21.26	30.02	7.14	4.61	6.37	6.00	10.77	15.99	3.12	19.02	20.70	3.57	
Yi-VL	10.71	17.02	3.24	16.43	19.68	4.31	2.06	6.92	0.59	10.89	11.33	1.77	15.43	17.75	3.05	
MiniCPM-V2.6	14.92	30.48	8.78	19.51	30.76	8.42	6.17	10.18	1.70	12.79	17.21	3.44	20.52	27.64	5.26	
Gemma-3	9.20	23.56	5.50	17.65	32.76	7.54	4.71	6.49	7.00	16.81	24.76	4.47	17.87	25.44	4.88	
							Medical M	odels								
LLaVA-Med	12.36	22.74	4.58	19.44	28.09	7.03	6.45	6.77	0.73	14.34	14.51	2.52	18.89	19.04	3.48	
MedVLM-R1	10.08	22.06	4.30	18.82	28.68	7.18	6.34	7.73	0.77	13.83	14.43	2.94	17.51	24.08	4.31	
HealthGPT-M3	6.64	14.12	2.81	13.28	22.95	6.40	7.20	7.29	1.00	12.39	14.20	2.74	19.12	20.81	4.39	
Eyexpert-3.8B	43.33	48.67	26.20	26.48	37.09	13.16	16.79	40.71	22.23	19.21	18.08	7.00	40.98	53.54	25.21	

Table 10: Other evaluation metrics for the open-ended question answering task in the main experiment.

Table 11: Other evaluation metrics for the report generation task in the main experiment.

	FA	-ICGA			СТ		UBM			
Model	BERTScoreF1	BLEU-1	BLEU-4	BERTScoreF1	BLEU-1	BLEU-4	BERTScoreF1	BLEU-1	BLEU-4	
				Generalist Mod	els					
LLaVA-1.5	81.12	6.06	0.23	82.27	18.44	0.84	81.01	7.57	0.12	
Qwen2.5-VL	84.54	26.81	0.76	84.32	16.28	2.95	81.30	8.01	0.59	
InternVL-2.5	82.21	7.84	0.50	83.15	17.56	0.70	81.98	9.46	0.36	
mPLUG-Owl3	81.12	4.10	0.12	81.89	18.42	0.38	81.52	9.06	0.34	
Yi-VL	80.83	7.70	0.34	83.03	19.33	1.61	80.95	7.45	0.17	
MiniCPM-V2.6	81.77	9.59	0.74	82.61	17.02	1.09	81.36	8.22	0.29	
				Medical Model	s					
LLaVA-Med	81.68	8.86	0.10	81.57	0.05	0.00	81.35	0.09	0.00	
MedVLM-R1	80.76	3.63	0.75	83.12	9.40	6	81.04	7.31	1	
HealthGPT-M3	83.20	10.91	0.59	85.01	27.91	1.82	82.29	11.27	0.29	
Eyexpert-3.8B	90.12	29.41	2.31	88.36	29.22	2.79	85.70	12.97	0.76	

Table 12: Comparative Experimental Results of EyecareGPT-3.8B Before and After Fine-tuning on Eyecare-100K.

Task	Dataset		Eyecare-100K									OmniMedVQA			
		FS	s. s	SL.	ОСТ	Fundus	FA-ICGA	UBN	л с	Т	OCT	Fundus	8		
Closed-QA	-	43.	90 66	5.67	62.48	18.28	78.31	64.7	6 76	.36	81.33	76.69	62.30		
	+ Eyecare-100	K 60.	87 77	7.03	89.76	75.10	91.43	81.6	6 85	.21	100	100	84.56		
Task	Dataset		ОСТ		Fundus		FA-ICGA			UBM		СТ			
		F1-1	Bio Roi	ıge-L I	F1-Bio	Rouge-L	F1-Bio	Rouge	e-L F1-	Bio	Rouge-L	F1-Bio	Rouge-L		
Onen OA	-	55.	53 20	0.10	69.80	23.43	51.79	16.4	7 82	.12	20.08	65.22	13.60		
Open-QA	+ Eyecare-100	K 95.	53 49	9.68	90.79	37.03	86.75	49.2	1 95	.78	47.61	83.90	36.11		
Task	Dataset		FA-ICGA				UBM					СТ			
		Acc ^{GPT}	ScoreAvg	F1-Rad	Rouge-	L Acc ^{GPT}	ScoreAvg	F1-Rad	Rouge-L	AccG	PT Score ^A	^{vg} F1-Rad	Rouge-L		
Report-Gen	-	19.21	75.35	12.78	15.19	4.51	63.41	9.36	12.19	10.7	1 63.93	14.94	18.92		
	+ Eyecare-100K	52.62	85.49	25.04	47.91	58.05	79.90	42.83	57.04	50.0	0 78.64	35.39	42.73		

Table 13: Comparative Experimental Results of EyecareGPT-7B Before and After Fine-tuning on Eyecare-100K.

Task	Dataset		Eyecare-100K									OmniMedVQA			
		F	s. s	SL.	ОСТ	Fundu	s FA-ICG	A UB	м	CT	OCT	Fundus	8		
Closed-QA	-	52.	17 70).33	68.82	77.36	74.71	44.7	78 5	4.93	81.93	77.36	66.93		
	+ Eyecare-100	K 61.	43 77	7.64	90.09	82.25	92.96	86.7	78 8	1.33	99.26	99.56	86.03		
Task	Dataset		ОСТ		Fundus		F	FA-ICGA		UBM		СТ			
		F1-]	Bio Rou	ige-L	F1-Bio	Rouge-	L F1-Bio	o Roug	e-L F1	-Bio	Rouge-L	F1-Bio	Rouge-L		
Onen OA	-	75.	84 25	5.91	80.24	25.12	55.01	17.0)1 8	3.14	23.66	73.17	20.28		
Open-QA	+ Eyecare-100	K 96.	26 50	0.10	90.88	38.13	87.86	51.2	24 9	5.60	47.26	87.27	36.70		
Task	Dataset		FA-ICGA				UBM					СТ			
		Acc ^{GPT}	ScoreAvg	F1-Rad	Rouge-	-L Acc ^G	PT Score ^{Avg}	F1-Rad	Rouge-I	Acc	PT Score ^A	^{vg} F1-Rad	Rouge-L		
Report-Gen	-	25.33	76.02	11.36	12.48	7.27	62.83	12.79	14.88	35.7	1 76.00	15.16	17.38		
	+ Eyecare-100K	53.91	85.97	26.04	48.32	60.0	6 80.56	42.98	58.43	52.4	3 80.71	36.19	43.54		

performance of the fine-tuned EyecareGPT-3.8B across all metrics for each task. Furthermore, the experimental results of EyecareGPT-7B before and after fine-tuning on Eyecare-100K are included in the Appendix, as shown in Table 13. This supplementary data allows for a more comprehensive evaluation of Eyecare-100K's significant value for ophthalmological AI research.

Findings from the results of the two tables reveal the following:

(i) Fine-tuning Significantly Improves Performance: Across nearly all tasks and the majority of datasets, both EyecareGPT-3.8B and EyecareGPT-7B demonstrate a substantial performance increase after fine-tuning on the Eyecare-100K dataset (+ Eyecare-100K). This highlights the effectiveness of the Eyecare-100K dataset in adapting these large language models for ophthalmology-specific tasks.

(ii) Larger Models Generally Perform Better: Comparing the rows before fine-tuning (those not including "+ Eyecare-100K"), EyecareGPT-7B generally exhibits higher initial performance than EyecareGPT-3.8B (Table 10) across most tasks and datasets. This aligns with the common trend that larger language models tend to possess superior zero-shot or few-shot capabilities.

(iii) Fine-tuning Significantly Enhances Large Model Performance: Despite the higher starting baseline of the 7B model, fine-tuning on Eyecare-100K results in similarly significant absolute gains for this larger model. In many instances, the performance level of the fine-tuned EyecareGPT-7B model considerably surpasses that of the fine-tuned 3.8B model.

C Case Study

In this section, we compare the generated answers of our proposed EyecareGPT with those of an open-source medical model (MedVLM-R1) and a closed-source general-purpose model (Claude 3.5). Figures 10, 11, and 12 illustrate the performance of the three models on UBM, CT, and Fundus modalities, respectively, and highlight the differences from the ground truth. Taking Figure 10 as an example, our answer is closer to the true answer, demonstrating EyecareGPT's strong understanding of fine-grained diagnostic questions. Figures 13, 14, and 15 present report generation examples for Fundus, UBM, and CT



Figure 10: A case of UMB 0pen-QA.

Fundus Open-QA	Ground Truth						
Question	This fundus image shows mild nonproliferative diabetic retinopathy with scattered microaneurysms, small dot hemorrhages, and early vascular changes, while maintaining good optic disc appearance.						
	MedVLM-R1 Answer						
	The early diabetic retinopathy changes visible in this fundus image include the presence of new blood vessels (vascular proliferation) and the presence of hemorrhages (bleeding into the retina).						
	Claude Answer The fundus image shows several microaneurysms (small red dots) scattered throughout the retina, which are characteristic early signs of diabetic retinopathy.						
What early diabetic retinopathy changes are visible in this fundus image?	Eyexpert Answer The image shows microaneurysms and dot-and-blot hemorrhages typical of early-stage nonproliferative diabetic retinopathy.						

Figure 11: A case of Fundus 0pen-QA.

modalities. These three figures show that our model can precisely respond to instructions for generating reports, producing well-structured and clear reports that accurately describe abnormal findings in the images.



Figure 12: A case of OTC 0pen-QA.



Please analyze the given medical image of the eve to generate an examination report of that image, including the type of examination. description of the image, and diagnostic recommendations. The report needs to be clear, accurate and use specialized terminology.

MedVLM-R1 Answer

FA Report Generation

Ground Truth

- Examination Type: FA (Fluorescein Fundus Angiography)
- Image Description:
- Punctate microaneurysms in the retinal vasculature · Capillary dilation and leakage observed
- Extensive vitreous hemorrhage significantly obscuring
- fluorescence Extensive hemorrhagic changes present
- Diagnostic Impression: Preliminary diagnosis: Proliferative Diabetic Retinopathy (PDR) with advanced vascular complications • Recommendation: Urgent ophthalmological intervention
- and comprehensive treatment plan required to manage severe retinal changes
- Patient Information: Age: 55 years
 Gender: Female Date of Examination: 2023-04-15 Time of Examination: 13:31:69 Location: Left Eye in Angiography (FA) Image Description: The image depicts a series of retinal angiograms taken from the left eye of a 55-year-old female patient. The examination was performed on April 15, 2023, at 13:31:69. The patient's left eye was examined using Fluorescein Angiography (FA). Image Analysis: 1. FA 0:18.61 (55°) • Type: Fluorescein Angiography • Location: Left Eye Time: 18.61 seconds
 Angle: 55° Image Quality: Clear and detailed
 Significant Features: Clear visualiz on of retinal vessels, including the macula and peripheral 2. FA 0:41.26 (55°) Claude Answer Type of Examination: Description: The images represent a series of fluorescein angiography scans of the retina and choroidal vasculature. The images are organized in three rows, each containing three sequential time-stamped fluorescent angiograms. 1. Vascular Network: Clear visualization of retinal and choroida blood vessels with varying degrees of fluorescence intensity. 2. Top two rows (first two rows) show more distinct vascular branching patterns. 3. Bottom row (third row) demonstrates more diffuse and less defined vascular structures. Notable Findings: Irregular fluorescence patterns in multiple images Possible signs of: Choroidal neovascularization Vascular abnormalities Diagnostic Recommendations: Comprehensive retinal examination Additional imaging studies (OCT, fundus photography)
 Detailed assessment of vascular integrity · Consider systemic health screening for underlying conditions Eyexpert Answer Examination Type: ngiography) Image Description:

Itiple hemorrhages obscuring Capillary dilation and widespread leaka
 Patchy areas of capillary non-perfusion **Diagnostic Impression:** Preliminary diagnosis: Diabetic Retinopathy, Stage 3 (Recommendation: U

Figure 13: A case of FA Report Generation.

Please analyze the given medical image of the eye and generate a report of the examination of that image, including the type of examination, description of the image, and diagnostic recommendations. The report needs to be clear, accurate and use specialized terminology. **UBM** Report Generation

Ground Truth
Examination Type:
Ultrasound Biomicroscopy (UBM)
Imaging Findings:
- Left eye central corneal thickness: 0.53 mm
- Anterior chamber depth: 2.75 mm
- Partial angle opening with some areas of closure
Diagnostic Impression:
- Normal to thin central corneal thickness
- Normal to thin central corneal thickness
- Relatively normal anterior chamber depth
- Mixed pattern of angle configuration suggesting potential risk for angleclosure glaucoma in soocie quadrants
- Recommendation: Clinical correlation, and comprehensive ophthalmic
evaluation advised for complete assessment of glaucoma risk

MedVLM-R1 Answer

Display the second sec

Figure 14: A case of UBM Report Generation.



CT Report Generation

Ground Truth



Ocular CT exam nation

Ocular CT examination **Image description:** • State of the eye: Bilateral eyes are symmetrical with normal size and morphology. Eye orbits are intact with uniform thickness. • Intraocular: Lens is present and normal with homogeneous density. No abnormal density shadows observed within the vitreous humor. • Orbital tissue: Retrobulbar space is Clare. Extraocular muscles and optic nerves show normal morphology and density. No abnormalities detected.

Other findings: Orbital wall bone structure appears normal.
 Diagnostic Recommendations:

Diagnostic Recommendations: • Preliminary diagnosis: No significant abnormalities detected on CT examination of both eyes. • Recommendations: No specific follow-up imaging is indicated.

Correlate with clinical findings if ocular symptoms persist.

MedVLM-R1 Answer

Please analyze the given medical image of the eye to generate an examination report of that image, including the type of examination, description of the image, and diagnostic recommendations. The

report needs to be clear, accurate

and use specialized terminology.

Image Description: The image depicts a set of nine axial CT scans of the head, arranged in a 3x3 grid. Each scan is a cross-sectional view of the brain, showing the skull, brain parenchyma, and surrounding structures. Type of Examination: The examination is a CT scan, which stands for Computed Tomography. CT scans are used to visualize internal structures of the body, including bones, soft tissues, and organs. Description of the Image: Description of the Image: The CT scans show a clear and detailed view of the brain/s anatomy. The skull is visible at the top of each scan, with the brain parenchyma filling the central and lower sections. The scans are arranged in a way that allows for a comprehensive view of the **brain**'s structure from different angles. **Diagnostic Recommendations:** 1. Skull and Bone Structures: <u>Skull and Bone Structures:</u>

 The skull appears intact with no signs of fractures or deformities. The bones show normal density and contour, indicating no significant pathology.
 <u>Brain Parenchyma</u>:
 The brain parenchyma appears homogenous in density, with no signs of hemorrhage, edema, or other abnormalities. The brain appears to be well-vascularized, with no signs of ischemia or infarction.

 abromances. The oran appears to be very vascularized, with no signs of information or inflatcont. • Soft issues • The soft tissues surrounding the brain appear normal, with no signs of inflammation or infection. • The meninges appear intact, with no signs of meningitis or other inflammatory processes. Claude Answer Type of Examination: (Computed Tomography) Scan of the Brain Imaging Technique: Axial and Coronal Non-Contrast CT Slices Axia and Coronal Non-Contrast CT Slices
Description:
- Multiple sequential axial and coronal CT images of the brain
- Multiple sequential axial and coronal CT images of the brain
- Images demonstrate a large, well-defined mass in the posterior fossa, consistent with a cerebellar tumor
- The lesion appears hyperdense and accupies a significant portion of the cerebellar hemisphere
- No obvious calcifications or hemorrhagic components are noted
- Surrounding brain parenchyma shows mild mass effect and slight midline shift
Ventricular system appears compressed but not severely hydrocephalic
Diagnostic Recommendations: 1. Urgent neurosurgical consultation. 2. MRI with contrast for detailed tumor
characterization. 3. Complete neurological examination. 4. Consider biopsy for definitive histopathological diagnosis.
5. Potential differential diagnoses include: - Cerebellar astrocytoma - Medulloblastoma - Metastatic
lesion - Hemanioblastoma lesion • Hemangioblastoma Clinical Implications: Requires immediate comprehensive diagnostic workup Potential surgical intervention may be necessary Further staging and treatment planning needed Recommended Next Steps: Comprehensive neurological assessment
 Advanced neuroimaging
 Multidisciplinary tumor board review Eyexpert Answer Examination type Image description: y. Orbital rims are intact. State of the eye: B Intraocular: Lens is present with regular morphology. Vitreous humor shows uniform density.
 Orbital tissue: Superior and inferior ophthalmic veins and ophthalmic artery show no significant dilation.
 Extraocular muscles and optic nerve demonstrate no definite abnormal changes. Muscle cone and extraconal spaces are clear. Lacrimal glands and lacrimal sac regions show no abnormal density. Eyelids show no definite abnorm changes. Other findings: No additional abnormalities detected.
 Diagnostic Recommendations:
 Preliminary diagnosis: No definite abnormalities detected in bilateral o
 Recommendations: Correlation with clinical findings is recommended.

