Modulated-bath AC calorimetry using modified commercial Peltier-elements

Rolf Lortz, Satoko Abe, Yuxing Wang, Frédéric Bouquet[†], Ulrich Tutsch, and Alain Junod Département de Physique de la Matière Condensée, Université de Genève, 24 Quai Ernest-Ansermet, CH-1211 Geneva 4, Switzerland.

We developed a new type of AC microcalorimeter based on a modulated-bath technique for measuring the specific heat of small microgram samples in the temperature range from 30-300 K, and tested it in magnetic fields up to 14 T. The device is built from a modified commercial Peltier element. The temperature of its top plate can be modulated periodically by Peltier effect, so that the oscillation is symmetrical about the temperature of the main bath. This avoids the problem of DC offsets which plagues conventional AC calorimeters. The sample is attached to a thin thermocouple cross, acting as a weak thermal link to a platform. The absence of a heater reduces the background heat capacity ('addenda') to a minimum. As an illustrative example of the performance of our device, the specific heat in fields up to 14 T of a small single crystal of the high-temperature superconductor Bi_{2.12}Sr_{1.71}Ca_{1.22}Cu_{1.95}O_y has been determined.

[†] Present address: Laboratoire de Physique des Solides, CNRS Université Paris-Sud, 91405 Orsay Cedex, (France).

I. INTRODUCTION

As new materials are often available only in small quantities, the development of devices allowing high-resolution measurements of thermal properties of tiny samples with masses of the order of several micrograms down to a few hundreds of nanograms is of particular interest. The temperature dependence of the specific heat is one of the most basic and informative properties, however calorimetric measurements of samples with such tiny masses require special techniques¹⁻⁷. In this article we present a new modulated-bath AC calorimeter for samples with masses of a few micrograms built from a commercial Peltier element, which presents the advantage of keeping the average temperature of the sample close to that of the thermal bath.

II. METHOD

Besides standard, adiabatic methods where the sample is thermally isolated from its surroundings, dynamical methods with a sample of heat capacity *C* being connected to a thermal bath through a weak thermal link of conductance *k* appear more promising for the measurement of samples with small masses. Alternating current (AC) methods working at a frequency $a \gg \tau^{-1} = k/C$ have the advantage that a thermocouple can be used as a thermal link and thermometer at the same time^{1,2,8}, which strongly reduces the addenda. In addition AC calorimetry allows signal averaging over many heating cycles. This leads to a high resolution for even the tiniest samples. Modulated-bath AC calorimeters^{1,2} using a sine wave heating function have been proven to be extremely useful for this purpose. The sample is attached to a thermocouple cross made of thin wires, while the heater is spatially separated from the sample. In this method the heat is flowing through the thermocouple legs towards the sample and may be more quantitatively controlled than e.g. with optical heating⁶⁻⁸. Due to the absence of a sample

platform the addenda are reduced to a minimum. Heat can be supplied to the sample by an electrically insulated heater on the legs of the thermocouples, measuring the temperature difference between sample and the thermal bath. The use of such a Joule heater as a source of the temperature modulation has the disadvantage that the heat is produced with the square of the supplied heating current. This causes a DC offset which plagues conventional AC calorimeters, since the condition $\omega\tau$ >>1 also means that the AC component of the thermocouple voltage carrying the useful information is much smaller than the DC offset. Alternatively, the temperature modulation can be created using the Peltier effect³. In the latter case heat is supplied on one half of the period of the alternating current, while the sample is cooled on the second half of the period with equal power, thus creating a temperature modulation on the sample symmetrically around the main bath temperature.

III. CONSTRUCTION OF THE CALORIMETER

The basic building block of the calorimeter is a commercially produced element, employing the Peltier effect for cooling various electronic devices⁹. The same elements were used as described in an earlier paper¹⁰, where they served as sensitive heat-flow meters to measure the specific heat and the isothermal magnetocaloric effect of large samples. Here we use them in an inverse way by applying the Peltier effect to modulate the sample temperature to extract the specific heat using an AC technique.

For the construction of the calorimeter the base plate of the Peltier element is soldered with In-Sb eutectic (melting temperature 117° C) onto a copper platform, representing a main thermal bath at a temperature T_B. To avoid damage to the Peltier elements due to differential thermal expansion, a loop of thin (50-100 µm) uncoated copper wire acts as a buffer between the

base plate of the Peltier element and the copper platform. The main bath temperature can be controlled using the temperature regulation of a variable temperature insert in a magnet cryostat with a field up to 16 T. The calorimeter is mounted in high vacuum to reduce heat losses due to thermal conduction and to avoid anomalies due to the condensation/evaporation of residual gases.

The top plate of the Peltier element (3.2 mm x 3.2 mm, see Fig. 1a for a schematic diagram of the element) can be considered as a platform representing a modulated bath. The temperature modulation is driven by the output voltage of a lock-in amplifier connected to the leads of the Peltier element. The top plate of the element is made of sapphire, coated with a thin layer of copper. In the middle of the top plate a cavity of approximately 2 mm diameter and ~ 0.3 mm depth has been ground with a small diamond grinder. The remaining copper coating has been cut with a diamond saw, so that four electrically insulated patches of copper in the corners and an additional patch on one of the sides of the top plate are left (Fig. 1b). A Chromel-Constantan thermocouple (th₁ in Fig. 1b), made of spot-welded 25 μ m \emptyset wire, is either soldered with indium or glued with a silver $epoxy^{11}$ with high thermal conductance to the additional patch on the top plate of the Peltier element. It enables to monitor the oscillating part of the temperature difference between the top plate and the main bath (T_{P-B}) . On the remaining copper patches in the corners of the plate the four legs of a Chromel-Constantan thermocouple cross - made of the same wire as th_1 , but flattened after spot welding (th_2 in Fig. 1b) - have been thermally anchored in the same way. Two legs thermally anchored to the copper patches are continued by electrical leads made of a pair of phosphor-bronze wires with the same diameter. They allow to measure the difference between the sample temperature and the modulated-bath temperature (T_{S-P}) . The second pair of legs are extended using the same thermocouple material to heat-sink chips on the main bath, to

measure the difference between the oscillating sample temperature and the (constant) temperature of the main thermal bath (T_{S-B}).

The sample is glued onto the thermocouple cross using a drop of highly diluted VGE7031 varnish¹². To do so, a drop of varnish is deposited into the hole on the top plate using a syringe so that the liquid fills the space under the thermocouple cross. The surface tension of the liquid helps to gently push the sample onto the thermocouple cross, where it remains stuck after the varnish dried. To remove the sample it is sufficient to put one or two drops of ethanol into the hole to wash away the varnish.

For the specific-heat measurement an AC voltage signal of typically 0.5-2 V is set at the output of the lock-in amplifier connected to the voltage leads of the Peltier element during a temperature sweep (typically 0.5-1 K/min) of the variable temperature insert. (Note that the current is determined essentially by the output resistance of the current source, 50 Ω , and that of the connecting wires, about 30 Ω ; the resistance of the Peltier element itself is 1-2 Ω). The temperature of the main bath T_B is measured using a calibrated Cernox sensor¹³, mounted just besides the Peltier element on the copper block. The signals of the three thermocouples (T_{P-B} , T_{S-B} and T_{S-P}) are amplified 10⁴ times by a picovolt DC amplifier¹⁴, when the frequency is low enough, and fed into a digital lock-in amplifier¹⁵.

IV. PERFORMANCE

The configuration of the thermal anchoring of the four legs of the thermocouples allows us to measure the specific heat in two different ways:

A) Modulated-bath AC calorimetry:

If the pair of legs being anchored on the Peltier element and measuring the difference between the temperature of the top plate of the Peltier element and that of the sample (T_{S-P}) is used, the method is the standard modulated-bath AC calorimetry¹. The particularity of our system is that the Joule-heater is replaced by a Peltier element with the advantages mentioned above. The thermocouple representing T_{S-P} measures the difference between the oscillating component of the sample temperature relative to the main bath $T_{S-B}(t)$ and that of the modulated bath relative to the main bath $T_{P-B}(t)$: $T_{S-P}=T_{S-B}(t)-T_{P-B}(t)$, where (see Fig. 2a for a schematic diagram). The system follows the differential equation (1):

(1)

$$k(T_{P-B} - T_{S-B}) = C \frac{dT_{S-B}}{dt}$$
$$\tau \left(\frac{dT_{S-P}}{dt} + \frac{dT_{P-B}}{dt}\right) + T_{S-P} = 0$$

with the harmonic solutions (2) and (3):

$$\frac{|T_{S-P}|}{|T_{P-B}|} = \sqrt{\frac{1}{1 + \frac{1}{\omega^2 \tau^2}}}$$
(2)

$$C = \frac{k}{\omega} \sqrt{\frac{1}{\frac{|T_{P-B}|^2}{|T_{S-P}|^2} - 1}}$$
(3)

where k is the heat conductance of the thermocouples between sample and top plate of the Peltier element, ω the angular frequency of the temperature modulation, C the heat capacity of the sample (including contributions from the addenda), and $\tau = k/C$ the time constant. Equation (2) allows us to determine the frequency dependence of the amplitude of the temperature modulation of the sample normalised by T_P . Below the cutoff frequency $\omega_0 = 1/\tau$ the sample temperature T_S asymptotically approaches the temperature modulation of the Peltier element T_P and thus the signal on the thermocouple T_{S-P} is small. In this frequency range it is possible to extract the heat capacity using equation (3) from the signal of the first harmonic using the lock-in amplifier (and not the second harmonic as would be the case for AC calorimetry with Joule heating). This method allows us thus to measure the specific heat at frequencies below the cutoff frequency, while in standard AC calorimetry the specific heat enters the signal above the cutoff frequency. This might be of advantage in the presence of additional time constants due to low thermal conductivity of the sample or imperfect thermal coupling between the sample and the thermocouple. In addition, the low frequency range leaves the choice between DC preamplifiers and transformers open to optimize the signal/noise ratio. The drawback is that the sample temperature is modulated with a large amplitude close to that of the modulated bath, while the amplitude of the signal $|T_{S-P}|$ is only a fraction of that. Indeed, to ensure a good resolution, a sufficiently large amplitude $|T_{S-B}|$ has to be chosen, which may broaden specific-heat anomalies at phase transitions. The working frequency during a measurement should thus be chosen rather close to the cutoff frequency where $|T_{S-P}|$ is of the same order of magnitude as $|T_{S-B}|$.

B) Standard AC calorimetry:

If the other pair of legs of the thermocouple cross is used, the difference between the oscillating component of the sample temperature and the temperature of the main bath T_{S-B} is measured. This method resembles standard AC calorimetry¹⁶, with the only difference that the source of the temperature modulation is not located on the sample but on the legs of the thermocouples (see Fig. 2b for a schematic diagram). In this case the solution of the differential equation describing the system is given by equations (4) and (5):

$$\frac{|T_{S-B}|}{|T_{P-B}|} = \sqrt{\frac{1}{1+\omega^2 \tau^2}}$$
(4)

$$C = \frac{k}{\omega} \sqrt{\frac{\left|T_{P-B}\right|^2}{\left|T_{S-B}\right|^2} - 1}$$
(5)

Equation (4) represents the frequency dependence of the amplitude of the temperature modulation of the sample normalised by T_{P-B} . The signal $|T_{S-B}|$ is maximal at low frequencies where it asymptotically approaches $|T_{P-B}|$, whereas it decreases above the cutoff frequency. The heat capacity can be extracted at frequencies above the cutoff frequency using equation (5). The advantage of this method is that the temperature modulation of the sample can be measured at higher frequencies, allowing a faster temperature-sweep rate to be chosen during the experiment, while still averaging over many cycles of the modulation with the lock-in technique. But most importantly, the measured signal T_{S-B} represents directly the real amplitude of the modulation of the sample temperature. Thus, no resolution is lost when the signal is measured at frequencies far

above the cutoff frequency with a sufficiently low amplitude to avoid broadening of anomalies at phase transitions.

Fig. 3a shows the frequency dependence of the amplitude of T_{S-P} , T_{S-B} and T_{P-B} for a sample with a mass of 26 µg at T=90 K. The dependence of $|T_{P-B}|$ follows the behaviour of a system which can be described by formula (6) with C_P the heat capacity of the top plate of the Peltier element, k_P the heat conductance of the thermopiles between top and base plate, and $\tau_P=C_P/k_P$:

$$\left|T_{P-B}\right| = T_0 \sqrt{\frac{1}{1 + \omega^2 \tau_P^2}} \tag{6}$$

By fitting the data to this formula, we find $\tau_P=0.32$ s for the time constant of the Peltier element, representative for the frequency dependence of the heating/cooling power creating the temperature modulation. Fig. 3b shows the frequency dependence of $T_{S\cdot P}/T_{P\cdot B}$ and $T_{S\cdot B}/T_{P\cdot B}$ together with a fit using formula (2) and (4). The data has been measured without preamplifiers, as their response might add an additional pole in the frequency dependence. The theoretical curves are in good agreement with the data in the relevant frequency range below ~2 Hz. The frequency corresponding to the crossing point of the two fitting functions determines the main time constant (in our case for a 26 µg Bi_{2.12}Sr_{1.71}Ca_{1.22}Cu_{1.95}O_y sample at 90 K: $\tau=0.88$ s, slightly above that of the Peltier element). This ensures that sufficient power to create a temperature modulation of reasonable amplitude is available over the possible range of working frequencies (~0.05-3 Hz). The upturn in the modulation above 4 Hz in Fig. 3a is associated with the internal time constant of the Bi-Te material in the Peltier cell. Fig. 3c shows T_{S-P} , T_{S-B} and T_{P-B} as a function of amplitude of the driving AC voltage. The dependence is in good approximation linear, as expected for reversible Peltier heating $P(t)=\pi I(t)$, where π is the Peltier coefficient and I(t) the driving current. To test if there is any additional effect due to Joule heating $P(t)=RI^2(t)$ from the internal resistance of the Peltier element (R=1-2 Ω at room temperature), we also measured the response of T_{P-B} and T_{S-B} in the second harmonic. The signal at T=100 K, using a driving AC voltage of $U_{eff}=1.8$ V (i.e. about 20 mA) and a frequency of 0.555 Hz, was smaller than 5 mK so that the effect of Joule heating can be neglected. The typical DC offset found between the temperature of the sample and that of the main bath is smaller than 0.1 K for all temperatures, as long as the temperature sweep rate during the experiment is chosen smaller than 1 K/min. Thus the DC voltage does not saturate the input of the amplifiers, and the average temperature of the sample is essentially that of the main bath.

The design of the calorimeter allows thus to measure the heat capacity with two different methods, so that, dependent on the value of the cutoff frequency (directly related to the heat capacity of the sample and the heat conductance of the thermocouple legs), it can be considered if it is more favourable to measure at relatively higher or lower frequencies. This might be important as high sensitivity DC preamplifiers can only be used at low frequencies (typically < 1Hz) whereas transformers only work at high frequencies (typically >10 Hz).

V. CALIBRATION OF THE CALORIMETER

For the calculation of the absolute value of the specific heat the heat conductance k of the thermocouples and a possible contribution from the addenda (thermocouple cross below the sample, glue) have to be determined independently. A first estimation for k of a Chromel/Constantan thermocouple can be obtained from literature data¹. A precise method to

estimate the relevant length of the thermocouple legs is to measure them on an enlarged photograph of the calorimeter with the sample mounted. We use the lengths from the edge of the sample, covering the thermocouple cross, to the top plate of the Peltier element. For an accurate calibration we used two silver samples with different masses (85 μ g and 40 μ g) of high purity. Using the difference of the heat capacity of both samples allowed us to cancel the contribution of the addenda and to extract *k* by comparing the specific heat with literature data¹⁷. After determination of *k* the heat capacity of one sample can be used to extract the additive contribution from the addenda (see Fig. 4). We note that the amount of glue used to mount the sample might be difficult to estimate and remains a factor of uncertainty in the determination of the absolute value, especially for samples with very small masses.

The sensitivity of the thermocouples is usually field dependent, but cancels in the formulas of the specific heat as only quotients of the signals of two thermocouples of the same kind enter. We tested this by measuring the heat capacity of an Ag sample (which is not field dependent) in a slow temperature sweep upon changing the field stepwise up to 14 T. No anomalies due to the field changes appeared in the final specific-heat data. As an upper limit, we estimate an error in the order of 1 % due to magnetic fields smaller than 14 T.

VI. EXAMPLES OF MEASUREMENTS, RESOLUTION AND ACCURACY

As an illustrative example of the performance of our device, we have chosen a single crystal of the high-temperature superconductor $Bi_{2.12}Sr_{1.71}Ca_{1.22}Cu_{1.95}O_y$ (Bi-2212) close to optimally doping with a mass of 26 µg. Details of sample growth and treatment can be found elsewhere^{18,19}. Fig. 5a shows the temperature dependence of amplitude $1/|T_{ac}|$ of the thermocouple signals T_{S-B} and T_{P-B} , allowing to extract the specific heat from formula 5. We have

chosen a working frequency of v=0.89 Hz, nearly 5 times the cutoff frequency (0.18 Hz) at 90 K. Fig. 5b shows the corresponding temperature dependence of the amplitude $1/|T_{ac}|$ of the thermocouple signals $T_{S\cdot P}$ and $T_{P\cdot B}$ for the use in the modulated-bath mode to calculate the specific heat using formula (3). Here the working frequency has been chosen as v=0.1 Hz, which is about half the cutoff frequency. We decided to use the standard AC-calorimetry method, using thermocouple $T_{S\cdot B}$, as this allows us to work at a higher frequency and thus to perform the experiment during a rather fast temperature sweep (0.5 K/min). The working frequency of v=0.89 Hz was still low enough to use DC preamplifiers. To obtain maximum resolution, we smoothed the temperature dependence $|T_{P\cdot B}(T)|/|T_{S\cdot B}(T)|$. The result in zero field and in fields up to 14 T is presented in Fig. 6. The absolute values of the specific-heat curves, as well as the width and shape of the superconducting transition at $T_c=$ 90 K is in reasonably good agreement with previously published specific-heat data, obtained on a ~10'000 times larger single crystal of Bi-2212 (0.2545 g) of the same batch using an adiabatic calorimeter¹⁸.

We note that the size of the anomaly at the superconducting transition of Bi-2212 is very small compared to the large phonon background. Using the difference of specific heat in zero field and in 14 T (see insert in Fig. 6) to estimate the size of the specific-heat anomaly at T_c , we find that the anomaly is only ~0.5 % of the total specific-heat value, which can still be resolved with good resolution in this type of calorimeter. For the measurements we set the time constant of the lock-in amplifier to one second and recorded data during a temperature sweep of ~5 mK/s. To characterize the noise in the measurement we subtracted a smooth polynomial from the structureless 14 T data of Fig. 6. The noise level as characterised by the root mean square is of

the order of 2.6 x 10^{-5} J/(K²gat). The limit of resolution due to noise is thus 2.3 x 10^{-2} % of the total specific heat at *T*=90 K (see Fig. 7), 27 times smaller than the anomaly at *T_c*.

VII. CONCLUDING REMARKS

In the following we would like to compare the performance of our device to other types of calorimeters. The measurements show that our calorimeter design allows us to measure the specific heat of samples of tiny masses of a few tens of micrograms, or even down to a few micrograms, with a sensitivity being comparable to standard adiabatic methods. Clearly, when measuring samples with such tiny masses, the absolute value of specific heat is plagued by errors due to uncertainties in the contributions from addenda and heat conductivity of the thermal link to the surrounding thermal bath (arising mainly from uncertainties in the relevant lengths of the thermocouple legs). Best absolute values are obtained with adiabatic methods, but require samples with masses >1 $g^{17,20,21}$ or >0.1 g^{18} . Relaxation methods^{4,22-24} can be used for samples with masses down to $\sim 100 \,\mu$ g, the limit of miniaturization is given by the fact that a thermometer and thus a sample platform is needed. Alternating methods¹⁶ allow the use of thermocouples, so that the addenda can be reduced to a minimum when using chopped infrared light from a laser diode $^{6-8}$ or a modulated-bath technique 1,2 to modulate the sample temperature. This and the high sensitivity due to averaging over many cycles of the temperature modulation is perfectly suited for the construction of miniature calorimeters.

The calorimeter presented in this article is highly specialized for the use in the temperature range between 30 to 300 K and for sample masses from a few μ g up to a few hundreds of μ g, when the main interest is to study contributions from phase transitions. For the use at lower temperatures the calorimeter is not suitable, as the Peltier-heating/cooling power

goes to zero faster than linearly with temperature approaching zero. A similar design, but with a modulated bath based on Joule heating in combination with low-temperature thermocouples such as <u>Au</u>Fe/Chromel may be used at lower temperatures, but then the advantage of having no DC offset in the measured signal is lost. With AC calorimetry it is in addition possible to measure the field dependence of the specific heat during a field sweep at stabilized temperature with a high density of data.

We note that AC calorimetry may underestimate latent heats at first-order transitions when the amplitude of the temperature modulation is kept smaller than the temperature hysteresis. We also note that larger and smaller Peltier elements are commercially available⁹, which can be used in combination with thicker thermocouple wires to construct calorimeters for larger samples in the same way (see Fig. 1c).

- ¹ J.E. Graebner, Rev. Sci. Instrum. **60**, 1132 (1989).
- ² S. Blanchard, *Transitions de Phases et Dynamique des Vortex dans (K, Ba)BiO₃*, PhD-thesis, Laboratoire d'Etudes des Propriétés Electroniques des Solides CNRS Grenoble, France (2003).
- ³ D.H. Jung, I.K. Moon, and Y.H. Jeong, cond-mat/0110304 (2001).
- ⁴ S.G. Doettinger-Zech, M. Uhl, D.L. Sisson, and A. Kapitulnik, Rev. Sci. Instrum. **72**, 2398 (2001).
- ⁵ R. Bachmann, F.J. DiSalvo, Jr., T.H. Geballe, R.L. Greene, R.E. Howard, C.N. King, H.C. Kirsch, K.N. Lee, R.E. Schwall, H.-U.Thomas, and R.B. Zubeck, Rev. Sci. Instrum. 43, 205 (1972).
- ⁶ M.B. Salamon, S.E. Inderhees, J.P. Rice, B.G. Pazol, D.M. Ginsberg, and N. Goldenfeld, Phys. Rev. B 38, 885 (1988).
- ⁷ A. Carrington, C. Marcenat, and F. Bouquet, Phys. Rev. B 55, R8674 (1997).
- ⁸ M.B. Salamon, Phys. Rev. B **2**, 214 (1970).
- ⁹ Thermion Company, 65009 Odessa, Ukraine.
- ¹⁰ T. Plackowski, Y. Wang, and A. Junod, Rev. Sci. Instrum. **73**, 2755 (2002).
- ¹¹ Model Epo-Tek H31LV, Epoxy Technology, 14 Fortune Drive, Billerica, MA 01821, U.S.A.
- ¹² Model VGE7031, Insulating varnish and adhesive, LakeShore Cryotronics, Inc., 575 McCorkle,
 Blvd. Westerville, Ohio 43082, U.S.A.
- ¹³ Model CX-1050-SD, LakeShore Cryotronics, Inc., 575 McCorkle Blvd.-Westerville, Ohio 43082,
 U.S.A.
- ¹⁴ Model A20 DC Nanovolt amplifier, EM Electronics, EM The Rise, Brockenhurst, Hampshire, SO427SJ, England. These amplifiers are free from 1/f noise.
- ¹⁵ Model SR830 DSP lock-in amplifier, Stanford Research Systems, Inc., 1290-D Reamwood Avenue, Sunnyvale, CA 94089.
- ¹⁶ P.F. Sullivan and G. Seidel, Phys. Rev. **173**, 679 (1968).
- ¹⁷ D.L.Martin, Can. J. Phys **65**, 1104 (1987).

- ¹⁸ A. Junod K.-Q. Wang, T. Tsukamoto, G. Triscone, B. Revaz, E. Walker, and J. Muller, Physica C 229, 209 (1994).
- ¹⁹ T. Tsukamoto, G. Triscone, J.-Y. Genoud, K.-Q. Wang, E. Janod, A. Junod, and J. Muller, J. Alloys and Compounds **209**, 225 (1994).
- ²⁰ T. Atake, Y. Takagi, T. Nakamura, and Y. Saito, Phys. Rev. B **37**, 522 (1988).
- ²¹ N. E. Phillips, R.A. Fischer, S.E. Lacy, C. Marcenat, J.A. Olsen, W.K. Ham, A.M. Stacy, J.E. Gordon, and M.L. Tan, Physica B & C 148, 360 (1987).
- ²² D. Sanchez, A. Junod, J.-Y- Genoud, T. Graf, and J. Muller, Physica C 200, 1 (1992).
- ²³ R. Bachmann, H.C. Kirsch, and T.H. Geballe, Rev. Sci. Instrum. **41**, 546 (1970).
- ²⁴ R. Bachmann, H.C. Kirsch and T.H. Geballe, Solid State Commun. 9, 57 (1971).

Figure captions

FIG. 1. (a) Schematic diagram of the Peltier element used to modulate the sample temperature. (b) Schematic diagram of the Peltier AC calorimeter (view from top). (c) Platform holding two AC calorimeters with their "bath" modulated by Peltier effect. The smaller one, $3.2 \times 3.2 \text{ mm}$ (as used for the measurements presented in this article), uses $25 \text{ }\mu\text{m} \text{ }\emptyset$ thermocouples, whereas the larger one, $6.4 \times 6.4 \text{ }\text{mm}$, uses $50 \text{ }\mu\text{m} \text{ }\emptyset$ thermocouples.

FIG. 2. (a) Schematic diagram of a modulated-bath calorimeter, where the shaded area represents the top plate of the Peltier element which can be considered as a modulated thermal bath of temperature $T_B+T_P(t)$, with T_B the temperature of the main thermal bath. The heat conductance k is given by the legs of the thermocouple cross in our device.

(b) Schematic diagram in the usual AC mode. The heater in our device is the top plate of the Peltier element, coupled to the sample by the legs of the thermocouple cross.

FIG. 3. (a) Frequency dependence of the amplitude $|T_{ac}|$ of thermocouples T_{S-B} , T_{S-P} and T_{P-B} for a 26 µg Bi-2212 sample at 90 K (measured without A20 preamplifier). (b) Frequency dependence of the amplitude $|T_{ac}|$ of thermocouples T_{S-B} and T_{S-P} divided by that of T_{P-B} . The lines are oneparameter fits of the data using formulas (3) and (5). The crossing point of the two lines represents the cutoff frequency ν =0.18 Hz (measured without A20 preamplifier). The data at high frequencies of the thermocouple T_{S-P} is below the fit, possibly due to an asymmetry of heat flow through the T_{S-B} and T_{S-P} thermocouples. (c) Dependence of the amplitude $|T_{ac}|$ of thermocouples T_{S-B} , T_{S-P} and T_{P-B} on the amplitude of the AC voltage driving the temperature modulation of the top plate of the Peltier element for a 26 μ g Bi-2212 sample at 90 K and ν =0.89 Hz. We note that the signals are not in phase.

FIG. 4. Specific heat of a pure Ag sample (40 μ g) with addenda compared to literature data¹⁷. The difference allows the estimation of the addenda due to the thermocouple cross and the glue used to mount the sample.

FIG. 5. (a) Temperature dependence of the amplitude $1/|T_{ac}|$ of the signals of thermocouples T_{S-B} and T_{P-B} taken at $\nu=0.89$ Hz above the cutoff frequency at 90 K (~0.18 Hz) (standard AC calorimetry). (b) Temperature dependence of the amplitude $1/|T_{ac}|$ of the signals of thermocouples T_{S-P} and T_{P-B} taken at $\nu=0.1$ Hz slightly below the cutoff frequency at 90 K (modulated-bath AC calorimetry).

FIG. 6. Absolute value of the specific heat C/T of a single crystal of $Bi_{2.12}Sr_{1.71}Ca_{1.22}Cu_{1.95}O_y$ with a mass of 26 µg in fields of 0 T, 1 T, 2 T, 8 T and 14 T. The anomaly at 90 K is the superconducting transition. In the insert the smoothed 14 T data has been subtracted as a background.

FIG. 7. Noise in the measurement after subtracting a smooth polynomial from the structureless 14 T data in Fig. 6. The root-mean-square value of the noise is $2.6 \times 10^{-5} \text{ J/K}^2$ gat which is 0.023 % of the absolute specific-heat value.





(c)



(a)

(b)











