

GDF - A GENERAL DATAFORMAT FOR BIOSIGNALS

VERSION 2.00

Alois Schlögl¹
<alois.schloegl@tugraz.at>

¹Institute for Human-Computer Interfaces, University of Technology Graz
(2004-2006)

Keyword(s): biomedical data format, ISO 9000, quality management, general public license GPL, ISO 11073, IEEE 1073, CEN/TC251

ABSTRACT

Biomedical signals are stored in many different data formats. Most formats have been developed for a specific purpose of a specialized community for ECG research, EEG analysis, sleep research, etc. However, in the field of Brain-Computer interfaces (BCI), no special format is available. BCI research requires not only the raw EEG data, but also information of cues and additional timing information. This additional information is similar to event and marker information, used in other biomedical signal data; however each domain uses its own coding scheme.

In order to address the needs of BCI research and to unify the various needs of other biomedical signal processing fields, a so-called “General Data Format for biomedical signals” (GDF) is developed. The GDF format is described in detail and aims to be a truly “open standard”.

1. INTRODUCTION

Biosignals are currently stored in a wide variety of different data-formats. BioSig – an open source software library – supports approximately 50 different data formats [1]. Many data formats have been developed by commercial companies, research groups and standardization organizations. Furthermore, some scientific data formats like HDF, netCDF and CDF should be mentioned, however these do not support the needs for storing biomedical signals. An overview is provided at [2]. Recently, also the use of XML is discussed, for example XML is used in HL7 for storing annotated ECG data. It is beyond the scope of this paper to discuss the various data formats; however the BioSig project [1] supports more data formats (currently ca. 50) than any other software project in the field. Support for these file formats has been implemented by the author of this work.

Most of the data formats addressed just some specific needs, but none was flexible enough to address new needs in the field of EEG recordings and Brain-Computer Interface research. Therefore, version 1 of the General Data format (GDF) [3] has been developed, and was successfully implemented and used in BCI research. It provides a common coding scheme for events, and supports many useful features (different sampling rates and calibration values for different channels, an automated overflow detection, support of different data types, encoding of filter settings etc.), that are only partly implemented in other formats. Within the recent years, new requirements became apparent. In this report, Version 2 of the GDF is described. Briefly, it addresses the need for

- Subject specific information (gender, age, impairment, etc)
- Recording location, identification of recording software, etc.
- Possibilities for storing the electrode positions in spatial coordinates, electrode impedances, etc.
- More efficient encoding of date and time, physical dimensions, filter information
- Non-equidistant (sparse) sampling

The structure of GDF v2.0 is an extension of the EDF [4], GDF1.x [3] and EDF+ [5], all of them have a fixed header with 256 bytes and a variable header with 256 times “number_of_channels”. Moreover, GDF v2.0 supports storing events (annotations and markers) and non-equidistant (sparsely) sampled data. Based in the version identification field (first 8 bytes), it is easy to write software supporting all formats simultaneously. This is also demonstrated by the BioSig project [1], which provides the prototype implementation for Octave/Matlab as well as C/C++. Support for other programming languages (like Python and Java) is under consideration.

2. SPECIFICATION

The general data structure is defined in Table 7 (see Annex), other tables define specific fields in more detail. The issues were considered during the design of the data format.

1. Data is stored in little Endian format. However, BioSig [1] supports also big-endian platforms by converting the data internally.
2. The Version field is of type char[8] and is stored at the beginning of the file. The first three letters are always 'GDF'. During the testing period of GDF writing functions, 'GDF 1.99' should be used. The final version will contain the value "GDF 2.00". This field is used to provide upwards compatibility with past and future versions of GDF.

3. The patient identification field is 66 bytes long (instead of 80 bytes in EDF [4], EDF+ [5] and GDF 1.x [3]). This fact is justified, because in EDF and GDF 1.x rarely all 80 bytes are used. In EDF+ [5] 14 bytes of this field are used for encoding Gender and Birthday; GDF 2.0 stores this information elsewhere (the gender is encoded in byte 88 and the Birthday is binary encoded in bytes 177-184). Therefore, the PID field is reduced by 14 bytes. The remaining 66 bytes contain the patient identification code, the patient name, the patient classification. The use of any remaining bytes is not specified. Each field is separated by and must not contain a space (ASCII 32) character. Empty fields are encoded by the single letter “X”; For reasons of privacy protection, we recommend not to include the real patient name.
4. The 85th byte (see Table 1) can be used to describe smoking, alcohol abuse, drug abuse and/or medication.

Table 1: Encoding of Byte 85

Byte 85	bit 7	bit 6	bit 5	bit 4	bit 3	bit 2	bit 1	bit 0
Smoking								
unknown							0	0
NO							0	1
YES							1	0
Alcohol abuse								
unknown					0	0		
NO					0	1		
YES					1	0		
Drug abuse								
unknown			0	0				
NO			0	1				
YES			1	0				
Medication								
unknown	0	0						
NO	0	1						
YES	1	0						

5. The weight field contains the body weight in kilograms [kg]; a value of 0 indicates an unknown weight, a value of 255 indicates a weight larger then 254 kg.
6. The height field contains the body height in centimeters [cm]; a value of 0 indicates an unknown height, a value of 255 indicates a height larger then 254 kg.
7. The gender, handedness and visual impairment of the subject is stored the 88th byte according to table 2.
8. The field recording identification (RID, starting at byte 89) has been reduced from a length of 80 to 68 bytes. This is justified because EDF+ uses 12 bytes to store the recording date; in GDF this information is stored in bytes 169-172. If byte 156 is zero, bytes 153-156 are used by the RFC1876 header (see next issue) and must not be used by the RID field.
9. The recording location (place of recording) is stored according the RFC1876 using 4 uint32 numbers (bytes 153 to 168 in the GDF header, see Table 3). If the version number (byte 156) is not zero, bytes 153-156 are used by the field “recording identification” and do not belong to the location information.

Table 2: Encoding of Byte 88.

Byte 88	bit 5	bit 4	bit 3	bit 2	bit 1	bit 0
Gender						
unknown					0	0
male					0	1
Female					1	0
Handedness						
Unknown/undefined			0	0		
Right			0	1		
left			1	0		
equal			1	1		
Visual impairment						
Unknown/undefined	0	0				
No impairment	0	1				
Visual impairment	1	0				
Impairment but Corrected	1	1				

Table 3: Encoding of the recording location

GDF Header bytes		LSB			MSB
153-156		Vertical Precision	Horizontal Precision	size	Version
157-160	Latitude [1/3600000 degrees]				
161-164	Longitude [1/3600000 degrees]				
165-168	Altitude [in cm]				

10. The date/time format of the begin of data recording is changed in Version 2.0 from an 16 byte ASCII format to an 8-byte (64 bit) binary format. The 4 least significant bytes indicate the fraction of a day (2^{32} corresponds to $3\,600 \times 24 = 86\,400$ s or 1 day) and provides a time resolution of approx 20 μ s. The 4 most significant bytes count the number of the days since 1. Jan 0000. Jan 1st, 1970 is day number 719 529. Conversion of this integer into a double (with 52 bits mantissa) does not result in round-off errors for the next 800 years. Only in the year 2800, the last bit will be lost which results in a time resolution of 40 μ s instead of 20 μ s for another 2800 years. Accordingly, this date/time format sufficiently accurate for storing the starting time and birthday information. The starting time of the recording is stored n bytes 169-176, the birthday of the subject is stored in bytes 177-184 for obtaining the subjects age at recording time. However, for privacy protection, the birthday could have been modified up to 1

year. The date-time conversion can be implemented in the following way:

```
/* conversions between t_time and gdf_time */

gdf_time = (int64) round((t_time/(3600*24) + 719529) * 2^32);
t_time = ((gdftime / 2.0^32) - 719529) * (3600*24);
```

11. The length of the header is defined in the field "number of blocks (each 256 bytes) in the header record". The number of blocks must be at least (1+NS). The number of data blocks beyond (NS+1) can be used to extend the header size. Currently, this feature is discouraged but could be used in future versions of the GDF format.
12. Equipment Provider Identification is defined as 8 byte code. Application for codes should be addressed to the author.
13. Technician ID and Lab ID have been removed in GDF v2.0. Instead, the IP address of the recording equipment is included in the header. The remaining bytes are reserved for future needs. Accordingly, the fixed header contains 30 reserved bytes in V2.0.
14. The size of the subject's head (circumference, distance nasion – inion, left to right mastoid) is stored in three 16-bit integer numbers in the bytes 207-212. The distances are stored in the units millimeter [mm]. Value 0 indicates an unknown value.
15. A "record" is a data segment with a fixed length as indicated in the field "Duration". The number of data records is of type int64 (8byte); the duration of data records in seconds is a rational number consisting of the numerator (uint32) and a denominator (uint32) (together 8byte). The number of channels is of type uint32 (4byte).
16. The field "NS" indicates the number of channels.

*Table 4: Examples of Physical Units. The full table is described in Annex A of the FEF standard and is available from http://biosig.cvs.sourceforge.net/*checkout*/biosig/biosig/doc/units.csv [1, 11(p. 62-75, Table A.6.3), 12].*

Code	Units	Description
0	Unknown/undefined	
512	-	Dimensionless
544	%	Per cent
736	degree	
768	rad	
2496	Hz	
3872	mmHg	
4256	V	
4384	K	Kelvin
6048	°C	degree Celsius

17. The "label"-field contains the description for each channel. It is recommended that the labels follow existing standards (e.g. [11,12]) whenever possible. The transducer type describes the sensor type. There exist already some standards for encoding of the channels and transducers (e.g. EN1064 [10] and [11,12]). Currently, the encoding scheme is not complete. Therefore, these fields are ASCII-encoded. However, this might change in future.
18. The Physical Dimensions are encoded with a 16 bit integer number according to The Annex of the File Exchange Format (FEF) for Vital signs [11,12]. The Physical Dimensions are encoded according to [11] (Table A.6.1, p.52-53, Table A.6.3, p.54-61) and [12] (Table A.4.3 p. 38- 44).

Some examples are listed in Table 4. The decimal factors are encoded as offsets (using the 5 least significant bits) according to Table 5. For example, the physical dimension “ μV ” (micro-Volt) is encoded as 4256 (Volt) + 19 (micro) = 4275.

Table 5: Table of Decimal Factors. The Table is also available from http://biosig.cvs.sourceforge.net/*checkout*/biosig/biosig/doc/DecimalFactors.txt [1, 11 (p. 60, Table A.6.1), 12].

Name	Magnitude	Code Offset
yotta	10^{+24}	10
zetta	10^{+21}	9
exa	10^{+18}	8
peta	10^{+15}	7
tera	10^{+12}	6
giga	10^{+9}	5
mega	10^{+6}	4
kilo	10^{+3}	3
hecto	10^{+2}	2
deca	10^{+1}	1
	10^{+0}	0
deci	10^{-1}	16
centi	10^{-2}	17
milli	10^{-3}	18
micro	10^{-6}	19
nano	10^{-9}	20
pico	10^{-12}	21
femto	10^{-15}	22
atto	10^{-18}	23
zepto	10^{-21}	24
yocto	10^{-24}	25

19. Physical minimum and maximum are stored as IEEE floating point numbers with double precision (8 bytes) instead of char[8] (=8*ASCII).
20. Digital minimum and maximum are of type "float64". This field as 52 bits significant, therefore all integer numbers up to 52 bits are stored in its full precision. In order to support quality control and support an automated overflow detection [6], the value of digital minimum and maximum should indicate the overflow (saturation) values of the recording system.
21. The field "pre-filtering" is in GDF Version 2.0 only 68 bytes (it was 80 bytes in EDF, BDF and GDF Version <1.9) and will become obsolete in future. Its purpose is replaced by the filter settings.
22. The filter settings are stored as floating point numbers in the variable header. Low pass, High pass and Notch are stored as a 32-bit floating point number. Unknown values are indicated by NaN's; Notch Off is indicated by a number smaller than zero (typically -1).
23. The "Samples Per Record"-field is uint32 (4bytes)*NS. The "Samples Per Record" is the sampling rate of the corresponding channel multiplied by the duration of a record. A value of zero (samples per record), indicates a channel with sparse (non-equidistant) sampling. In this case, the actual sampling values are stored in the event table.
24. Various data types are defined (see Table 6). The type information is stored in the variable header as uint32 * NS (NS = number of channels) after "number of samples per record".

Table 6: Data types, memory requirement, coding scheme. Similar data types are defined in [13].

data type (GDF_TYP)	size per sample [bytes]	code	Range [min .. max], remarks
int8, INT-I8	1	1	-128 .. 127
uint8, INT-U8	1	2	0 .. 255
int16, INT-I16	2	3	-32 768 .. 32 767
uint16, INT-U16	2	4	0 .. 65 535
int32, INT-I32	4	5	-2 147 483 648 .. 2 147 483 647
uint32, INT-U32	4	6	0 .. 4 294 967 295
int64	8	7	$-2^{63} \dots 2^{63}-1$
uint64	8	8	$0 \dots 2^{64}-1$
float32,	4	16	According to [14]
float64	8	17	According to [14]
float128	16	18	According to [15]
int24, INT-I24	24	279	-8 388 608 .. 8 388 607
uint24, INT-U24	24	535	0 .. 16 777 215

25. The X-Y-Z coordinates of the EEG electrodes are included in the variable header. The coordinates of the references and ground electrodes are included in the fixed header from bytes 213-224 and 225-236, respectively. All positions should be stored in a coordinate system according to [9].
26. The electrode impedance will be stored in compressed form in 1 byte. For this purpose, the digital value will correspond to the impedance according to the following formulas.

$$\begin{aligned} \text{DIGVAL} &= \text{round}(\log_2(Z[\text{Ohm}]) * 8) \\ Z[\text{Ohm}] &= 2^{(\text{DIGVAL}/8)} \end{aligned}$$

This allows to store the value of the impedance with less than +5 % relative error. DIGVAL is 255, if the impedance is undefined or larger than $2^{(254/8)}=3.9 \text{ GOhm}$.

27. Beyond the fixed and variable header, it is possible to define a free header section. It has size $\text{HeadLen} - (\text{NS}+1)*256$. Currently, this section is ignored; this might change in future versions.
28. The data section contains the samples of the equidistant sampled channels. The samples of each channel and record are stored consecutively. First all samples of the first channel from the first record are stored, then all samples of channel 2 of the first block, until all channels of the first block are stored. Then, the samples of the second block are stored in the same order.
29. The table of events is stored after the data section. The starting position of the event table (event table position ETP) can be calculated in the following way.

$$\text{ETP} = \text{Length_of_header} + \text{Number_of_Records} * \text{Bytes_per_Record}.$$

The value of the first byte (mode of event table) can be 1 or 3. A value of 1 indicates that the event table contains event-type and position. A value of 3 indicates that the position, type,

associated channel and the duration of the event is stored.

30. Bytes 2 to 4 of the event table represent a 24-bit integer value indicating the number of events.
31. The next four bytes contain a floating point number representing the sample rate associated with the event position. Dividing the event position (and duration) by this value yields the position (duration) in seconds.
32. Then the position of all events is saved in 32-bit integers, followed by the event type as 16-bit integers. The encoding of the various event types is defined in Table 8 (see Annex). If the mode of the event table has value 3, it follows the channel number associated with each event (a value of 0 indicates the event refers to all channels); and the duration or the sample value. If the channel number refers to a non-equidistant (sparsely) sampled channel and the event type is 0x7fff, instead of the duration the value of a the sparsely sampled channel is represented. The the scaling (Physical and digital minimum and maximum) applies to the sample value.

3. RECOMMENDATIONS AND REMARKS:

GDF v.1 defined also data types of bit lengths e.g. bit 1, bit7, bit 4 etc. However, this has not been much of use. Furthermore, it complicated software development. Therefore, the use of bit-types is discouraged in GDF v.2.

The event table is defined only, if the number of records is known. Consequently, no eventtable can be stored in an ongoing recording, because the number of records is usually not known (NRec = -1). The eventtable can be only written to the GDF file, once the recording length (i.e. data size, number of records) is known (often only after the end of the data acquisition).

Non-equidistant sampled data can be stored, too. For this purpose, an additional channel header must be defined, and the corresponding sampling rate must be set to zero. The actual samples are stored in the event table, and the channel information must contain the corresponding channel number, then the 32 bit duration field is used for the actual sample value. The data type must correspond to the value of GDFTYP of the respective channel, and must not contain more than 32 bits. Preferably, the data type should be of type uint32.

The format definition of GDF is nearly as simple as the definition of EDF [4], as can be seen in Table 5. The use of binary encodings enables a more compressed representation; accordingly, more information can be stored within the header information. This enables a higher accuracy (e.g. in date and time information) and additional information can be stored without extending the header size.

The proposed format specification was successfully implemented in C/C++ as well as an M-file which can be used with Octave (>2.1.36) and Matlab® (>5.1). The software implementation requires only minor changes to upgrade from EDF to GDF. The software is available "online" (see BioSig – an open source software library for biomedical signal processing [1]) and is "free" under the terms of the "General Public License" (GPL) [7]. BioSig demonstrates also how different data formats can be simultaneously supported, without additional workload for the user.

4. DISCUSSION AND CONCLUSION

It is reasonable to say that GDF provides a superset of features from many other data formats. Accordingly, GDF 2.0 is (upwards) compatible to most other data formats; this means most other data can be converted to GDF without loss of any information. Routines for reading and writing of

GDF files in Octave and Matlab as well as in C are implemented in the open source package BioSig [1].

5. ACKNOWLEDGMENT

I thank (in alphabetical order) Oliver Filz, Arnulf Heller, Peter Jacobi, Claudia Keinrath, Herbert Ramoser, Stefan Sauermann, Reinhold Scherer and Darren Weber for their contributions to the discussion.

6. REFERENCES

- [1] A. Schlögl, The BioSig project. 2003-2006. online available at: <http://biosig.sf.net/> .
- [2] Overview of different data formats: <http://hci.tugraz.at/schloegl/matlab/eeg/>
- [3] A. Schlögl, O. Filz, H. Ramoser, G. Pfurtscheller, GDF - A general dataformat for biosignals, Technical Report, 2004. available online at:
http://hci.tugraz.at/schloegl/~schloegl/matlab/eeg/gdf4/TR_GDF.pdf
- [4] B. Kemp, A. Värri, A.C. Rosa, K.D. Nielsen and J. Gade (1992): A simple format for exchange of digitized polygraphic recordings. *Electroenceph. clin. Neurophysiol.*, 82: 391-393.
- [5] Kemp B, Olivan J. European data format 'plus' (EDF+), an EDF alike standard format for the exchange of physiological data. *Clin Neurophysiol.* 2003 Sep;114(9):1755-61.
- [6] A. Schlögl, B. Kemp, T. Penzel, D. Kunz, S.-L. Himanen, A. Värri, G. Dorffner, G. Pfurtscheller. Quality Control of polysomnographic Sleep Data by Histogram and Entropy Analysis. *Electroenceph. clin. Neurophysiol.* 1999, Dec; 110(12): 2165 - 2170.
- [7] The General Purpose License <http://www.fsf.org/>
- [8] RFC1876 How latitude and longitude are stored in a [DNS](http://www.faq.s.org/rfcs/rfc1876.html) record. available online:
<http://www.faq.s.org/rfcs/rfc1876.html>
- [9] J.A. Malmivuo. Consistent System of Rectangular and Spherical Coordinates for Electrocardiology and Magnetocardiology. online available:
<http://butler.cc.tut.fi/~malmivuo/bem/bembook/aa/aa.htm>.
- [10] SCP-ECG, Health informatics - Standard Communication Protocol - Computer-assisted electrocardiography, CEN TC 251/N02-015 prEN1064:2002. CEN/TC251 Secretariat: SIS - Swedish Standards Institute.
- [11] ISO/IEEE 11073-10101, Health informatics - point-of-care medical device communication Part 10101: nomenclature.
- [12] CEN/TC251 PT40, File Exchange Format of Vital Signal Annex A (Normative).
- [13] ISO/IEEE 11073-20101, Health informatics - point-of-care medical device communication Part 20101: Application profiles - Base standard. Appendix A, Table A.1. p. 19.

[14] IEEE Standard for Binary Floating-Point Arithmetic (ANSI/IEEE Std 754-1985), also known as IEC 60559:1989, Binary floating-point arithmetic for microprocessor systems.

[15] Revising ANSI/IEEE Std 754-1985, <http://754r.ucbtest.org>

Table 7: Overview of the header definition for GDF 2.0. The “remark” column indicates which fields have changed from Ver 1.25 to Ver 2.0.

FIXED HEADER	Remark 1.x -> 2.0	Position Start:End [bytes]	Bytes	Type
Version identification (GDF 2.00)	updated	0	8	char[8]
Patient identification (P-id)	smaller	8	66	char[80]
reserved		74	10	uint8
Smoking / Alcohol abuse / drug abuse / medication	Added	84	1	bit[4x2]
Weight [in kg]	Added	85	1	uint8
Height [in cm]	Added	86	1	uint8
Gender [bits 0..1]	Added	87	2/8	bit[2]
Handedness [bits 2..3]	Added	87	2/8	bit[2]
Visual impairment [bits 4..5]	Added	87	2/8	bit[2]
Recording identification (R-id)	smaller	88	64	char[64]
Recording Location (Lat, Long, Alt)	Added	152	16	uint32[4]
Startdate and -time of recording	changed	168	8	uint32[2]
Birthday	Added	176	8	uint32[2]
Header Length (number of 256-byte blocks)	smaller	184	2	uint16
reserved		186	6	byte[6]
Equipment Provider identification (EP-id)		192	8	uint64
Technician ID	removed		8	
Lab ID	removed		8	
IP address	Added	200	6	byte[6]
Headsize [in mm]	Added	206	6	uint16[3]
Position Reference Electrode [X,Y,Z]	Added	212	12	float32[3]
Position Ground Electrode [X,Y,Z]	Added	224	12	float32[3]
reserved	removed	236	20	byte
number of data records (-1 if unknown)		236	8	int64
Duration of a data record, as a rational number in seconds (first the numerator, secondly the denominator.		244	8	uint32[2]
NS: number of signals (channels)		252	2	uint16
reserved		254	2	uint16
VARIABLE HEADER		256		
Label		256	NS*16	char[16]*NS
Transducer type		256+16*NS	NS*80	char[80]*NS
Physical dimension	obsolete	256+96*NS	NS*6	char[8]*NS
Physical dimension Code	NEW	256+102*NS	NS*2	uint16*NS
Physical minimum		256+104*NS	NS*8	float64*NS
Physical maximum		256+112*NS	NS*8	float64*NS
digital minimum	changed	256+120*NS	NS*8	float64*NS
digital maximum	changed	256+128*NS	NS*8	float64*NS
Pre-filtering	obsolete	256+136*NS	NS*68	char[68]*NS
Lowpass	Added	256+204*NS	NS*4	float32*NS
Highpass	added	256+208*NS	NS*4	float32*NS
Notch	added	256+212*NS	NS*4	float32*NS
Samples Per Record (0 indicates non-equidistant sampling)		256+216*NS	NS*4	uint32*NS
Channel TYPE		256+220*NS	NS*4	uint32*NS
Electrode Positions XYZ	added	256+224*NS	NS*12	float32*NS*3
Electrode Impedance	added	256+236*NS	NS*1	uint8*NS
reserved	smaller	256+237*NS	NS*19	char[32]*NS

FREE HEADER [optional]				
Any information in this section is ignored. This might change in future versions.		$(NS+1) * 256 :$ Header Length-1	Header Length – $(NS+1) * 256$	ASCII (char)
DATA RECORD		$256*(NS+1)$		
nr samples from channel [1] of TYPE[1]				Type[1]
nr samples from channel [2] of TYPE[2]				Type[2]
nr samples from channel [3] of TYPE[3]				Type[3]
nr samples from channel [NS] of TYPE[NS]				Type[NS]
EVENT TABLE	NEW	ETP		
mode of event table can be 1 or 3.		ETP + 0	1	uint8
Number of events NEV	Changed	ETP + 1	3	uint8[3]
Samplerate associated with Event positions.	Changed	ETP + 4	4	float32
Position [in samples]		ETP + 8	NEV*4	uint32
Type		ETP+ 8+NEV*4	NEV*2	uint16
Channel [optional]	Only if mode is 3	ETP+ 8+NEV*6	NEV*2	uint16
Duration [in samples, optional] or Value of non-equidistant sampling	Only if mode is 3	ETP+ 8+NEV*8	NEV*4	uint32 float32(*)

*Table 8: Table of event codes. The most recent version of this table will be available from [1] http://biosig.cvs.sourceforge.net/*checkout*/biosig/biosig/doc/eventcodes.txt .*

```

### Table of event codes.
# This file is part of the biosig project http://biosig.sf.net/
# Copyright (C) 2004 Alois Schloegl <a.schloegl@ieee.org>
# $Revision: 1.3 $
# $Id: eventcodes.txt,v 1.3 2004/06/17 17:08:54 schloegl Exp $
#
### table of event codes: lines starting with # are omitted
### add 0x8000 to indicate end of event
#
### 0x010_      EEG artifacts
0x0101 artifact:EOG
0x0102 artifact:ECG
0x0103 artifact:EMG/Muscle
0x0104 artifact:Movement
0x0105 artifact:Failing Electrode
0x0106 artifact:Sweat
0x0107 artifact:50/60 Hz mains interference
0x0108 artifact:breathing
0x0109 artifact:pulse
### 0x011_      EEG patterns
0x0111 eeg:Sleep spindles
0x0112 eeg:K-complexes
0x0113 eeg:Saw-tooth waves
### 0x03__      Trigger, cues, classlabels,
0x0300 Trigger, start of Trial (unspecific)
0x0301 Left - cue onset (BCI experiment)
0x0302 Right - cue onset (BCI experiment)
0x0303 Foot - cue onset (BCI experiment)
0x0304 Tongue - cue onset (BCI experiment)
0x0306 Down - cue onset (BCI experiment)
0x030C Up - cue onset (BCI experiment)
0x030D Feedback (continuous) - onset (BCI experiment)
0x030E Feedback (discrete) - onset (BCI experiment)
0x0311 Beep (accustic stimulus, BCI experiment)
0x0312 Cross on screen (BCI experiment)
0x03ff Rejection of whole trial
### 0x040_      Sleep-related Respiratory Events
0x0401 Obstructive Apnea/Hypopnea Event (OAHE)
0x0402 Respiratory Effort Related Arousal (RERA)
0x0403 Central Apnea/Hypopnea Event (CAHE)
0x0404 Cheyne-Stokes Breathing (CSB)
0x0405 Sleep Hypoventilation
### 0x041_      Sleep stages according to Rechtschaffen&Kales
0x0410 Wake
0x0411 Stage 1
0x0412 Stage 2
0x0413 Stage 3
0x0414 Stage 4
0x0415 REM
### 0x050_      ECG events
0x0501 ecg:Fiducial point of QRS complex
0x0502 ecg:P-wave
0x0503 ecg:Q-point
0x0504 ecg:R-point
0x0505 ecg:S-point
0x0506 ecg:T-point
0x0507 ecg:U-wave
### 0x_____    OTHER
0x0000 No event
0x7FFF non-equidistant sampled value

```